

Titer determination of strong bases

Description

This application report describes the general procedure for the titer determination of aqueous strong bases like Sodium hydroxide and Potassium hydroxide. The procedure is also usable for strong bases like Potassium hydroxide and Tetrabutylammonium hydroxide in water-soluble organic solvents like alcohols. This procedure is usable for base solutions up to 3 mol/l.

The titer is a dimensionless number about 1 for correcting the indicated concentration. In the software of the titration devices and application reports from SI Analytics®, the term "Titer" describes the exact concentration in mol/l and not the dimensionless factor.

Instruments

Titration	TL 7000 or higher
Exchangeable Unit	WA 20
Electrode	N 62 or A 7780 1M-DIN-ID or similar
Cable	L 1 A (only for electrodes with plug head)
Stirrer	Magnetic stirrer TM 235 or similar
Lab accessoires	Glas beaker 150 ml
	Magnetic stirrer bar 30 mm

Reagents

1	the base from which the titer is to be determined
2	Distilled water
3	Potassiumhydrogenphthalate (KHPht) – certified reference material, volumetric standard
4	KCl solution 3 mol/l
5	Soda lime (a mixture of solid NaOH and Ca(OH) ₂)
All reagents should be in analytical grade or better.	

Titration procedure

Reagents

The KHPht volumetric standard is dried as described in the corresponding certificate of analysis, mostly 3h at 105°C.

The base solutions can absorb CO₂ from the atmosphere and thus alter their content. The reagent bottles with alkaline reagents are closed with CO₂ absorption tubes. These contain soda lime, a mixture of solid NaOH and Ca(OH)₂. This must be exchanged regularly. An indicator in soda lime indicates the exchange too late.

For very accurate titrations with very low concentrated bases they can be used under inert gas like nitrogen.

Cleaning and storage of the electrode

Use distilled water for cleaning the electrode. For storage use KCl solution 3 mol/l or electrolyte solution L 911.

Sample preparation

The amount of volumetric standard depends on the size of the burette and the concentration of the base. The amount should be chosen so that about half of the burette volume is consumed. The most common is the 20 ml burette. The required quantity of TRIS can be estimated according to this rule of thumb:


$$W [g] = 2 * Concentration[mol/l]$$

To determine the titer of a 0.1 mol/l base, about 0.2 g KHPht volumetric standard are weighed into a 150 ml beaker with an accuracy of 0.1 mg and filled up to 80 ml with distilled water. When the KHPht is completely dissolved, it is titrated with the base to an EQ.

If the specified assay of the volumetric standard is significantly different from 100%, the weight for calculating the concentration must be corrected:

$$W = \frac{Weight * specified assay \%}{100}$$

Titration parameter



Default method	Titre NaOH		
Method type	Automatic titration		
Modus	Dynamic		
Measured value	pH		
Measuring speed / drift	Normal	Minimum holding time	2 s
		Maximum holding time	15 s
		Measuring time	2 s
		Drift	20 mV/min
Initial waiting time	0 s		
Dynamic	Steep	Max step size	1.0 ml
		Slope max ml	15
		Min. step size	0,02 ml
		Slope min. ml	230
Damping	None	Titration direction	Increase
Pretitration	off	Delay time	0 s
End value	12		
EQ	on	Slope value	700
Max. titration volume	20 ml		
Dosing speed	100%	Filling speed	30 s

Calculation:

$$T [\text{mol/l}] = \frac{W * F2}{(EQ - B) * M * F1}$$

B	0	Blank value
W	man	Weight of the sample [g]
F2	1000	Conversion factor ml - l
EQ1		Consumption of titrant until first Equivalence point
M	204,22	Molecular mass
F1	1	Conversion factor

We recommend to write the exact concentration T to the Exchangable Unit (WA) automatically.

Any questions? Please contact the application team:

Xylem Analytics Germany Sales GmbH & Co. KG, SI Analytics
Hattenbergstraße 10
D-55122 Mainz, Germany
Telefon: + 49 6131 66 5126
Fax: + 49 6131 66 5101
E-Mail: titration@si-analytics.com

SI Analytics
a xylem brand

Xylem Analytics Germany Sales GmbH & Co. KG · Hattenbergstr. 10 · D-55122 Mainz · Germany
Telefon: +49 6131.66. 5111 · E-Mail: Info.si-analytics@Xyleminc.com · www.si-analytics.com

Alle Namen sind eingetragene Handelsnamen oder Warenzeichen der Xylem Inc. oder eines seiner Tochterunternehmen. Technische Änderungen vorbehalten.
© 2018 Xylem Analytics Germany Sales GmbH & Co. KG.