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8157916_GA_ProzesspHredox_211208_D

Instruction manual for process electrodes disposed to measure pH-values and redox (ORP) potentials

The electrode shown is an example type

Also available with integrated temperature sensor

Xylem Analytics Germany GmbH
Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49(0)6131.66.5119
Fax: +49(0)6131.66.5001
E-mail: Service-Sensors.si-analytics@xylem.com
www.XylemAnalytics.com

Preparation and general

The process electrodes are supplied ready for use. If there is a irrigation cap over membrane/sensor and diaphragm, it must be removed for measurement and calibration. It contains irrigation fluid. The electrode is ready to perform a measurement. Electrodes stored dry must be irrigated for 24 hours in irrigation fluid and then checked. Lacking electrolyte solution has to be refilled.

Electrodes with liquid electrolyte should always be connected to a vessel with electrolyte solution. During measurements without pressure it is sufficient to place the vessel 1 m above the location of measuring point. With this hydrostatic pressure it is guaranteed that enough electrolyte solution is flowing out through the diaphragm in order to protect the electrode. If the measurement is performed under pressure the vessel should be charged with an excess pressure of 0.5 bar above the measuring pressure. Refilling is not required for low-maintenance electrodes with gel filling or the REFERID® system. Irrigation with electrolyte solution is however particularly important for these electrodes.

Measuring the pH value

Please also note the usage instructions of the measuring device during calibration and measurement. In order to minimize inaccuracies in the measurement results, electrodes used under extreme conditions or at the limits of the specified usage range should be calibrated more frequently. For accurate calibration, we recommend the use of our hot steam sterilized, certified buffer ampoules to DIN 19 266. Always use fresh buffer solution.

Measuring the Redox voltage

In the case of metal single-rod measurement chains, an Ag/AgCl reference system is used as standard. A calibration of these electrodes isn't performed. Redox test solutions are available in order to check them.

Storage and maintenance

Electrodes should be stored at a temperature of between 0° and 40 °C. Depending on storage conditions (temperature and air humidity), the irrigation fluid in the cap may dry out prematurely. In this case, the electrode must be irrigated for at least 24 hours in potassium chloride solution $c(\text{KCl}) = 3 \text{ mol/l}$ before it is ready for use.

In the case of combination electrodes (pH or redox) as well as reference electrodes, the reference electrolyte must be refilled or replaced occasionally. Crystals in the electrolyte area can be dissolved by warming in a water bath. The electrolyte solution should then be replaced, and the correct potassium chloride concentration checked.

Cleaning

Dirt and contamination at the membrane/sensor and diaphragm will lead to measurement inaccuracies:

- **Coatings and deposits** can be removed with diluted mineral assets (e.g. hydrochloric acid 1:1), **Organic** contamination with the aid of suitable solvents, **Grease** with tenside solution and **Proteins** with hydrochloric pepsin solution (cleaning solution L 510).
- After cleaning, rinse off the electrodes with distilled water, do not rub dry.
- Ceramic diaphragms blocked from outside can be restored to working order by carefully rubbing down with fine sandpaper or a diamond file. **Take care not to scratch the pH glass membrane!**
- Platinum diaphragms must not be subjected to any mechanical treatment. Rinsing (e.g. vacuuming) may be followed by chemical cleaning (e.g. with diluted hydrochloric acid).
- Ground diaphragms are cleaned and prepared for measurement by gently raising and then replacing the ground sleeve on the core. The refilling hole must be open during this process. **Caution:** More electrolyte flows out during this process to ensure full coverage of the ground surface.
- The glass membrane can be cleaned by rubbing off with a fluff-free cloth soaked in ethanol.

Quality

Every electrode must meet the strict quality requirements of final testing. The durability depends mainly on the usage conditions. Extreme conditions include for example high or frequently fluctuating temperatures, strong acids and caustic solutions, protein and heavily contaminated solutions as sulphides, bromides and iodides. Hydrofluoric acid and hot phosphoric acid corrode glass.

Further information

Further information can be found in our meters and electrodes catalog.
Subject to technical amendment.

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