**GEBRAUCHSANLEITUNG**  
Originalversion  
**OPERATING MANUAL**  
**MODE D'EMPLOI**  
**MANUAL DE INSTRUCCIONES**

# TITRONIC<sup>®</sup> 300

KOLBENBÜRETTE | PISTON BURETTE |  
BURETTE A PISTON | BURETA DE EMBOLO

**SI Analytics**

a **xylem** brand

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## **Operating Manual ..... Page 59 ... 114**

### **Important notes:**

The operating manual is part of the product. Before initial operation, please carefully read and observe the operating manual and keep it. For safety reasons the product may only be used for the purposes described in these present operating manual. Please also consider the operating manuals for the devices to be connected.

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties. A potentially more recent version of this manual is available on our internet website. The German version is the original version and binding in all specifications!

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## TABLE OF CONTENT

<b>1</b>	<b>Technical Specifications of the TITRONIC® 300 piston burette.....</b>	<b>61</b>
1.1	Notes to the operating manual .....	61
1.2	Intended Use .....	61
1.3	Technical Specifications.....	62
1.3.1	Piston burette TITRONIC® 300 .....	62
1.4	Warning and safety information .....	64
1.4.1	General.....	64
1.4.2	Chemical and biological safety.....	65
1.4.3	Flammable liquids .....	65
<b>2</b>	<b>Installation and Commissioning.....</b>	<b>66</b>
2.1	Unpacking and setting up.....	66
2.2	Back panel of the piston burette TITRONIC® 300.....	66
2.3	Connection and installation of the piston burette and the magnetic stirrer TM 50.....	67
2.4	Setting the language .....	68
2.6	Dosing unit and Accessories.....	69
2.7	Installing the burette tip .....	70
2.7.1	Initial Filling or Rinsing of the Entire Interchangeable Unit .....	71
<b>3</b>	<b>Working with the piston burette TITRONIC® 300 .....</b>	<b>73</b>
3.1	Front Keyboard.....	73
3.2	Display.....	73
3.3	Manual controller.....	74
3.4	External PC Keyboard.....	74
3.5	Menu Structure .....	75
3.6	Main Menu.....	77
3.6.1	Manual Titration.....	77
3.6.2	Dosage .....	80
<b>4</b>	<b>Method parameters.....</b>	<b>82</b>
4.1	Method editing and new method .....	82
4.2	Default method .....	83
4.3	Copy Method .....	83
4.4	Delete Method .....	84
4.5	Change Method Parameters .....	84
4.5.1	Method type.....	84
4.5.2	Result .....	85
4.5.3	Dosing parameter.....	91
4.5.4	Sample identification .....	92
4.5.5	Documentation .....	93
<b>5</b>	<b>System settings .....</b>	<b>94</b>
5.1	Setting the language .....	94
5.2	Dosing unit - Reagents.....	94
5.2.1	Replacing the dosing unit.....	95
5.2.2	Reagent data.....	98
5.3	Globale Memory .....	99
5.4	RS-232 Settings .....	99
5.5	Date and Time.....	101
5.6	RESET .....	101
5.7	Printer .....	102
5.8	Device Information .....	102
5.9	System Tone .....	103
5.10	Data exchange .....	103
5.11	Software Update.....	105
<b>6</b>	<b>Communication via RS-232 and USB-B interface .....</b>	<b>107</b>
6.1	General Information.....	107
6.2	Chaining multiple devices - "Daisy Chain Concept" .....	107
6.3	Instruction Set for RS-Communication.....	107

<b>7</b>	<b>Connection of Analytical Balances and Printers .....</b>	<b>109</b>
7.1	Connection of Analytical Balances .....	109
7.2	Balance data editor .....	110
7.3	Printers .....	111
<b>8</b>	<b>Maintenance and Care of the piston burette .....</b>	<b>112</b>
<b>9</b>	<b>Guarantee .....</b>	<b>113</b>
<b>10</b>	<b>Storage and transportation.....</b>	<b>113</b>
<b>11</b>	<b>Recycling and Disposal .....</b>	<b>113</b>
<b>12</b>	<b>EC – Declaration of Conformity.....</b>	<b>113</b>

---

**Copyright**

© 2021, Xylem Analytics Germany GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization.

Germany. Printed in Germany.

# 1 Technical Specifications of the TITRONIC<sup>®</sup> 300 piston burette

## 1.1 Notes to the operating manual

The provided operating manual will allow you the proper and safe handling of the product. For maximum security, observe the safety and warning instructions in the operating manual!

-  Warning of a general danger:  
Non-compliance results (can result) in injury or material damage.
-  Important information for device use.
-  Refers to another part of the operating manual.

The menu screens shown in this operating manual serve as an example and may differ from what you see!

## 1.2 Intended Use

The TITRONIC<sup>®</sup> 300 is a piston burette and suitable for the following applications:

- Manual titrations with or without calculation of the result
- Dosing
- It can be used as a dosing or titrating burette with the TitroLine<sup>®</sup> 7000, 7750, 7800 titrator
- It can be used as a dosing or titrating burette in combination with the TitriSoft control software from version 3.1

Each method allows for the setting of a variety of dosing and filling rates.

Up to 3 user methods can be memorised in the device.

### Solutions to be used:

Virtually, any liquids and solutions with a viscosity of  $\leq 10 \text{ mm}^2/\text{s}$  such as concentrated sulphuric acid may be used.

 However, one has to avoid the use of chemicals that may attack glass, PTFE or FEP or that are explosive, such as hydrofluoric acid, sodium azide or bromine! Suspensions containing high solids percentages may clog or even damage the dosing system.

 **Do not use the device in hazardous locations!**

 **General:**

The safety guidelines that are applicable to the handling of chemicals have to be observed under all circumstances. This applies in particular to inflammable and/or etching liquids.

## 1.3 Technical Specifications

### 1.3.1 Piston burette TITRONIC® 300

Translation of the legally binding German version

(Release: 07. July 2020)

EMC compatibility according to the Council Directive: 2014/30/EU;  
 applied harmonized standards: EN 61326-1  
 Low-voltage directive according to the Council Directive 2014/35/EU;  
 Testing basis EN 61 010-1: for laboratory equipment  
 RoHS Council Directive 2011/65/EU  
 FCC Part 15B and ICES 003

**Country of origin:** Germany, Made in Germany

**The following solvents/titration reagents are allowed to be used:**

- All common titration solutions.
- As reagent water and all non-aggressive non-organic and organic fluids are allowed.
- If using combustible fluids fire please adhere to the Guidelines for Explosion Protection and Prevention of the chemical industry.
- For fluids with higher viscosity ( $\geq 5 \text{ mm}^2/\text{s}$ ), lower boiling point or affinity to outgas, the filling and dosage speed can be adjusted.
- Fluids with viscosity over  $20 \text{ mm}^2/\text{s}$  cannot be dosed.

**Display:** 3.5 inches -1/4 VGA TFT display with 320x240 pixels

**Power supply:** by external multi-range power supply from 100 – 240 V, 50/60 Hz  
 Input voltage: 12 Volt DC, 2500 mA  
 Power consumption 30 W  
 Corresponds to protection class III:  
 Protection class for dust and humidity IP 50 according to DIN 40 050

**⚠ Only use the power supply TZ 1853 or a power supply approved by the manufacturer!**

**RS-232-C Interface:**

Daisy Chain function available

Data bits: adjustable, 7 or **8** Bit (default: 8 Bit)

Stop bit: adjustable, **1** or 2 Bit (default: 1 Bit)

Start bit: static **1** Bit

Parity: adjustable: even / odd / **none**

Baud rate: adjustable: 1200, 2400, **4800**, 9600, 19200 (Default 4800 baud)

Address: adjustable, (0 to 15, default: 01)

RS-232-1 for computer, input Daisy Chain

RS-232-2 devices of SI Analytics®:  
 - Titrator TitroLine® 7000 / 7500 / 7500 KF / 7750 / 7800  
 - Sample Changer TW alpha plus, TW 7400  
 - Piston burette TITRONIC® 300 and 500, TITRONIC® 110 *plus*, TITRONIC® *universal*,  
 - Balances of the types Mettler, Sartorius, Kern, Ohaus, (for more, please contact us)  
 - Exit Daisy-Chain

**USB Interface:** 1 x USB-type-A and 1 x USB-type-B

USB-type A for connecting of USB keyboard, - printer, - manual controller, - data media (e.g. USB stick) and USB-Hub

USB-type B for connecting a PC

**Stirrer/pump:** Plug connection with integrated low-voltage supply (15 V  $\approx$ ) for the TM 50 Stirrer incorporated in the bottom of the casing of the TITRONIC<sup>®</sup> 300 piston burette

**Housing:**

Material: Polypropylene

Front keyboard: polyester coated

Dimensions: 13.5 x 31 x 20.5 cm (W x H x D), height incl. interchangeable unit

Weight: approx. 2.0 kg for basic unit

**Ambient conditions:**

 **Do not use the device in hazardous locations!**

Climate: Ambient temperature: + 10 ... + 40 °C for operation and storage  
Humidity according to EN 61 010, Part 1:  
Max. relative humidity 80 % for temperatures up to 31 °C,  
linear decrease down to 50 % relative humidity at a temperature of 40 °C

Altitude: Device: No restrictions  
Power supply: up to 5000 m

Pollution degree:  
Pollution degree IP 20, indoor use only.

**Dosing units:**

Cylinder: 20 ml and 50 ml, borosilicate glass 3.3 (DURAN<sup>®</sup>)

Valve: volume neutral cone valve made from fluorocarbon polymers (PTFE), TZ 3000

Hoses: FEP hose set, blue

Dosing accuracy:  
after DIN EN ISO 8655, part 3:  
Accuracy: 0.15 %  
Precision: 0.05 %

## 1.4 Warning and safety information

### 1.4.1 General

The device corresponds to protection class III.

It was manufactured and tested according to DIN EN 61 010, Part 1, “**Protective Measures for electronic measurement devices**” and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information contained in the present operating instructions. Development and production is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard.

 For reasons of safety, the device must only be used for the range of application described in the present operating manual. Nonobservance of the intended proper use of the device may result in personal injury or damage to property.

 For reasons of safety, the device and the power supply must be opened by authorised persons only; this means, for instance, that work on electrical equipment must only be performed by qualified specialists. **In case of nonobservance of these provisions the titrator and the power supply may constitute a danger: electrical accidents of persons or fire hazard!** Moreover, in the case of unauthorised intervention in the titrator or the power supply, as well as in the case of negligently or deliberately caused damage, the warranty will become void.

 Prior to switching the device on it has to be ensured that the operating voltage matches the mains voltage. The operating voltage is indicated on the specification plate (underside of the device and backside of the power supply). **Nonobservance of this provision may result in damage to the titrator and the power supply, or in personal injury or damage to property!**

 **If it has to be assumed that safe operation is impossible, the device has to be put out of operation and secured against inadvertent putting to operation.** In this case please switch the device off, pull plug of the mains cable out of the power supply, and remove the device from the place of work.

Examples for the assumption that a safe operation is no longer possible,

- if the package is damaged,
- if the device shows visible damages,
- if the power supply shows visible damages,
- if the device does not function properly,
- if liquid has penetrated into the casing.
- if the unit has been altered technologically or if unauthorized personnel tried or succeeded to open the device as attempt to repair it.

In case that the user operates such a device, all thereof resulting risks are on the user!

 The device must not be stored or operated in humid rooms.

 **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the environment or that they are supervised by specialized persons, respectively.

 For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

 The device is equipped with integrated circuits (EPROMs). X rays or other high energy radiation may penetrate through the device's casing and delete the program.

 For working with liquids, not being common titration solvents, especially the chemical resistance of the construction materials of the device have to be considered (see  1.3 Technical Specifications).

 For the use of liquids with high vapour pressure or (mixture of) substances not being mentioned in  1.3 Technical Specifications) as allowed substances, the safe and proper operation of the device has to be guaranteed by the user. When the piston moves upwards within the cylinder, a microfilm of dosing liquid or titration solution will always remain adhered to the inner wall of the cylinder, but this has no influence on the dosing accuracy. This small residue of liquid, however, may evaporate and thus penetrate into the zone

underneath the piston, and if non-admitted liquids are being used, the materials of the may be dissolved or corroded (see  8 Maintenance and Care of the )

### 1.4.2 Chemical and biological safety

 The device is not intended for use with potentially biohazardous substances.

 **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the devices or that they are supervised by specialized persons, respectively.

 When using biohazardous substances, the regulations for handling the substances used must be observed. In such cases, the use is the sole responsibility of the user.

 For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

 Dispose of all used solutions in accordance with national regulations and laws. Select the type of protective equipment according to the concentration and quantity of the hazardous substance at the respective workplace.

### 1.4.3 Flammable liquids

When handling flammable liquids, make sure that there is no naked flame in the vicinity of the equipment. Adequate ventilation must be provided. Only small quantities of flammable liquids should be kept in the workplace.

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## 2 Installation and Commissioning


### 2.1 Unpacking and setting up

The device has been put together especially for you (basic unit + corresponding modules and accessories), there may be differences with respect to the delivery and the accessories described in this chapter. The scope of delivery, please refer to the attached packing list. For any questions please contact us directly (see backside of this operating manual).

The device itself as well as all related accessory and peripheral parts have been carefully checked at the factory to ensure their correct function and size. Please ensure that the small accessories are also removed in full from the packaging.

The device may be placed on any flat surface.

### 2.2 Back panel of the piston burette TITRONIC® 300

**Fig. 1**

The TITRONIC® 300 is equipped with the following connections:

Two RS-232 ports, 4-channel (Mini-DIN):

- 1) RS2 for connection of a weighing balance and other devices from SI Analytics®
- 2) RS1 for connection to the PC
- 3) USB-A interfaces for connecting USB devices
- 4) USB-B interface for connection to a PC
- 5) Socket "in": Connection of the external power supply TZ 1853
- 6) On/Off switch

### 2.3 Connection and installation of the piston burette and the magnetic stirrer TM 50


The low voltage cable of the power supply TZ 1853 has to be plugged in to the 12 V socket “in”, on the back panel of the device (Fig. 2). Then plug the power supply into the plug socket.

**Fig. 2**

**⚠** Place the power supply easily accessible in order to be able to remove the device anytime easily from the power circuit.

The stirrer connects to the right side at the bottom and is locked in position by pushing it backwards (Fig. 3). In this way the power supply of the TM 50 stirrer is automatically established.

**Fig. 3**

The stand rod TZ 1748 is screwed into the thread and the titration clamp Z 305 may now be mounted on the stand rod (Fig. 4). Instead of the magnetic stirrer TM 50, also the titration stand without stirring function (TZ 3886) can be connected.

**Fig. 4**

## 2.4 Setting the language

The ex-factory default language setting is English.

After the device is switched on and the start-up process is complete, the main menu appears (Fig. 5).


Fig. 5

Using <SYS> or <MODE>, you navigate to the system settings («System settings»). The very first menu is to be used for setting the language (Fig. 6).


Fig. 6

Use <ENTER>/<OK> to call the menu. Select the language using the arrow keys <↑↓>. Confirm with <ENTER>/<OK>.



Fig. 7

The selected language will appear immediately (Fig. 7). Pressing <ESC> twice will return the user to the main menu.

## 2.6 Dosing unit and Accessories

**Fig. 8**

- 7) TZ 2003 - drying tube
- 8) TZ 3282 - dosing hose without dosing tip and holding bracket, or  
TZ 3620 - dosing hose with dosing tip and holding bracket
- 9) TZ 3802 - threaded cap with borehole GL 45, incl. adapter with 2 openings for drying tube and suction hose
- 10) TZ 1748 - Stand rod
- 11) Z 305 - Titration clamp
- 12) TZ 3803 - 1 litre reagent bottle, brown
- 13) TZ 3656 - Titration tip unit
- 14) TZ 3801 - Valve cover lid and TZ 3000 - 3/2-way valve
- 15) TZ 3130 - 20 ml dosing unit or  
TZ 3160 - 50 ml dosing unit
- 16) TZ 3283 - Connection hose
- 17) TZ 3281 - Suction hose

## 2.7 Installing the burette tip

The burette tip consists of the elements shaft with threaded clamping joint, hose and slip-on tip (



Fig. 9).

**Fig. 9**


Burette tip - Sequence of assembly:

1. Cut of hose end evenly.
2. Slip parts of the threaded clamping joint on to the hose.
3. Guide hose through shaft.
4. Press the free hose end over the ring marking until it reaches the stop of the tip.
5. Push the tip with pressed in hose onto the shaft.
6. Hold tip firmly, and screw threaded clamping joint to the shaft

### 2.7.1 Initial Filling or Rinsing of the Entire Interchangeable Unit


**⚠** While the initial filling or rinsing programme is being run, please place a sufficiently dimensioned waste vessel under the titration tip.

Initial filling of the interchangeable unit is done using the «**rinsing**» program.

**Fig. 10**


On the main menu (Fig. 10), press <**MODE**> to navigate to the methods/system. Pressing <↑> twice will take you to the «**Rinsing**» selection immediately (Fig. 11).

**Fig. 11**

Confirm the selection by pressing <**ENTER**>/<**OK**>. At this point you can select the number of rinsing cycles (Fig. 12).

**i** Initial filling requires a minimum of two rinsing cycles!

**Fig. 12**

Confirm the selection by pressing <**ENTER**>/<**OK**>. At this point you can select the number of rinsing cycles (Fig. 13 - Fig. 16).


Fig. 13

Fig. 14

Fig. 15


Fig. 16

**i** You can stop the rinsing operation at any time by pressing <STOP> and then resume rinsing with <START>. When the rinsing is finished, you can get back to the start menu by pushing 2 x <ESC>.

## 3 Working with the piston burette TITRONIC® 300

### 3.1 Front Keyboard

Fig. 17

**i** Apart from alphanumeric input (a-z, A-Z, 0-9) and a few other functions, almost all functions can be performed using the front keyboard (Fig. 17).

<Mode>:	Methods selection, rinsing, system settings
<EDIT>:	Changing the current method, new method, copy and delete method
<ESC>:	<ESC> will take you back to the previous menu level
<START>:	Start and Stop of a current method
<FILL>:	Filling the unit
<↑>:	Arrow-key up: Choosing of a menu or changing a number
<↓>:	Arrow-key down: Choosing of a menu or changing a number
<→>:	Arrow-key right: Changing position of the cursor, Enter

The individual functions are described in detail in [3.4 External PC Keyboard](#).

### 3.2 Display

The display (Fig. 18) consists of a graphical LCD display with a resolution of 320 x 240 pixels.


Fig. 18

### 3.3 Manual controller

The manual controller (Fig. 19) is needed for manual titration. It can also be used for starting dosage or other methods.

Fig. 19

Mode	Black key	Grey Key
Manual titration	Start of titration, single-step and continuous titration	Filling Stop of titration including evaluation
Dosage through Dosage method	Start dosage	Filling

### 3.4 External PC Keyboard

Keys	Function
<ESC>	<ESC> will take the user to the previous level on the menu
<F1>/<START>	Start of a selected method
<F2>/<STOP>	Stop of the current method
<F3>/<EDIT>	Change of the current method, new method, copy method
<F4>/<FILL>	Fill the interchangeable unit
<F5>/ 	Display and modification of the balance data. With <Shift + F5> display and modification of the global memories
<F6>/<MODE>	Selection of method, rinsing, system settings
<F7>/<SYS>	System settings (language selection, time/date ...)
<F8>/<CAL>	Without function on TITRONIC® 300
<F9>/+/-	Change of sign
<F10>/<DOS>	Start dosing menu
Num/ Scroll Lock/ Lock	Without function
Prt Sc Sys Rq	Without function
<↑> <↓> <←> <→>	Selection of individual menus and numeric values
0...9	Input of numeric values
<ENTER>	Confirmation of input parameters
< ←Backspace >	Deletion of one input digit / an input character to the left of the flashing cursor
Letters, ASCII-symbols	Alphanumeric input possible, Uppercase and lowercase possible
All other keys	Do not have any function


### 3.5 Menu Structure

**i** The menu screens shown in this manual serve as an example and may differ from what you see!

There are 4 selection menus:


- Start or main menu
- Method parameters
- Method selection
- System settings.

After power-up, the main menu is always the first menu to appear.  
The method displayed will always be the last method that was used (Fig. 20).

**Fig. 20**

Pressing **<START>** will result in the immediate execution of the method shown.  
**<EDIT>** will take you to the method parameters (Fig. 21).

**Fig. 21**

At this point you can

- modify the current method
- create a new method
- call and memorise standard methods
- copy or delete an existing method

Use **<↓>** and **<↑>** to select the submenus.  
Confirm your selection with **<ENTER>/<OK>**.  
**<ESC>** will take you back to the main menu.

<MODE>/F6 leads you to the select method menu (Fig. 22).


Fig. 22

Existing methods can be selected by pressing <↓> and <↑> and confirming the selection with <ENTER>/<OK>. Once the selection made, you will return to the main menu with the newly selected method. If no method is selected <ESC> will also take you back to the main menu.

To navigate directly to the system settings (Fig. 23 and Fig. 24) you can use the <SYS>/F7 key; you can also navigate there through the method selection menu.


Fig. 23


Fig. 24

### 3.6 Main Menu

After power-up, the main menu is always the first menu to appear. The method displayed will always be the last method that was used (Fig. 25).


Fig. 25

#### 3.6.1 Manual Titration

**i** Manual titration is impossible without the manual controller.


Fig. 26

<START> or pressing the black key on the manual controller will start the manual titration method.

Depending on the method settings, you will be prompted for the sample identification (Fig. 27) and the weighed-in quantity (Fig. 28). You can use an external PC keyboard for entering a 20-digit alphanumeric sample ID.


Fig. 27


Fig. 28

The balance data can be entered using the front keyboard or an external keyboard. The input is to be confirmed with <ENTER>/<OK>.

In the case of an automatic acceptance of the balance data, the weighed-in quantities will be read in from a memory. If the memory does not contain any balance data, a message will appear (Fig. 29).


Fig. 29

Pressing the Print key will transfer the balance data, too.

**i** Titration will then begin directly after the transfer of the balance data without any further confirmation being necessary.

Following the input of the sample description and/or the weight/volume the display changes (Fig. 30).



Fig. 30

You can control the metering rate with the black key of the manual controller (Fig. 30).

- a) A single depression of the key will cause a step up to the first level. Depending on the size of the interchangeable unit, this corresponds to 0.005 ml (20 ml) or 0.025 ml (50 ml).  
Displayed are 0.003 or 0.006 ml from a single step.
- b) If one keeps the black key depressed on the first level, titration will be continued at a low rate.
- c) If you press the black key fully down (2<sup>nd</sup> level) titration will proceed at a higher rate.

The rate of the second level can be set in five stages using the <↓↑> arrow keys.

**i** These stages can also be changed during manual titration (Fig. 31).

**Fig. 31**


Stage 5 corresponds to maximum titration speed. Speed is reduced by 50% each time.

**Example:**

20 ml dosing unit

Stage 5	100 %	approx.	40 ml/min
Stage 4	50 %	approx.	20 ml/min
Stage 3	25 %	approx.	10 ml/min
Stage 2	12.5 %	approx.	5 ml/min
Stage 1	6.8 %	approx.	2.5 ml/min

Even if the titration is completed, press <**STOP**> or approx. for 1 sec. the grey key of the manual controller. The titration result will be calculated and displayed (Fig. 32) and optionally printed on the connected printer or USB stick in PDF format.

**Fig. 32**

<**ESC**> will take you back to the start menu way to start the next titration immediately. Filling of the interchangeable unit occurs automatically.

### 3.6.2 Dosage

#### 3.6.2.1 Dosing operation with dosing method

Use <START> or the black key of the manual controller to start a dosage method (Fig. 33 and Fig. 34).

**Main menu**

**0.000 ml**  
2.000 ml

dosing NaOH  
Method parameter  
Select method / system  
20 ml

START  
EDIT  
MODE

07/20/14 19:51

Fig. 33

**Dose**  
dosing NaOH

**1.323 ml**  
2.000 ml

STOP      ESC

20 ml      07/20/14 19:52

Fig. 34

The dosed volume will be briefly displayed (Fig. 35) before the display returns to the main menu (Fig. 36).

**Dose**  
dosing NaOH

**2.000 ml**  
2.000 ml

ESC

20 ml      07/20/14 19:58

Fig. 35

**Main menu**

**0.000 ml**  
2.000 ml

dosing NaOH  
Method parameter  
Select method / system  
20 ml

START  
EDIT  
MODE

07/20/14 19:51

Fig. 36


The next dosage operation can be started immediately.

**i** The next dosage operation can be started immediately. Filling of the unit following dosage will not occur automatically, unless the maximum cylinder volume has been reached or the automatic filling option ("always") was activated:

The unit can be filled at any time using **<FILL>**.  
**<ESC>** will take you back to the main menu.


### 3.6.2.2 Dosing operation without dosing method

A dosing operation can also be performed without any dosing method with **<DOS>/<F10>** of the external keyboard (Fig. 37).

**Fig. 37**

The volume will be dosed with **<ENTER>/<OK>** (Fig. 38).

**Fig. 38**

The next volume can be carried out immediately with **<ENTER>/<OK>**.

Filling of the unit following dosage will not occur automatically here, unless the maximum cylinder volume has been reached.

The unit can be filled at any time using **<FILL>**.  
**<ESC>** will take you back to the main menu.

## 4 Method parameters

From the main menu, <EDIT> will take you to the method parameters (Fig. 39).


Fig. 39

### 4.1 Method editing and new method

If you select «**edit method**» or «**new method**» you will be taken to the modification or new creation of a method.

Selecting «**new method**» will always lead to the prompt for the input of a method name (Fig. 40). This prompt will not appear in the case of the modification of an already created method.


Fig. 40

The method name can contain up to 21 characters. Special characters are also possible.


**i** If no keyboard is connected, the method name being displayed **has to be** adopted.

Numbering of methods will occur automatically. Press <ENTER>/<OK> to confirm the input. The method name can be changed at any time.

Please continue at this point with  4.5 Change Method Parameters.


## 4.2 Default method

The «**Default methods**» item of the device contains a series of ready-made standard methods which can be conveniently selected (Fig. 41).

**Fig. 41**

Once the selection made, you are directly prompted for the input of the method name (Fig. 42).


**Fig. 42**

The standard name may be adopted or modified. Subsequently, you will be taken to «**Change method parameters**».

Please continue at this point with  4.5 Change Method Parameters.

## 4.3 Copy Method

Methods can be copied or stored with a new name (Fig. 43). If you select this function, the current method will be copied and you can include a new name.

**Fig. 43**

 A new name with the suffix [1] is assigned automatically in order to avoid the existence of two methods having the same name. Subsequently, you will be taken to <**Change method parameters**>. Please continue at this point with  4.5 Change Method Parameters.

## 4.4 Delete Method

In this function you will be prompted to know whether the current method is actually to be deleted (Fig. 44). You have to reply «Yes» in explicit terms and also confirm this reply with <ENTER>/<OK>.


Fig. 44

## 4.5 Change Method Parameters

The input or modification of the method name was already described in [Fig. 4.1](#) and [4.3](#)


Fig. 45

### 4.5.1 Method type

On the «Method type» you can select whether you wish to perform a manual titration or a dosage or whether you wish to prepare a solution (Fig. 46).


Fig. 46

**i** The selection of the Method type will have an influence on the further parameterization of the method. For instance, if you select the dosing mode, neither a selection of a formula will be available.

#### 4.5.2 Result

The «**Result**» menu offers the following possible settings (Fig. 47).


Fig. 47

The «**Result text**» may contain up to 21 alphanumeric characters including special characters (Fig. 48).


Fig. 48

Please confirm your input with <ENTER>/<OK>.

##### 4.5.2.1 Calculation Formula

The appropriate calculation formula is selected on the «**Formula selection**» submenu (Fig. 49).


Fig. 49

The following calculation formulae are available for manual titration:

Titration formula	Additional information
$(ml - B) * T * M * F1 / (W * F2)$	Formula for calculating the concentration of a sample taking into account a blank value in terms of ml.
$(B - ml) * T * M * F1 / (W * F2)$	Formula for calculating the concentration of a sample taking into account a blank value in terms of ml. Reverse titration (examples. CSB, saponification number)
$(B * F3 - ml * F1) * T * M / (W * F2)$	Formula for calculating the concentration of a sample taking into account a blank value, including a multiplicative factor. Reverse titration.
$(W * F2) / (ml - B) * M * F1$	Formula for calculating a titer (T) of a titration solution.
ml	Used to calculate the consumption in ml

The abbreviations used here have the following meaning:

- ml: Titration solution consumption in ml
- B: Blank value in ml. Mostly determined by way of titration
- T: Titer of the titration solution (e.g. 0.09986)
- M: Mol; mol- or equivalence weight of the sample (e.g. NaCl 58.44)
- F1: Factor 1, Conversion factor
- F2: Factor 2, Conversion factor
- W: "Weight", weighed-in quantity in g or volume in ml

**i** After selecting a formula, please confirm your selection with <ENTER>/<OK>.

The values of the individual parameters of the selected calculation formula can now be input one by one (Fig. 50 and Fig. 51).


Fig. 50


Fig. 51

#### 4.5.2.2 Sample weight and volume (sample quantity)

The Sample Quantity (W) item (Fig. 52) is used to select whether one is wishing to use a sample weight or a sample volume for titration or solution preparation (Fig. 53).

**Formula parameter**  
(EQ1-B)\*T\*M\*F1/(W\*F2)

T (Titre)	1.00000000 ▲
M (Mol)	36.46000
F1 (Factor 1)	1.0000
<b>W (Amount)</b>	<b>1.0000g ▼</b>

Selection   
 Enter   
 Back

20 ml NaOH 0.1 mol/L 09/13/11 12:45

Fig. 52

**Formula parameter**  
Amount

- Weight manual**
- Weight automatic
- Fixed weight
- Manuel Volume ▼

Selection   
 Enter   
 Back

20 ml NaOH 0.1 mol/L 09/13/11 12:44

Fig. 53

You have the following options:

- «**Manual sample weight**»: The sample weight is enquired by a prompt at the start of the method and manually input.
- «**Automatic sample weight**»: The sample weight is automatically transferred by a connected balance.
- «**Fixed sample weight**»: A fixed sample weight is input in g. This weight will then automatically be used for each start of the method.
- «**Manual sample volume**»: The sample volume in ml is prompted at the start of the method and manually input.
- «**Fixed sample volume**»: A fixed sample volume is input in ml. This volume will then automatically be used for each test of the method.

### 4.5.2.3 Formula unit

The formula unit can be selected in the «Unit» submenu (Fig. 54).


Fig. 54

Once the selection made (e.g. «%»), the unit will also be displayed as piece of information (Fig. 55).


Fig. 55

By pressing the «INS» (Insert) key on the external keyboard, you can also add new units.

### 4.5.2.4 Decimal digits

To conclude, it is possible to determine the number of decimal digits from 0 - 6. The standard setting is 2 (Fig. 56).



Fig. 56

#### 4.5.2.5 Global Memories


If a titration result is to be used again later, such as the factor or titer of a solution or a blind value, this can be saved automatically. The creation of a global memory is only possible if an external keypad is used. The creation of a global memory is possible in the «**system settings**». By pressing <SHIFT> + <F5> on the external keypad, you can access the «**Global Memories**» (Fig. 57).

**Fig. 57**


Using <F3> it is possible to add a global memory (Fig. 58).

**i** The name of the memory can be changed in reference to the application. \_\_\_\_\_


**Fig. 58**

The titrator proposes a memory name, such as **M01** (M01- M10) (Fig. 59).


**Fig. 59**

M01 can be accepted or provided with a designation such as blank value or titer (Fig. 60)


**Fig. 60**

This simplifies later the allocation of the global memory in another method (Fig. 61).


**Fig. 61**

The blind value which was possibly titrated in advance, is always taken into consideration automatically (Fig. 62).

**Fig. 62**

**Example:** The blank value of a chloride titration is defined with the support of an extra method. The result in ml is thereby automatically written into global memory M01 by using the name "Blanc value" (Fig. 63). The blank value is then automatically deducted from the titrant consumption within the chloride method.

**Fig. 63**

### 4.5.3 Dosing parameter

The dosing parameters (dosing speed, filling speed and max. dosing/titration volume) are determined for each method (Fig. 64). This applies to both types of methods such as manual titration (Fig. 65) and dosage (Fig. 66).

```

Edit method parameter
HCl titration
Method name
Method type          man
Result
Dosing parameter
Sample ID            man
Documentation        GLP
  ^ v |      OK |      ESC |
20 ml                07/21/14 12:40
  
```

Fig. 64

```

Edit dosing parameter
HCl titration
Dosing speed        100.00 %
Dosing speed        40.0 ml/min
Filling speed        20 s
Max. titration volume 50.000 ml
  ^ v |      OK |      ESC |
20 ml                07/20/14 22:50
  
```

Fig. 65

```

Edit dosing parameter
Dosieren
Dosing volume       0.000 ml
Dosing speed        100.00 %
Dosing speed        40.0 ml/min
Filling speed        30 s
Automatic filling    always
  ^ v |      OK |      ESC |
20 ml 0.1 NaOH      07/30/19 16:28
  
```

Fig. 66

The dosing speed can be set in % from 1 to 100 %. 100 % is the maximum dosing speed:


Dosing unit	Max. dosing speed [ml/min]
20 ml	40
50 ml	100

The filling speed can be set in terms of seconds from 20 to 240. The standard setting of this value is 30 seconds. For diluted aqueous solutions the filling speed can be six to 20 seconds. For non-aqueous solutions the filling speed should be set to the 30 seconds. In the case of highly viscous solutions such as concentrated sulphuric acid the filling speed should be further reduced down to 40 - 60 seconds.

The (maximum) the living volume or titration volume can be set to 9999.999.

#### 4.5.4 Sample identification


In the manual titration it is possible to input a sample identification (Fig. 67). The possible input includes **manual**, **automatic** or **no** sample description at all.

**Fig. 67**

For a sample description of the «**manual**», a prompt for the sample description will always be displayed at the start of the method (See also  3.6 Main Menu).

For an «**automatic**» sample description there will be selected a master description (e.g. Fig. 68 in the current case this is water), which will then automatically be numbered starting on 01.

**Fig. 68**

After a new power-up, numbering will resume with 01.

#### 4.5.5 Documentation

Three different format settings are available for documentation (Fig. 69) on a printer or USB device: «**short**», and «**GLP**» (Fig. 70).

```

Edit method parameter
HCl titration
Method name
Method type                man
Result
Dosing parameter
Sample ID                  auto
Documentation            GLP
  ^ v |      OK |      ESC |
20 ml                      07/21/14 12:45

```

Fig. 69

```

Documentation
HCl titration
Short
GLP
Only Display
  ^ v |      OK |      ESC |
20 ml                      07/21/14 12:46

```

Fig. 70

Method type	Short documentation	Standard documentation	GLP-documentation
Manual titration	Method name, date, time, sample description, sample weight/sample volume, results and calculation formula	N/A	Same as 'Short documentation' + plus method contents
Dosing	Only method printout, Method name, date, time, dosing parameters	N/A	N/A

## 5 System settings


Fig. 71

From the main menu (Fig. 71) you can access the system settings with <SYS> (Fig. 72).


Fig. 72

### 5.1 Setting the language

Setting the national language was already described in [2.4](#).

### 5.2 Dosing unit - Reagents

You can set up the attachment size in the menu (20 or 50 ml), perform an attachment change and enter reagent data, which are put out into the GLP documentation during manual titration (Fig. 73).


Fig. 73

### 5.2.1 Replacing the dosing unit

As a rule, the need for replacing the dosing unit occurs only rarely. The dosing unit has to be replaced, if such a replacement becomes necessary as a result of a defect, or of an inspection of the titration unit.

The dosing unit is equipped with lateral ribs around its circumference, with one of these ribs being in double design. This double rib serves as a mark for the correct placement of the dosing unit (Fig. 81).

With <ENTER>/<OK> confirm the «Unit size» and select «Dosing unit exchange» (Fig. 74).

**⚠ The exchange procedure starts directly with any additional warning!**  
Please take care that the titration tip is placed in a beaker or in the reagent bottle.


Fig. 74

The piston is raised to about 85% (Fig. 75).


Fig. 75

Then you will be prompted to unlock the attachment (Fig. 76).



Fig. 76

Now unlock the dosing unit (Fig. 77).

**Fig. 77**

Confirm with **<ENTER>/<OK>** after the attachment has been unlocked. Now, the attachment will be raised all the way up (Fig. 78).

**System settings**  
Dosing unit exchange

**Dosing unit moves up**

20 ml ESC 10/02/14 13:02

**Fig. 78**

You can now change the attachment (Fig. 79).


**System settings**  
Dosing unit exchange

**Please remove dosing unit  
and attach a new one**

20 ml OK 10/02/14 13:02

**Fig. 79**


Pull the attachment off toward the top (Fig. 80).

**Fig. 80**

Attach the new dosing attachment in the same manner (Fig. 81).


**i** The two struts of the UV protection must match up with the marking on the housing!

**Fig. 81**


Confirm with **<ENTER>/<OK>**

If you changed the attachment size, you can select the size here now (Fig. 82).


**Fig. 82**

If you want to change reagents, you can reset the data completely (Fig. 83).


**Fig. 83**

The attachment then moves down again (Fig. 84) and you will be asked to lock the attachment (Fig. 85).

**Fig. 84**

**Fig. 85**

### 5.2.2 Reagent data

The following reagents can be entered:

- Unit size 20 or 50 ml (selectable)
- Reagent name (default: empty)
- Concentration (default: 1.000000)
- Concentration determined on (default: empty)
- Expire date (default: empty)
- Opened/Produced on: (default: empty)
- Test according to ISO 8655: (default: empty)
- Batch ID: (default: empty)
- Last modification (default: current Datum).

### 5.3 Globale Memory

The handling with the global memories is already described in [4.5.2.5 Global Memories](#).

### 5.4 RS-232 Settings

The «RS232 settings» item can be used to determine the device address of the TITRONIC® 300 and set the parameters of the two RS-232 interfaces independent from each other (Fig. 86).


Fig. 86

The device address can be set from 0 - 15. Address 1 is the default setting (Fig. 87).


Fig. 87

The baud rate is preset to 4800 (Fig. 88).




Fig. 88

It may be set to 1200 - 19200 (Fig. 89).


**Fig. 89**

The parity can be selected amongst «No», «Even» and «Odd». «No» is the default setting (Fig. 90).

**Fig. 90**

You may select between 7 and 8 data bits. 8 bits is the default setting (Fig. 91).

**Fig. 91**

The stop bits can be set to 1, 1.5 and 2. The default setting is 1 (Fig. 92).


Fig. 92

**i** The RS-232 parameters can be set to the factory settings.

### 5.5 Date and Time

The factory time setting is Central European Time. This setting may be changed, where necessary (Fig. 93).


Fig. 93

### 5.6 RESET

RESET will reset all settings to the factory setting.

**i** All methods will also be deleted! So please print the methods or export/copy them to a connected USB storage medium (this will be possible with a higher update!).

The RESET has to be confirmed separately once again (Fig. 94).


Fig. 94

## 5.7 Printer

For connecting printers (Fig. 95) please refer to  7.3 Printers.


Fig. 95

## 5.8 Device Information

This point contains information about the device (Fig. 96).


Fig. 96

### 5.9 System Tone

The system tone can be switched on and off (Fig. 97).


Fig. 97

### 5.10 Data exchange

All methods with all parameter settings and global memories can be stored and restored on a connected USB-memory. It is also possible to transfer the settings from one titrator to another one. The backup will be started with «Settings backup» (Fig. 98).


Fig. 98

“Backup settings” is displayed during the backup in blue (Fig. 99).


Fig. 99

After a Reset or a maintenance case it is possible to restore the backup with «Restores settings» (Fig. 100)


Fig. 100

The backup folder on the USB-memory Stick starts with the backup date (Fig. 101).


Fig. 101

Confirm the selection with <ENTER>/<OK>.

“Settings are being restored” is displayed during the restoring process of the backup in blue (Fig. 102).




Fig. 102

### 5.11 Software Update

**Fig. 103**


An update of the device software (Fig. 103) requires a USB stick containing a new version. For this operation, the two files that are needed have to be located in the root directory of the USB device (Fig. 104).

**Fig. 104**


Plug the USB device into a free USB-A port, wait for some seconds, and then select the Software Update function. The valid software updates will be shown on the display.

In the present case (Fig. 105) this is Version "15\_50" from week 50 and year 2015.


**Fig. 105**

After starting the update using <ENTER>/<OK>, next thing to appear is the following graphic (Fig. 106),

**Fig. 106**

which will change after a few seconds to the following display (Fig. 107).


**Fig. 107**

Upon completion of the update (approx. 2 - 3 minutes), the device will shut down the software completely and proceed to a new start.

**i** In the course of an update, the methods will not be deleted! You can continue to use them.

If no valid update file is stored on the USB stick, a message will appear (Fig. 108):

**Fig. 108**

## 6 Communication via RS-232 and USB-B interface

### 6.1 General Information

The TITRONIC® 300 has two serial RS-232-C interfaces to communicate data with other devices. By means of these two interfaces it is possible to operate several devices on one computer (PC) interface. In addition to that, the TITRONIC® 300 also has an alternatively USB-B interface, which can only be used to connect a PC. RS-232-C-1 establishes the connection to a connected computer or to the previous device of the “Daisy Chain”. At the RS-232-C-2 it is possible to connect additional devices (Daisy Chain Concept).

PIN assignment of the RS-232-C interfaces:

PIN No.	Meaning / Description
1	T x D Data output
2	R x D Data input
3	Digital mass

### 6.2 Chaining multiple devices - “Daisy Chain Concept”

In order to activate several devices in a chain individually, each device must have an own device address. For this it is at first necessary to establish a connection from the computer to the RS-232-C interface 1 of the first device in the chain by means of a RS-232-C data cable, e.g. Type No. TZ 3097. With the additional RS-232-C data cable, Type No. TZ 3094, the RS-232-C- interface 2 of the first device is connected with the RS-232-C-interface 1 of the second device. At interface 2 of the second device it is possible to connect an additional device.

The TITRONIC® 300 can also be connected via USB cable TZ 3840 (type A (M) - type B (M), 1.8m) to a USB interface of a PC. To accomplish this connection, a driver has to be installed on the PC. Then the USB-B interface takes over the function of the RS-232-1 interface.

The address always consists of two characters: e.g. address 1 of the two ASCII- characters <0> and <1>. The addresses can be set from **00** to **15**, i.e. 16 possibilities. It must be ensured that the devices in a chain have different addresses. If a device is addressed with its address, this device will process this command without sending it to another device. The reply to the computer has also an own address. The addresses are allocated as described in  5.4 RS-232 Settings.

The TITRONIC® 300 receives commands from a PC at the interface **1** (USB- B) if the computer knows the address. It also sends the answer via this interface. If the address of the incoming command does not match the device address, the complete command will be forwarded to interface **2**. Interface 2 is connected to interface 1 of another device. This device checks the address as well and reacts to the command as the first TITRONIC® 300 did before.

All information (data strings) which arrive at interface 2 of the TITRONIC® 300 will immediately be send to the computer via interface 1 (or USB-B interface). Thus, the computer receives the data of all devices. In practice it is possible to connect up to 16 devices to one computer- (PC-) interface.

### 6.3 Instruction Set for RS-Communication

The commands consist of three parts:

Adresse, two-digit aa	e.g. <b>01</b>
Command	e.g. <b>DA</b>
Variable, if necessary	e.g. <b>14</b>
an end of command	<b>&lt;CR&gt; &lt;LF&gt;</b>

 **Every command must be completed with the ASCII - sign <CR> and <LF>** (Carriage Return and Line Feed). Only if the respective action has ended the answers will be returned to the computer.

Example:

The command to dose 12.5 ml shall be sent to the TITRONIC® 300 with the address 2.

The command consists of the characters:

**02DA12.5<CR LF>** in detail:

02	=	Device address
DA	=	Dosage command with filling and zero points of the display
12.5	=	Volume in ml to be dosed
<CR LF>	=	Control character as command end

Command	Description	Reply
aaAA	automatic allocation of device address	aaY
aaMC1...XX	choosing a method	aaY
aaBF	„filling burette“. Aufsatz wird gefüllt.	aaY
aaBV	output of dosed volume in ml	aa0.200
aaDA	dose volume without filling, with adding the volume	aaY
aaDB	dose volume without filling, reset of the volume	aaY
aaDO	dose volume with filling, without adding the volume	aaY
aaGDM	dosing speed in ml/min	aaY
aaGF	filling time in seconds (min is 20, default 30)	aaY
aaES	„ESC“ function one step backwards	aaY
aaEX	„EXIT“ function.back to main menu	aaY
aaGDM	dosing speed in ml/min (0.01 – 100 ml/min)	aaY
aaGF	filling time in sec (adjustable 20 – 999 seconds)	aaY
aaGS	output serial no. Of device	aaGS08154711
aaLR	output report (short report)	aaY
aaLI	output method content	
aaLO	output documentation (as configured)	
aaRH	request of identification	aalident: TITRONIC® 300
aaRC	send last command	aa“last command“
aaRS	report status	aaStatus:“text
	possible answers are:	
	„STATUS:READY“ for ready	
	„STATUS:dosing“ dosing	
	„STATUS:filling“ filling	
	„ERROR:busy“ if no interchangeable unit has been attached	
aaSM	start selected method	aaY
aaSEEPROM	EEPROM reset to factory defaults	aaY
aaSR	stop the actual function	aaY
aaSYS5	adjust language to „German“	aaY
aaSYS1	adjust language to „English“	aaY
aaSYS2	adjust language to „French“	aaY
aaSYS3	adjust language to „Spanish“	aaY
aaVE	Version number of the software	aaVersion:

## 7 Connection of Analytical Balances and Printers

### 7.1 Connection of Analytical Balances

As it often happens that the sample is weighed in on an analytical balance, it makes sense to connect this balance to the TITRONIC® 300. The balance must have a RS-232-C-interface and the connection cable must be configured accordingly. For the following types of balances there are already assembled connection cables:

Balance	TZ-Number
Sartorius (all type with 25-pole RS-232), partially Kern	TZ 3092
Mettler, AB-S, AG, PG, Sartorius with USB-Port	TZ 3099
Precisa XT-Serie	TZ 3183
Kern with 9-pole RS-232	TZ 3180

For all other types of balances it is possible to obtain an already assembled connection cable (on demand). For this we need detailed information about the RS-232-C-interface of the balance used.

The connection cable is to be connected to the RS-232-C-interface 2 of the TITRONIC® 300. This side of the connection cables always consists of a 4-pole mini-plug. The other side of the cable can, depending on the type of balance, be a 25-pole plug (Sartorius), a 9-pole plug (Mettler AB-S) or a 15-pole specialised plug (Mettler AT) etc.

In order to allow the balance data to be sent to the TITRONIC® 300, the data transmission parameters of the titrator and the balance must correspond to each other. Additionally, it is necessary to carry out some more standard settings on the side of the balances:

- The balance is to send the balance data via RS-232-C only by means of a print command
- The balance is to send the balance data only after the display standstill
- The balance should never be set to “automatic sending” and/or “send continuously”
- “Handshake” on the balance must be set to “off”, or even “Software Handshake” or “Pause”

No special characters such as **S** or **St** are allowed to be used as prefix in the balance data of the balance data string. In such a case it might be possible that the TITRONIC® 300 cannot process the balance data correctly.

After you have connected the balance with the appropriate cable and have adjusted all settings in the balance software, and possibly in the TITRONIC® 300, you can now test the data transfer of the balance very easily. Start the one method. Confirm the sample designation. Then, the display asks you:

- a) to press the print-button at the balance  
→ Parameters to “weighted sample automatically”
- b) to enter the weighted sample → then the parameters are still set to “weighted sample manually”

Put an object onto the balance and press the print button.

After the standstill of the balance display there will be beep and the transmitted balance data appear:

- a) the display changes automatically into the measuring display.
- b) the weighted sample must again be confirmed with <ENTER>/<OK>.

If no balance data is available, the «No balance data found» message will appear (Fig. 112).


Fig. 112

### 7.3 Printers

The results, calibration data and methods can be printed on the following media

- HP PCL compatible printer (A4)
- Seiko DPU S445 (Thermo paper 112 mm width)
- On the USB stick in PDF- and CSV-format

To connect the printers to the burette please use the USB socket.

When printing, please check whether the correct printer is connected. It is not possible to print "HP" printer layouts on another thermal printer or vice versa. The printer settings should always be checked and adjusted after changing the printer (Fig. 113).


Fig. 113

**i** Only one printer should be connected for one titrator because automatic printer recognition is not activated. «Print PDF» is the default setting.

## 8 Maintenance and Care of the piston burette

 The preservation of the proper functioning of the device requires testing and maintenance work to be performed on a regular basis. Regular inspections are essential prerequisites for the correctness of the volume and the proper functioning.

The accuracy of the volume is determined by all chemicals-carrying components (piston, cylinder, valve, titration tip and hoses). These parts are subject to wear and tear. The piston and cylinder are subject to particular strain, hence they require special attention.

### Heavy strain:

Use of e.g. concentrated solutions, reagents and chemicals (> 0.5 mol/L); chemicals attacking glass, such as fluorides, phosphates, alkali solutions; solutions with a tendency to crystallising out; Fe (III) chloride solutions; oxidising and corroding solutions such as iodine, potassium permanganate, Cer (III), Karl-Fischer titration agent, HCl; solutions with a viscosity of > 5 mm<sup>2</sup>/s; frequent, or even daily use.

### Normal strain:

Use of solutions, reagents and chemicals (up to 0.5 mol/l) which do not attack glass, crystalize out or corrode.

### Interrupted use:

If the dosing system is not in use for more than two weeks, we recommend emptying and cleaning the dosing unit [6]. This applies in particular under the operating conditions referred to in the „Heavy strain” section. If this recommendation is not adhered to, the piston of the valve may become leaking, this may result in damage to the piston burette.

 If the liquid is left within the system, you will also have to reckon with corrosion and an alteration of the solutions used over time, which includes e.g. crystallisation. Considering that as of the state of the art there are no plastic hoses available for the use in titration equipment which would be perfectly free of diffusion phenomena, particular attention is to be paid to the range of the hose lines.

### We recommend the following inspection and maintenance work:

	Heavy strain	Normal strain
Simple cleaning: • Wiping off splashed chemicals from the outer surface [1]	Whenever required in operation	Whenever required in operation
Sight check: • Check for leakage in the area of the dosing system. [2] • Is the piston tight? [3] • Is the valve tight? [4] • Titration to clear? [5]	Weekly, when putting back into operation	Monthly, when putting back into operation
Basic cleaning of the dosing system: • All parts of the dosing system to be cleaned separately. [6]	Every three months	Whenever necessary
Technical inspection: • Check for air bubbles in the dosing system. [7] • Visual inspection • Check of the electrical connections. [8]	Semi-annually, when putting back into operation	Semi-annually, when putting back into operation
Verification of the volume according to ISO 8655 • Perform basic cleaning • Inspection according to ISO 8655 Part 6 or Part 7. [9]	Semi-annually	Annually

 Depending on the respective application, there may be different specifications for the entirety of the inspection and maintenance work to be performed. The individual intervals may be extended if no complaints occur, but they will have to be shortened again as soon as any problem has arisen

The inspection of the metrological reliability including maintenance work is offered as a service (including a manufacturer's certificate, if so ordered). In this case the titration device is to be sent in. Please contact the service (see backside of this manual).

## Detailed description of the inspection and maintenance work

- [1] Wipe off using a soft cloth (and some water with a normal household detergent).
- [2] Leaking connections can be identified by moisture or crystals at the threaded connections of the hoses, at the sealing lips of the piston inside the dosing cylinder or at the valve.
- [3] If any liquid becomes visible below the first sealing lip, it has to be checked at short timely intervals whether any liquid will build up under the second sealing lip, too. In this case both the piston and the glass cylinder have to be replaced immediately. It is easily possible that in operation small liquid droplets build up under the first sealing lip, but they may also disappear again. This phenomenon alone is no reason for replacement.
- [4] The valve has to be removed from its housing for inspection. In this process, the hoses remain connected to the valve. Please check for moisture underneath the valve. When reinserting the valve, please make sure that the small cam at the rotating axis is fitted into the corresponding groove again.
- [5] The titration tip must be free of sedimentation or crystals which might obstruct the dosing process or falsify the results.
- [6] Remove the cylinder, take the valve out of the valve housing, unscrew the hoses and then rinse all parts carefully with distilled water. For the assembly of the cylinder, hoses and other parts of the interchangeable unit, please refer to the operating instructions.
- [7] Dose one burette volume, then refill. Air bubbles will gather at the tip of the cylinder and in the titration hose where they can be detected easily. If bubbles become visible, please re-tighten all connections finger tight, and then repeat dosing. If air bubbles still remain within the system, [6] please check the valve and replace the hose connections. The air bubbles may also occur at the interface between the sealing lip of the piston and the cylinder. If a reduction of the filling speed will not do, the dosing unit has to be replaced.
- [8] Check the electrical plug contacts for corrosion and mechanical damage. Defective parts have to be repaired or replaced by new parts.
- [9] Please refer to the application "Burette inspection according to ISO 8655 Part 6".

## 9 Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts, as lobes, cylinders, valves and pipes including the thread connections and the titration tips. The breach of glass parts is also excluded. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

## 10 Storage and transportation

If the TITRONIC® 300 or the interchangeable units have to be stored over some time, or to be dislocated, the use of the original packing will be the best protection of the devices. However, in many cases this packing will not be available anymore, so that one will have to compose an equivalent packaging system. Sealing the lower section in a foil is hereby recommended. The devices should be stored in a room with a temperature between +10 and +40°C, and the (relative) humidity of the air should not exceed 70 %.

 If the interchangeable have to be stored over some time, or to be dislocated, the fluids inside the system, especially aggressive solution have to be removed.

## 11 Recycling and Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

The TITRONIC® 300 and her packaging are manufactured as far as possible from materials which can be disposed of environmental-friendly and recycled in a technically appropriate manner. If you have any question regarding disposal, please contact the service (see backside of this manual).

 The main printed board carries a lithium battery (type CR 2430). Batteries should not to be disposed of with the normal domestic waste. They will be taken back and recycled or disposed of properly by the manufacturer at no cost.

## 12 EC – Declaration of Conformity

The corresponding declaration of conformity of the device can be found on our homepage. It will also be made available to you on request.

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
**Xylem Analytics Germany Sales GmbH & Co.KG**
© 2021 Xylem, Inc. Version 211029 D 825 288 0