**GEBRAUCHSANLEITUNG**  
Originalversion  
**OPERATING MANUAL**  
**MODE D'EMPLOI**  
**MANUAL DE INSTRUCCIONES**

# ViscoPump III

VISKOSITÄTSMESSGERÄT  
VISCOSITY MEASURING DEVICE  
L'APPAREIL DE MESURE DE LA VISCOSITE  
EQUIPO MEDIDOR DE VISCOSIDAD

**SI Analytics**

a **xylem** brand

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## **Operating Manual ..... Page 25 ... 46**

### **Important notes:**

The operating manual is part of the product. Before initial operation, please carefully read and observe the operating manual and keep it. For safety reasons the product may only be used for the purposes described in these present operating manual. Please also consider the operating manuals for the devices to be connected.

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties. A potentially more recent version of this manual is available on our internet website. The German version is the original version and binding in all specifications!

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## TABLE OF CONTENTS

<b>1</b>	<b>Specifications of the of the ViscoPump III.....</b>	<b>27</b>
1.1	Notes to the operating manual .....	27
1.2	Intended Use .....	27
1.3	Technical Specifications.....	27
1.3.1	ViscoPump III .....	27
1.4	Warnings and safety information.....	29
1.4.1	General.....	29
1.4.2	Chemical and biological safety.....	30
1.4.3	Flammable liquids .....	30
1.5	Functioning of the device .....	30
1.6	Capillary viscometry .....	30
1.7	Measurement principles .....	31
<b>2</b>	<b>Initial operation .....</b>	<b>32</b>
2.1	Unpacking .....	32
2.2	Installation and dismantling of the module .....	32
2.2.1	Installation .....	32
2.2.2	Dismantling.....	33
2.3	Connecting the devices .....	34
2.3.1	Suitable viscometer types, racks, and measurement stands.....	34
2.3.2	Deployable hose combinations .....	36
2.4	Connecting the viscometers and other devices .....	37
2.4.1	Viscometers using light-barrier sensing .....	37
2.4.2	TC viscometers with thermistor sensors .....	38
2.4.3	Connection of VZ 7215 absorption traps.....	38
2.4.4	Connection of the VZ 8552 overflow guard.....	38
2.5	Troubleshooting.....	39
<b>3</b>	<b>Data transfer.....</b>	<b>40</b>
3.1	RS-232-C Interfaces.....	40
3.2	Interface configuration.....	40
<b>4</b>	<b>Software-Update .....</b>	<b>41</b>
<b>5</b>	<b>Maintenance and care of the device and the viscometers .....</b>	<b>42</b>
5.1	Maintenance work to be carried out.....	42
5.2	Maintenance and care of the VZ 7215 absorbent bottle .....	42
5.3	Periods without operation.....	43
5.4	Reproducibility of results .....	43
5.5	Viscometers within quality assurance systems.....	43
<b>6</b>	<b>Guarantee .....</b>	<b>44</b>
<b>7</b>	<b>Storage and transportation.....</b>	<b>44</b>
<b>8</b>	<b>Recycling and Disposal.....</b>	<b>45</b>
<b>9</b>	<b>EC – Declaration of Conformity .....</b>	<b>45</b>

---

**Copyright**

© 2022, Xylem Analytics Germany GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization.  
Germany. Printed in Germany.

# 1 Specifications of the of the ViscoPump III

## 1.1 Notes to the operating manual

The provided operating manual will allow you the proper and safe handling of the product. For maximum security, observe the safety and warning instructions in the operating manual!

-  Warning of a general danger:  
Non-compliance results (can result) in injury or material damage.
-  Important information for device use.
-  Refers to another part of the operating manual.

The menu screens shown in this operating manual serve as an example and may differ from what you see!

## 1.2 Intended Use

The ViscoPump III module is a functional unit used to determine absolute and relative viscosity of liquid media. The module is controlled exclusively through the use of a special program to be run either on a computer or on a suitable control unit (such as the AVS<sup>®</sup> 470 Viscosity Measuring Unit). The module forms an integral part a viscosity measuring device and is designed for being built into the respective device or controller unit.

## 1.3 Technical Specifications

### 1.3.1 ViscoPump III

Translation of the legally binding German version

(Release: 18. Oct 2021)

 **The data given below is only valid if the ViscoPump III module is built into a certified device (e.g. AVS<sup>®</sup> 370, AVS<sup>®</sup> 470, AVS<sup>®</sup> Pro)!**

EMC compatibility according to the Council Directive: 2014/30/EU;  
applied harmonized standards: EN 61326-1

Low-voltage directive according to the Council Directive 2014/35/EU;  
Testing basis EN 61 010-1: for laboratory equipment  
RoHS Council Directive 2011/65/EU  
FCC Part 15B and ICES 003

**Country of origin:** Germany, Made in Germany

**Measurement parameters:** Flow time in seconds [s]

**Capture of measurement value:**

Flow time: Optoelectronic or thermo-resistive capture of the meniscus passage through the measurement planes of the viscometers

**Measurement ranges:**

Viscosity: “pressing” mode 0.35 ... 1800 mm<sup>2</sup>/s (cSt) measuring temperature of approx. 20 ... 25 °C  
“sucking” mode 0.35 ... 5800 mm<sup>2</sup>/s (cSt) measuring temperature of approx. 20 ... 25 °C

Time: 5 s to 9999.99 s, Resolution 0.01 s

Pumping pressure: Fully automatically controlled “sucking” action to approx. -160 mbar  
Fully automatically controlled “pressing” action to approx. +160 mbar

**Measurement precision:**

Precision (reproducibility and comparability) in accordance with DIN 51562, Part 1

Time measurement:  $\pm 0.01 \text{ s} \pm 1 \text{ digit}$ , but not more accurate than 0.01 %  
measurement uncertainty in the determination of absolute, kinematic viscosity furthermore depends on the uncertainty of the numerical value of the viscometer constants and the measurement conditions, especially as concerns the measurement temperature.

**Power supply:** via AVS<sup>®</sup> measuring systems AVS<sup>®</sup> 370, AVS<sup>®</sup> 470 und AVS<sup>®</sup> Pro  
 Input: +24 Volt DC, 800 mA  
 Power Input: 19 W

**Data transfer parameters:**

Data interface: bidirectional, serial interface according to EIA RS-232-C  
 Data format: word length 7 bits, 2 stop bits, 4800 baud, no parity

**Connections Back panel of the device:**

96-channel VG bar for contact making with the module plugged in a controller unit, including voltage supply and 2 serial interfaces

Gas exhaust pipe: Stainless-steel tube used to expel or suck the gases delivered by the rotary pump

**Front panel (to be connected to the front panel of ViscoPump III modules)**

Pneumatic connectors: Pressure/suction venting, to be connected to viscometer.

Overflow guard for VZ 8552 suction line:  
 Capacitive sensor: Round-plug connector, 4-channel,  
 Screwed cap according to DIN 45321

**Type VZ 8561 ViscoPump III for optoelectronic sampling**

Measurements stand AVS<sup>®</sup>/S:  
 Round-plug connector with bayonet catch,  
 DIN 5 channels for AVS<sup>®</sup>/S (measurement stand), 5-channel DIN socket

**Type VZ 8562 ViscoPump III for TC sampling**

TC viscometer: 4-channel DIN socket

**Operating indicator:** Green LED (voltage indicator)

**Various add-on devices:**

EIA-RS-232-C interface with 9-channel Submin D plug

**Casing:** No casing, module in plug-in design

Dimensions: approx. 50 x 173 x 140 mm (W x H x D) overall dimensions

Weight: approx. 0.8 kg

Front panel: 2 mm anodized aluminum

**Ambient conditions:**

 **Do not use the device in hazardous locations!**

Climate: Ambient temperature: + 10 ... + 40 °C for operation and storage  
 Humidity according to EN 61 010, Part 1:  
 Max. relative humidity 80 % for temperatures up to 31 °C,  
 linear decrease down to 50 % relative humidity at a temperature of 40 °C

Pollution degree: Pollution degree IP 20, indoor use only.

## 1.4 Warnings and safety information

### 1.4.1 General

The device corresponds to protection class III.

It was manufactured and tested according to DIN EN 61 010, Part 1, “**Protective Measures for electronic measurement devices**” and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information contained in the present operating instructions. Development and production is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard.

 For reasons of safety, the device must only be used for the range of application described in the present operating manual. Nonobservance of the intended proper use of the device may result in personal injury or damage to property.

 **The ViscoPump III modules control the entire measurement process!**

(among other things the temperature pre-adaptation process of the samples in the viscometers, the process of pumping the liquid up into the storage containers of the viscometers, the measurement of the flow times etc.).

If the module is incorporated into an AVS<sup>®</sup> 370 or AVS<sup>®</sup> 470 Viscosity Measuring Unit, please verify after start-up that the proper operating mode is set, i.e. “sucking” or “pressing” action; these operating modes can be seen on the corresponding LED on the front panel of the AVS<sup>®</sup> 470, or in the context menu of the WinVisco software of the AVS<sup>®</sup> 370 Viscosity Measuring Unit, respectively. Please refer to the user manual of the device concerned to make the proper setting for the application, and make sure to use the suitable hose set.

 For reasons of safety, the device must be opened by authorized persons only! This means, for instance, that work on electrical equipment must only be performed by qualified specialists. In the case of unauthorized intervention in the device, as well as in the case of negligently or deliberately caused damage, the warranty will become void.

 **If it has to be assumed that safe operation is impossible, the ViscoPump III module has to be put out of the device and secured against inadvertent putting to operation!**

(In this case please switch the basic unit off, pull plug of the mains cable out of the mains socket, remove the device and return it to our service).

Examples for the assumption that a safe operation is no longer possible,

- the package is damaged
- the unit shows visible damages,
- the unit does not function properly
- liquid has penetrated into the casing. the package is damaged

In case that the user operates such a device, all thereof resulting risks are on the user!

 The device must not be stored or operated in humid rooms.

 **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. Special care is essential with acidic, flammable, toxic or otherwise health hazardous samples/chemicals, which may be measured during capillary viscometry or which are used as a rinsing solvent. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the environment or that they are supervised by specialized persons, respectively.

 When handling the substances used, the chemical resistance of the materials of the device must be taken into account.

 For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

 Please read the separate operating instructions of the transparent thermostats as well.

### 1.4.2 Chemical and biological safety

**i** The device is not intended for use with potentially biohazardous substances.

**⚠** **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the devices or that they are supervised by specialized persons, respectively.

**⚠** When using biohazardous substances, the regulations for handling the substances used must be observed. In such cases, the use is the sole responsibility of the user.

**⚠** For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

**⚠** Dispose of all used solutions in accordance with national regulations and laws. Select the type of protective equipment according to the concentration and quantity of the hazardous substance at the respective workplace.

### 1.4.3 Flammable liquids

When handling flammable liquids, make sure that there is no naked flame in the vicinity of the equipment. Adequate ventilation must be provided. Only small quantities of flammable liquids should be kept in the workplace.

## 1.5 Functioning of the device

The ViscoPump III Module is used to perform flow-time measurements in capillary viscometers.

The available capillary viscometers enable viscosity measurements 0.35 of approx. 5,000 mm<sup>2</sup>/s (cSt) to be carried out at a measuring temperature of approx. 20 ... 25 °C.

The ViscoPump III Module can be fitted out in two ways to be used for meniscus sensing. Connecting TC viscometers to the ViscoPump III VZ 8562 module will also enable the measurement of black and opaque liquids. As an alternative option it is possible to use viscometers in combination with the optoelectronic ViscoPump III VZ 8561 module viscometers for meniscus sensing with the measurement stand, e.g. the AVS<sup>®</sup>/S.

Time recording extends up to 9999.99 seconds with a resolution of 0.01 s. Serial data transmission is used to transfer the time readings to the respective controller unit.

Prior to the measurement as such, the liquid to be measured is sucked upwards inside the capillary viscometer through two measurement planes (N2 and N1) which are designed as light barriers or thermistor sensors, depending on the viscometer type (see Fig.1 and Fig. 2).

When using Ubbelohde viscometer, the design of the program ensures that the suspended spherical level will form prior to the start of the measurement.

## 1.6 Capillary viscometry

Capillary viscometry is the most accurate method for the determination of the viscosity of liquids with a Newtonian flowing behaviour. The measurement as such consists in a time measurement. The time measured is that which a specific quantity of liquid requires to pass through a capillary having a defined width and length. Conventionally, this process is watched with the human eye, and the flow time is measured manually using a stop watch.

In the case of the AVS<sup>®</sup> 370 Viscosity Measuring Unit, as with all viscometry measuring devices from SI Analytics<sup>®</sup>, the liquid meniscus is captured on the measurement planes, either in an optoelectronic manner by means of light barriers, or else on a thermo-resistivity basis by thermistors.

## 1.7 Measurement principles

### a) ViscoPump VZ 8561: Optoelectronic sensing of the liquid meniscus

The near-infrared light which is generated in LEDs located in the upper section of the measurement stand is conducted through a glass-fiber light-conductor cable onto the measurement planes. The light shines through the viscometer before it arrives at another light-conductor cable located on the opposite side; inside this second cable, the light is conducted to a receiver in the upper section of the measurement stand.

While the liquid meniscus passes through the measurement planes, the lens-like effect of the meniscus causes a short-term darkening of the light beam, followed by a magnification. This process generates a measurement signal which can be evaluated accurately.


Fig.1 Viscometer for optoelectronic measurements

### b) ViscoPump VZ 8562: Viscometer with thermistor sensors (TC viscometer)

In the case of TC viscometers, glass-coated thermistors serving as sensors are inserted on the level of the measurement planes. While the meniscus passes through the measurement planes, the differences in the thermal conductivity properties of air and liquid lead to a change in the heat balance. The thermistors of the TC viscometers are inserted hermetically tightly into the glass coating of the viscometer.


Fig. 2 Viscometer for optoelectronic measurements

## 2 Initial operation

### 2.1 Unpacking

The ViscoPump III module is delivered in a folding packaging with a plastic-foam protection. Please take the module out of this packaging and make sure that all items specified on the packing list are present.

### 2.2 Installation and dismantling of the module

**⚠ In the work mode "Press", in case of malfunction, the sample can emerge from the capillary tube and could enter the thermostat bath!** Therefore, especially during the pressing work mode, a correct connection of the measuring tripod to the ViscoPump III and a properly used viscometer are very important for the safe operation. With health hazardous and/or corrosive samples (e.g. acidic, toxic or hazardous to the environment), the work mode "suction" is the preferred choice.

**⚠ In the work mode "suction", the sample can be pumped into the ViscoPump in case of a malfunction!** In order to avoid this and the associated damage to the ViscoPump III, you must always use a suitable tube kit including the safety bottles VZ 7021 or VZ 7022 during the work mode "suction". In order to review a possible immersion of the sample into the safety bottle, a capacitive sensor VZ 8552 is highly recommended.

#### 2.2.1 Installation

**⚠ Please switch off the basic / controller unit first! Prior to installing or replacing this functional unit, please make sure that the mains plug is pulled off of the mains socket!**


The ViscoPump III module is installed by inserting it into the slot provided. In this process, please see that both the printed board and the metal shielding plate are inserted correctly into the corresponding guidance rails, and make sure that the socket bar locks in place. Subsequently, complete the installation by tightening the front-panel screws.

After switching on, the operating indicator (LED) flashes for approx. 4 seconds.

**⚠ The ViscoPump III modules control the entire measurement process, among other things the temperature pre-adaptation process of the samples in the viscometers, the process of pumping the liquid up into the storage containers of the viscometers, the measurement of the flow times etc.!**

If the module is incorporated into an AVS<sup>®</sup> 370 or AVS<sup>®</sup> 470 Viscosity Measuring Unit, please verify after start-up that the proper operating mode is set, i.e. "sucking" or "pressing" action; these operating modes can be seen on the corresponding LED on the front panel of the AVS<sup>®</sup> 470, or in the context menu of the WinVisco software of the AVS<sup>®</sup> 370 Viscosity Measuring Unit, respectively. Please refer to the operating instructions of the device concerned to make the proper setting for the application, and make sure to use the suitable hose set!

**⚠ The AVS<sup>®</sup> 370 can also be operated with ViscoPump II - modules, the older ViscoPump generation. A mixed configuration with ViscoPump II and III - modules is also possible. In the case of a mixed configuration, please note: **The ViscoPump II - modules may only be used in the first measuring positions, i.e. in the slots to the left of the ViscoPump III - modules** (see  Fig. 3). If the ViscoPump modules are placed in the wrong order, transmission errors will occur during RS communication!**

**Fig. 3 Measuring positions ViscoPump II and III – Module in the AVS<sup>®</sup> 370**

With all further connected devices with serial communication (e.g. burettes for rinsing the viscosity meters) the following must be adjusted as communication parameter:

Baud: 4800

Data bit: 7

Stop bits: 2

Parity: No

### 2.2.2 Dismantling

** Please switch off the basic/controller unit first! Prior to dismantling or replacing this functional unit, please make sure that the mains plug is pulled off of the mains socket!**


To dismantle or replace the ViscoPump III module, please proceed as follows:

- Remove the pneumatic and electrical connections from the front panel of the ViscoPump III module to be replaced.
- Loosen the screws located at the corners of the front panel.
- Use the upper and lower insertion handles to leverage the ViscoPump III module out of its rear-side plug connection.
- Pull the ViscoPump III module out of the AVS 470 Viscosity Measuring Unit.
- After inserting the new ViscoPump III module, please secure it again with the screws of the front panel. Re-establish the electrical and pneumatic connections.


## 2.3 Connecting the devices

### 2.3.1 Suitable viscometer types, racks, and measurement stands


Viscosimeter (Type)	Rack (Type no.)	Measuring stand (Type)
Ubbelohde (DIN) 532... 530... 501... 541... 545...	053 92	AVS <sup>®</sup> /S AVS <sup>®</sup> /SK
Ubbelohde (ASTM) 525... 526... 527... 545...	053 92	AVS <sup>®</sup> /S AVS <sup>®</sup> /SK
Mikro Ubbelohde 536... 537... 538...	053 92	AVS <sup>®</sup> /SAVS <sup>®</sup> /SK
Ubbelohde dilute-solution viscometers 531...	---	AVS <sup>®</sup> /SK-V
Cannon-Fenske-routine 513... 520...	---	AVS <sup>®</sup> /SK-CF
Mikro-Ostwald 516... 517... 518...	053 97	AVS <sup>®</sup> /S AVS <sup>®</sup> /SK
Ubbelohde (TC) 562... 563... 564... 567... 568... 569...	053 93	---

DIN-Ubbelohde-viscometer


Micro-Ubbelohde-viscometer

Cannon-Fenske-Routine-viscometer

Micro-Ostwald-viscometer

Ubbelohde-viscometer  
with TC sensorMicro-Ubbelohde-viscometer  
with TC sensor**Fig. 4 Suitable viscometer types**

### 2.3.2 Deployable hose combinations

Hose combinations Type no.	Description	Application
VZ 5505	Silicone hose kit, oppressive, for Ubbelohde viscometers (3 legs), and Cannon-Fenske and Ostwald viscometers	Standard, but please observe: the sample can leak from the capillary tube during a malfunction
VZ 5505 + VZ 8526	Silicone hose kit, suctioning, for Ubbelohde viscometers (3 legs)	Standard, safer than oppressive operation, as the sample cannot leak from the capillary tube. Unsuitable for volatile samples.
VZ 8523	PTFE hose set, suction, for Ubbelohde viscometer (3 legs)	For aggressive samples that attack silicone, such as sulphuric acid. Hose lengths are designed to fit the AVS <sup>®</sup> 370 on the support console VZ 8571. All hose sets for suctioning operation can be combined with the suction set and the sample fill set.
VZ 8524	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (3 legs)	For aggressive samples whose vapors are absorbed by soda lime filters to protect the ViscoPump. Depending on the sample, an active carbon filter VZ 7216 can be used instead of the soda lime filter VZ 7215. Hose lengths are designed to fit the AVS <sup>®</sup> 370 on the support console VZ 8571.
VZ 7218 + VZ 8535	PTFE hose kit, suctioning, for Ubbelohde viscometers (4 legs)	For Ubbelohde viscometers with additional 4 <sup>th</sup> pipe for filling and cleaning Hose lengths are designed to fit the AVS <sup>®</sup> 370 on the support console VZ 8571.
VZ 8530	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (4 legs)	For Ubbelohde viscometers with additional 4 <sup>th</sup> pipe for filling and cleaning Hose lengths are designed to fit the AVS <sup>®</sup> 370 on the support console VZ 8571. For aggressive samples whose vapors are absorbed by soda lime filters to protect the ViscoPump.
VZ 5606	For TC viscometers (3 and 4 legs) with screw connections: Silicone hose fittings with connecting cable. For oppressive operation	For TC viscometers (3 and 4 legs) with 4 <sup>th</sup> pipe for filling and cleaning Typical applications are measurements of oils.

**i** The hose combinations are to be selected on the basis of the required application.

**⚠** **When using Micro Ubbelohde viscometers with TV sensors, the ignition temperature of the media to be measured has to be taken into account!** It has to be higher than 250 °C.

**⚠** **In the work mode "Press", in case of malfunction, the sample can emerge from the capillary tube and could enter the thermostat bath!** Therefore, especially during the pressing work mode, a correct connection of the measuring tripod to the ViscoPump III and a properly used viscometer are very important for the safe operation. With health hazardous and/or corrosive samples (e.g. acidic, toxic or hazardous to the environment), the work mode "suction" is the preferred choice.

**⚠** **In the work mode "suction", the sample can be pumped into the ViscoPump in case of a malfunction!** In order to avoid this and the associated damage to the ViscoPump III, you must always use a suitable tube kit including the safety bottles VZ 7021 or VZ 7022 during the work mode "suction". In order to review a possible immersion of the sample into the safety bottle, a capacitive sensor VZ 8552 is highly recommended.

## 2.4 Connecting the viscometers and other devices

The ViscoPump III module allows the use of most various viscometer types:

DIN, ASTM, Ubbelohde and Micro Ubbelohde viscometers as well as Cannon-Fenske Routine, Micro TC and Micro Ostwald viscometers.

Owing to careful manufacture and quality-assurance procedures, all viscometers from SI Analytics® meet the highest accuracy standards.

The K viscometer constant is determined individually by way of a calibration of each glass capillary viscometer. Owing to the use of high-quality measurement and testing equipment and the application of national standard gauges, SI Analytics® guarantees an absolutely precisely reproducible calibration. For Ubbelohde viscometers having the same constant, the same correction seconds (Hagenbach correction) are valid.

**⚠ In the work mode "Press", in case of malfunction, the sample can emerge from the capillary tube and could enter the thermostat bath!** Therefore, especially during the pressing work mode, a correct connection of the measuring tripod to the ViscoPump III and a properly used viscometer are very important for the safe operation. With health hazardous and/or corrosive samples (e.g. acidic, toxic or hazardous to the environment), the work mode "suction" is the preferred choice.

**⚠ In the work mode "suction", the sample can be pumped into the ViscoPump in case of a malfunction!** In order to avoid this and the associated damage to the ViscoPump III, you must always use a suitable tube kit including the safety bottles VZ 7021 or VZ 7022 during the work mode "suction". In order to review a possible immersion of the sample into the safety bottle, a capacitive sensor VZ 8552 is highly recommended.

**⚠ It is also possible to connect or control other devices** (such as absorption traps, overflow guards etc.). Depending on the intended use of the ViscoPump III modul, it may be highly recommendable to connect these devices, please refer to the items below.

### 2.4.1 Viscometers using light-barrier sensing

Use the hose/cable combination to make an electrical and pneumatic connection between the ViscoPump III module VZ 8561 and the measurement stand. The plugs are firmly connected to the sockets by rotating the union sleeve. Please observe the color codes (red = sucking, black = venting) when screwing the threaded pneumatic connections into the ViscoPump III module. Please insert the selected capillary viscometer into the fixating rack as is shown in Fig. 5, then fill it. Insert the fixating rack together with the viscometer into the measurement stand (with the cut-out at the bottom sheet pointing forwards). The cut-out will latch into the lug provided. Pressing the viscometer slightly against the fixating rack will latch it into the holding spring located on the measurement stand.

**Fig. 5 Inserting or replacing a viscometer with light-barrier sensing**

### 2.4.2 TC viscometers with thermistor sensors

Fill the viscometer, then place it in the thermostat bath. Connect the ViscoPump III module VZ 8562 and the TC viscometer using the hose/cable combination which comes with the device. To do so, place the device in the holders, then attach the quadruple plug of the cable to the viscometer and the ViscoPump III module (first plug, then screw); subsequently, make the screwed connections in accordance with the numbers indicated on the hoses and the rack. In the case of “pressing” operation, the capillary tube remains open, for “sucking” operation the filling tube is to be left open. Please observe the color codes (red = suction, black = venting) when connecting the pneumatic screwed connections to the ViscoPump III module for TC viscometers.

### 2.4.3 Connection of VZ 7215 absorption traps

In the “sucking” mode (vacuum), volatile components can enter the ViscoPump III module. This is particularly problematic for corrosive solvents such as formic acid or dichloroacetic acid.

**⚠ For these cases, a hose fitting “suctioning” VZ 8524 must be used, which includes the absorption traps VZ 7215 and appropriate connecting hoses.**

The absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals. If sodium lime is used as an absorption agent with acidic solvents, the color condition of the indicator is to be checked on a daily basis. As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.

**⚠ If such a color shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a decoloration; this may then appear as “normal”, but it will definitely result in the destruction of the pneumatic system after some time.**

This situation is explicitly excluded from the warranty coverage!

For non-corrosive solvents and oils, which contain volatile constituents, absorption traps with activated carbon filling are available. When using activated carbon as an absorption agent (e.g. with solvents or used mineral oils), a replacement should be made at intervals between 1 and 2 weeks; this depends on the load factor which, in turn, is a function of the volatility of the materials.

### 2.4.4 Connection of the VZ 8552 overflow guard

We urgently recommend the connection of the VZ 8552 overflow guard (available as an option) for the suction-mode operation of the ViscoPump III module. The connection of the VZ 8552 overflow guard (capacitive sensor for the safety bottle) excludes over-pumping in suction mode (contamination of the ViscoPump III module). The holder on the safety bottle accommodates the capacitive sensor.

When using the ViscoPump III module VZ 8561 (meniscus sensing by light barriers) the holder for the safety bottle is to be attached to the measurement stand, e.g. the AVS<sup>®</sup>/S.

When using the ViscoPump III module VZ 8562 (thermo-resistive measurement) the holder for the safety bottle is to be attached to the “viscometer gallows” provided for the TC viscometer 5732.

Should any liquid be over-pumped into the safety bottle, the safety sensor will trigger a stop. After emptying the safety bottle, the lateral LED on the capacitive sensor will go out. You may continue with the measurements. The electrical connection of the VZ 8552 overflow guard is made using DIN plugs on the front side of the respective module of the ViscoPump III.

**⚠ The sensitivity of the capacitive sensor has to be adapted to the medium being used.** To do so, please use the enclosed screw driver to adjust the lateral set screw in such a manner that the capacitive sensor in the built-in condition (i.e. without medium) is just close from responding (i.e. the LED is on).

## 2.5 Troubleshooting

Check whether the Basic Measuring Device is switched on.

Trouble	Trouble shooting
<b>Air bubbles in viscometer</b>	<ul style="list-style-type: none"> <li>Is the filling quantity sufficient? Check, fill viscometer if required.</li> <li>Is the viscometer of properly connected?               <ul style="list-style-type: none"> <li>- in the case of <b>pressing mode</b>, please checked whether the filling tube is connected; if necessary connect properly.</li> <li>- for operation in <b>suction mode</b>, please check whether the capillary tube is connected; if necessary, connect properly.</li> </ul> </li> <li>please check whether the venting port is tightly connected; if necessary re-tighten its screwed connection.</li> </ul>
<b>Excessive pumping of measuring medium into the thermostat bath or into the safety bottle:</b>	<ul style="list-style-type: none"> <li>Is the viscometer properly connected?               <ul style="list-style-type: none"> <li>- for operation in pressing/suction mode?</li> </ul> </li> </ul> <p><b>When using AVS measuring tripods</b></p> <ul style="list-style-type: none"> <li>Check the electrical connections between the measuring tripod an the visco pump.               <ul style="list-style-type: none"> <li>- Correct: the green LED on the measuring tripod is illuminated.</li> </ul> </li> </ul> <p><b>When using TC viskometers:</b></p> <ul style="list-style-type: none"> <li>Is the viscometer properly connected?</li> </ul>
<b>ViscoPump III blinks (twice, continuously) and the sample in the viscometer is not pumped</b>	<ul style="list-style-type: none"> <li>Check the electrical connections between the measuring tripod and the module of the visco pump.               <ul style="list-style-type: none"> <li>- Correct: the green LED on the measuring tripod is illuminated.</li> </ul> </li> </ul>
<b>Problems with the controlling or serial interface</b>	<p><b>⚠ Attention</b></p> <p>With a simultaneous use of ViscoPump III (VP III) and the precursor devices ViscoPump (VZ 7001/ VZ 7002,VP I) and ViscoPump II (VZ 8511/ VZ 8512, VP II) at a viscosity measurement system with several measurement positions (AVS<sup>®</sup> 370 or AVS<sup>®</sup> Pro), the following must be observed for the installation:</p> <p><b>The VP III must be used on measurement positions following the VP I /VP II.</b></p> <p><b>Correspondingly, the VP III must be installed in the inserts to the right of the precursor devices!</b></p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px; text-align: center;"> <p>VP II    VP III</p> </div> <div style="border: 1px solid black; padding: 5px; text-align: center;"> <p>VP III    VP II</p> </div> </div> <p><b>Transfer errors will occur in the RS communication with a wrong sequence of the ViscoPump modules!</b></p> <p>(With all further connected devices with serial communication (e.g. burettes for rinsing the viscosity meters) the following must be adjusted as communication parameter:</p> <p>Baud:            4800            Data bit:        7            Stop bits:      2            Parity:          No</p>

## 3 Data transfer

### 3.1 RS-232-C Interfaces

The ViscoPump III module has three RS-232 C interfaces.

The interfaces on the internal data bus (96-channel VG bar) are used for communication with the respective controller unit or computer. The interface on the front panel serves to connect optional devices such as thermostats etc. The connected are treated as sub-devices of the ViscoPump III module.

### 3.2 Interface configuration

The interface parameters cannot be changed. All transmission parameters are firmly set to the following values:

Parity:	None
Stop bits:	2
Data bits:	7

Changes to the transmission parameters of the RS-232-C interface are only possible on a software basis. Setting is done through the program being used or through the software of the controlling device. Please refer to the operating instructions of the respective device or program being used for a description on how to make the required settings. The parameter settings of the connected peripherals must be identical, otherwise communication via this interface will be impossible.

All other devices to be connected must have the same parameter settings.

The other possible settings include:

<b>RS parameter</b>	<b>bit:</b>	<b>Stop:</b>	<b>parity:</b>
<b>baud:</b> 2400, 4800, 9600	7	2	No
2400, 4800, 9600	8	1/2	No/Even/Odd


## 4 Software-Update

The internal software (firmware) of the ViscoPump III can be updated if desired by the user by means of a USB stick. The PCB has a USB host port for this, see (see  Fig. 6).

Based on the design, a mini USB stick featuring a maximum length of 30 mm is required (including the plug). (Example: Intenso® Micro Line).

How to perform the software update:

1. Copy the firmware file (e.g. "VP3\_FW\_1732.bin") directly onto a FAT-formatted USB stick – without subdirectory (root). Do not rename the file. There must not be any other update file present in the root directory of the USB stick!
2. Switch off the viscosity meter (AVS® 370, AVS® 470, AVS® Pro) and pull the power plug. Pull the ViscoPump III out of its slot. Please observe the notes in Section 2.2.
3. Insert the USB stick into the USB host socket of the removed ViscoPump III.

**Fig. 6 ViscoPump III with inserted mini USB stick**

4. Return the ViscoPump III into its slot, plug in the power plug and switch on the power supply.
5. After the USB stick is detected, the loading of the update file will be indicated by the rapidly flashing of the status LED on the front panel. The rapid blinking will alternate with several short breaks.
6. A slow blinking sequence after approx. 10 seconds signals the successful update.

In case of an error, the status LED will blink very rapidly. The update will have to be performed once more. For this, shut off the power supply briefly in order to restart the ViscoPump III.

7. Switch off the power supply, pull the power plug, remove the ViscoPump III from its slot, pull out the USB stick, replace the ViscoPump III into its slot, insert the power plug and restart the ViscoPump III by switching on the viscosity meter.

## 5 Maintenance and care of the device and the viscometers

Maintaining the proper functioning requires certain inspection and maintenance work.

Maintenance and service work includes:

- Visual check
- Interface function, ViscoPump III, and rinsing burettes
- Once per quarter, the electrical contacts have to be inspected for corrosion, if the viscosity measuring unit is used in premises with an occasional occurrence of corrosive matters in their atmosphere.

### Maintenance intervals

Normal operation	intervals for carrying out all work is 6 months
Under particular strain	the max. intervals for carrying out all maintenance work are 4 weeks
In case of disturbances	If any disturbance, malfunction, or other defect becomes obvious, the work has to be carried out immediately.

### 5.1 Maintenance work to be carried out

- Check the hoses, screwed connections for signs of visible damage, contamination, and leaks.
- Check the electrical plug contacts for corrosion and mechanical damage (on the AVS<sup>®</sup> 370 Viscosity Measuring Unit and on the cables).
- If necessary, the exterior of the casing of the viscosity measuring unit can be cleaned with a piece of cloth soaked with a household cleaning agent. The lower and rear sections have to be dry-treated.
- In no case must liquid penetrate into the interior of the lower section.
- Defective parts must be repaired or replaced with new ones.
- Defective glass parts must always be replaced.

### 5.2 Maintenance and care of the VZ 7215 absorbent bottle

The VZ 7215 absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals.

**⚠ The measurements will not work if there is an incorrect connection.** There is a risk of the sample leaking from the viscometer or being suctioned into the ViscoPump III.

**⚠ If sodium lime is used acidic solvents as an absorption agent,** the color condition of the indicator is to be checked on a daily basis.

- As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.
- If such a color shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a de-coloration; this may then appear as "normal", but it will definitely result in the destruction of the pneumatic system after some time.

This situation is explicitly excluded from the warranty coverage!

**⚠ When using activated carbon as an absorption agent** (e.g. with solvents or used mineral oils), a replacement should be made at intervals between 1 and 2 weeks; this depends on the load factor which, in turn, is a function of the volatility of the materials.

### 5.3 Periods without operation

If the capillary viscometers are not used over a long period of time, the liquids contained in the system, in particular aggressive solutions, have to be drained. If the liquid is left in the system, one has to reckon that the solutions used will alter in the course of time and attack the glass, in particular the capillaries.

**⚠ Cleaning agents should be matched to the previous samples or impurities.**

In many cases, an aqueous cleaning agent (glass cleaners, detergents) or organic solvents (such as acetone or hydrocarbons) are sufficient.

**⚠ Strong oxidizing cleaning agents such as chromic acid may only be used by trained personnel and must be suitably disposed for safety and environmental reasons - the current guidelines for handling hazardous materials must be observed.**

**⚠ In the last rinse cycle, the viscometer should be rinsed with a suitable solvent with a low boiling point (such as acetone), and dried by an air flow, which is preferably generated by negative pressure (for example, water jet pump). The viscometer is dry and dust-free by this treatment and can thus be used for manual and automatic measurements.**

### 5.4 Reproducibility of results

The measurement or analysis results depend on a variety of factors. Please check the plausibility of the measurement results or analysis results at regular intervals, and carry out the required reliability tests. In this regard, please adhere to the usual validation procedures and especially to the "Viscometers within quality assurance systems" chapter.

### 5.5 Viscometers within quality assurance systems

Recommendations for companies that have introduced a quality management system (OM system) according to DIN EN ISO 9001: In this quality assurance system, an inspection of the measuring equipment is planned. The intervals and required accuracy can be defined by each company according to its own requirements. The standard DIN/ISO 10 012, Part 1 serves as a guideline in this matter. We recommend regular inspection of the viscometers in defined intervals.

**Inspection of the viscometer constants:**

#### a) Calibration using comparative measurements with reference measuring standards

Comparative measurements must be performed with a viscometer (reference measuring standard) which was tested at the PTB (Federal German Physical-Technical Institute) and provided with a constant. During this comparative measurement, the viscometer to be inspected and the PTB - tested viscometer were placed simultaneously in the same thermostat bath. The test liquid tested, the viscosity of which must not be known exactly, is filled into both viscometers, tempered and the flow-through time then measured. The constants of the viscometers to be inspected are then calculated according to the following equation:

$$K = \frac{K_{PTB} \times t_{PTB}}{t}$$

K constant of the tested viscometer

$K_{PTB}$  constant of the standard reference viscometer

t low time (HC) of the tested viscometer (corrected by Hagenbach-Couette)

$t_{PTB}$  flow time (HC) of the standard reference viscometer (corrected by Hagenbach-Couette)

Within the quality management system in accordance with DIN EN ISO 9001, traceability of the measuring equipment to national measuring standards is demanded. This traceability can be achieved by inspecting the comparative viscometers (reference measuring standards) at regular intervals at the PTB. The time intervals are defined according to the specifications made in the quality assurance quality management system of the user.

### **b) Calibration of the capillary viscometer with normal oils of the PTB**

During this calibration, a normal oil from the PTB with known viscosity is used as a reference measuring standard. The measurement is performed by means of flow-through measurement of the PTB normal oil in the viscometer to be inspected in a temperature bath, the temperature of which must correspond precisely to the test temperature of the PTB. In this case, it is extremely important to make sure that the temperature is absolutely correct. In case of temperature variation, this will always result in a constant for the viscometer that deviates from the constant applied. A temperature variation of 0.01 K, for instance, will result in a measuring error of up to 0.01 %. The calibration of the deviating temperature into the viscometer constant is not permitted.

### **c) Inspection by Xylem Analytics Germany with a quality certificate in accordance with DIN 55 350-18, 4.2.2**

The inspection at the manufacturer is carried out by means of comparative measurements using viscometers as reference measuring standards that were tested at the PTB (corresponds to Item 1).

#### **i General Information on the stability of viscometer constants**

Each inspection (even with a certificate) can guarantee the technical measuring direction only for a limited period of time. The constants of viscometers made of the borosilicate glass DURAN<sup>®</sup>, however, can remain unchanged for long periods of time if the viscometers are kept away from altered influences. Especially extreme changes can be expected, for instance, during the use of liquids that attack glass, in particular hot caustic soda hydrated (NaOH) orduring glass-blowing repairs (even for apparently insignificant repairs).

Liquids whose components adhere to the glass wall also cause errors. In such cases, regular cleaning is required whereby the corrosive action cleaning agent on the glass must be eliminated. For this reason, we recommend that the user should write up a special processing instruction for all important measurements and include them in his quality management manual in accordance with DIN EN ISO 9001. In all cases the user is responsible for the correctness of his measuring and testing equipment and is not released from his responsibility for quality (cp. DIN 55 350, Part 18).

## **6 Guarantee**

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts. The breach of glass parts is also excluded. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

## **7 Storage and transportation**

If the ViscoPump III has to be stored over some time, or to be dislocated, the use of the original packing will be the best protection of the devices. However, in many cases this packing will not be available anymore, so that one will have to compose an equivalent packaging system. Sealing the lower section in a foil is hereby recommended.

The devices should be stored in a room with a temperature between +10 and +40 °C, and the (relative) humidity of the air should not exceed 70 %.

 If the interchangeable have to be stored over some time, or to be dislocated, the fluids inside the system, especially aggressive solution have to be removed.

## 8 Recycling and Disposal

Please observe the applicable local or national regulations concerning the disposal of “waste electrical and electronic equipment”.

The ViscoPump III and his packaging are manufactured as far as possible from materials which can be disposed of environmental-friendly and recycled in a technically appropriate manner. If you have any question regarding disposal, please contact the service (see backside of this manual).

## 9 EC – Declaration of Conformity

The corresponding declaration of conformity of the device can be found on our homepage. It will also be made available to you on request.

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
**Xylem Analytics Germany Sales GmbH & Co.KG**
© 2022 Xylem, Inc. Version 220125 M 825 299 6