GEBRAUCHSANLEITUNG
Originalversion
OPERATING MANUAL
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES

AVS[®] 470

VISKOSITÄTSMESSGERÄT
VISCOSITY MEASURING DEVICE
L'APPAREIL DE MESURE DE LA VISCOSITE
EQUIPO MEDIDOR DE VISCOSIDAD

SI Analytics

a **xylem** brand

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Operating Manual Page 33 ... 62

Important notes:

The operating manual is part of the product. Before initial operation, please carefully read and observe the operating manual and keep it. For safety reasons the product may only be used for the purposes described in these present operating manual. Please also consider the operating manuals for the devices to be connected.

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties. A potentially more recent version of this manual is available on our internet website. The German version is the original version and binding in all specifications!

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TABLE OF CONTENT

1	Technical Specifications of the AVS[®] 370	35
1.1	Notes to the operating manual	35
1.2	Intended Use	35
1.3	Technical Specifications	35
1.3.1	AVS [®] 470	35
1.4	Warning and safety information	37
1.4.1	General	37
1.4.2	Chemical and biological safety	38
1.4.3	Flammable liquids	38
1.1	Functioning of the device	38
1.2	Capillary viscometry	39
1.3	Measurements principles	39
2	Installation and Commissioning	40
2.1	Unpacking and setting up	40
2.2	Connecting the devices	40
2.2.1	Connecting cables	40
2.2.2	Deployable hose combinations	41
2.2.3	Suitable viscometer types, racks, and measurement stands	42
2.3	Connecting the viscometers and other devices	44
2.3.1	TC-Viscometer with thermistor sensors	44
2.3.2	Viscometer with light-barrier sending	44
2.3.3	Connection of VZ 7215 absorption traps	45
2.3.4	Connection of the VZ 8552 overflow guard	45
2.3.5	Flow coolers	46
2.3.6	ViscoPump III module	46
2.3.7	System enhancements	46
2.3.8	Troubleshooting	47
2.4	Initialization and software upgrade of the AVS [®] 470 measuring unit	48
2.4.1	Initialization	48
2.4.2	Software upgrade	48
2.5	Description of the front-panel elements	48
2.6	Measurement with the device	49
2.7	Completing a measurement	49
3	Data transfer	50
3.1	RS-232-C Interfaces	50
3.2	Interface configuration	50
4	The Operating Software of the AVS[®] 470	51
4.1	Introduction	51
4.2	Hardware requirements	51
4.3	Operation	51
4.3.1	Operation using the PS2 keyboard	51
4.3.2	Operation using the front membrane keyboard	51
4.4	General information	52
4.4.1	Selection of menu items	52
4.4.2	Selection of the operating mode:	52
4.4.3	ViscoPump parameters:	52
4.5	Notes on programming	53
4.6	Functional description	53
5	Maintenance and care of the device and the viscometers	59
5.1	Maintenance work to be carried out	59
5.2	Maintenance and care of the absorber traps VZ 7215	59
5.3	Periods without operation	59
5.4	Reproducibility of results	60
5.5	Viscometers within quality assurance systems	60

6	Guarantee.....	61
7	Storage and transportation.....	61
8	Recycling and Disposal	61
9	EC – Declaration of Conformity.....	62
10	Appendix viscometer type list.....	62
11	Appendix spare parts.....	62

Copyright

© 2022, Xylem Analytics Germany GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization.
Germany. Printed in Germany.

1 Technical Specifications of the AVS[®] 370

1.1 Notes to the operating manual

The provided operating manual will allow you the proper and safe handling of the product. For maximum security, observe the safety and warning instructions in the operating manual!

-  Warning of a general danger:
Non-compliance results (can result) in injury or material damage.
-  Important information for device use.
-  Refers to another part of the operating manual.

The menu screens shown in this operating manual serve as an example and may differ from what you see!

1.2 Intended Use

The AVS[®] 470 is a measuring unit for determining absolute and relative viscosity of liquid media. This device must only be operated by skilled personnel.

It is operated using either the built-in membrane keyboard at the front or the TZ 2835 PS2 keyboard. Calculation of the results on the basis of the determined values is done by the built-in computing unit. The readings can be shown on the display or documented on the optional TZ 3460 printer.

1.3 Technical Specifications

1.3.1 AVS[®] 470

Translation of the legally binding German version

(Release: 18. Oct 2021)

EMC compatibility according to the Council Directive: 2014/30/EU;
applied harmonized standards: EN 61326-1
Low-voltage directive according to the Council Directive 2014/35/EU;
Testing basis EN 61 010-1: for laboratory equipment
RoHS Council Directive 2011/65/EU
FCC Part 15B and ICES 003

Country of origin: Germany, Made in Germany

Display: LCD-Type (70 x 40 mm)

Measurement parameters: Flow time in seconds [s]

Capture of measurement value:

Flow time: Optoelectronic or thermo-resistive capture of the meniscus passage through the measurement planes of the viscometers

Optional parameters: to be selected at the AVS[®] 470 measuring unit

Method: Absolute or relative viscosity

Viscometers: Ubbelohde viscometers (DIN, ASTM, micro); micro Ostwald; Cannon-Fenske routine; TC Ubbelohde viscometers and dilute-solution viscometers

ViscoPump: Pump parameters (ramp, pressure, suck over N1)

Time for temperature adaption: 0...20 min, to be selected in increments of 1 min

Number of measurements: 1...10 for each sample

Measurement ranges:

Viscosity: "pressing" mode 0.35 ... 1800 mm²/s (cSt) measuring temperature of approx. 20 ... 25 °C
"sucking" mode 0.35 ... 5800 mm²/s (cSt) measuring temperature of approx. 20 ... 25 °C

Pumping pressure: Fully automatically controlled "sucking" action to approx. -160 mbar
Fully automatically controlled "pressing" action to approx. +160 mbar

Measurement precision:

Precision (reproducibility and comparability) in accordance with DIN 51562, Part 1

Time measurement: $\pm 0.01 \text{ s} \pm 1 \text{ digit}$, but not more accurate than 0.01 %
 measurement uncertainty in the determination of absolute, kinematic viscosity furthermore depends on the uncertainty of the numerical value of the viscometer constants and the /measurement conditions, especially as concerns the measurement temperature.

Connections Back panel of the device:

Data in- and output: 1 serial RS-232-C interface: 9-channel sub-miniature D sockets:
 for use with a serial data printer
 1 device connector, round-plug connector with 7 channel DIN bayonet catch for connecting an automatic viscometer cleaning unit

Bath backlighting: Circular connectors with bayonet lock DIN 4 pin, 24 V, 350 mA

Mains connector: Device plug with safety interrupter according to VDE 0625, IEC 320/C14, EN 60320/C14, DIN 49 457 B

To be connected to the front panel of ViscoPump III modules:

Pneumatic connectors:

Aerate pressure/suction to be connected to viscometer.
 Overflow protection for suction hose VZ 8552

Capacitive sensor: DIN round-plug connector, 4-channel (ViscoPump III module)
 Screwed cap according to DIN 45321

ViscoPump III module:

DIN round-plug connector with bayonet catch
 - for AVS[®]/S (measurement stand), 5-channel or
 - for TC viscometer, 4-channel

Power supply:

Corresponds to Protection Class 1 according to DIN 57 411, Part 1 / VDE 0411, Part 1
 Mains connection: 100 - 240 V, 50 ... 60 Hz
 Mains fuse: Fine-wire fuse 5 x 20 mm, 250 V~, 4 A, time-lag design
 Power consumption: 100 VA

Casing:

Steel/aluminium casing with chemically resistant two-component coating, stackable

Dimensions: approx. 255 x 204 x 320 mm (W x H x D)

Weight: approx. 5.34 kg with 1 ViscoPump III module

Ambient conditions:

 **Do not use the device in hazardous locations!**

Climate: Ambient temperature: + 10 ... + 40 °C for operation and storage
 Humidity according to EN 61 010, Part 1:
 Max. relative humidity 80 % for temperatures up to 31 °C,
 linear decrease down to 50 % relative humidity at a temperature of 40 °C

Pollution degree: Pollution degree IP 20, indoor use only.

1.4 Warning and safety information

1.4.1 General

The device corresponds to protection class I.

It was manufactured and tested according to DIN EN 61 010, Part 1, “**Protective Measures for electronic measurement devices**” and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information contained in the present operating instructions. Development and production is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard.

 For reasons of safety, the device must only be used for the range of application described in the present operating manual. Nonobservance of the intended proper use of the device may result in personal injury or damage to property.

 For reasons of safety, the device and the power supply must be opened by authorised persons only; this means, for instance, that work on electrical equipment must only be performed by qualified specialists. **In case of nonobservance of these provisions the device and the power supply may constitute a danger: electrical accidents of persons or fire hazard!** Moreover, in the case of unauthorised intervention in the device or the power supply, as well as in the case of negligently or deliberately caused damage, the warranty will become void.

 Prior to switching the device on it has to be ensured that the operating voltage matches the mains voltage. The operating voltage is indicated on the specification plate (backside of the device). Only insert the power plug into a receptacle with ground contact. Any interruption of the conductor inside or outside of the device or the loosening of a protective conductor connector is not permitted and can lead to hazardous risk situations of the device. Only use fuses of the indicated type and rated amperage as a replacement. The use of repaired fuses or shorting the fuse box is prohibited. **Nonobservance of this provision may result in damage to the device and the power supply, or in personal injury or damage to property!**

 **If it has to be assumed that safe operation is impossible, the device has to be put out of operation and secured against inadvertent putting to operation.** In this case please switch the device off, pull plug of the mains cable out of the power supply, and remove the device from the place of work.

Examples for the assumption that a safe operation is no longer possible,

- if the package is damaged,
- if the device shows visible damages,
- if the power supply shows visible damages,
- if the device does not function properly,
- if liquid has penetrated into the casing.
- if the unit has been altered technologically or if unauthorized personnel tried or succeeded to open the device as attempt to repair it.

In case that the user operates such a device, all thereof resulting risks are on the user!

 The device must not be stored or operated in humid rooms.

 **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the environment or that they are supervised by specialized persons, respectively.

 When handling the substances used, the chemical resistance of the materials of the device must be taken into account.

 For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

 Please also note the corresponding Operating instructions for the devices to be connected!

1.4.2 Chemical and biological safety

i The device is not intended for use with potentially biohazardous substances.

⚠ **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the unit are experts in the handling of substances used in the devices or that they are supervised by specialized persons, respectively.

⚠ When using biohazardous substances, the regulations for handling the substances used must be observed. In such cases, the use is the sole responsibility of the user.

⚠ For all work with chemicals: **Always wear protective glasses!** Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

⚠ Dispose of all used solutions in accordance with national regulations and laws. Select the type of protective equipment according to the concentration and quantity of the hazardous substance at the respective workplace.

1.4.3 Flammable liquids

When handling flammable liquids, make sure that there is no naked flame in the vicinity of the equipment. Adequate ventilation must be provided. Only small quantities of flammable liquids should be kept in the workplace.

1.1 Functioning of the device

The AVS[®] 470 is used to perform flow-time measurements in capillary viscometers. It is operated via a Personal Computer.

The AVS[®] 470 is equipped with two options for meniscus scanning by using a corresponding ViscoPump III unit.

The available capillary viscometers enable viscosity measurements 0.35 of approx. 5,000 mm²/s (cSt) to be carried out at a measuring temperature of approx. 20 ... 25 °C.

Connecting TC viscometers to the ViscoPump III VZ 8562 module will also enable the measurement of black and opaque liquids. As an alternative option it is possible to use viscometers in combination with the light-optical ViscoPump III VZ 8561 module viscometers for meniscus sensing with the measurement stand, e.g. the AVS[®]/S.

Prior to the measurement as such, the liquid to be measured is sucked upwards inside the capillary viscometer through two measurement planes (N2 and N1) which are designed as light barriers or thermistor sensors, depending on the viscometer type (see  Fig.1 and Fig. 2).

The pumping pressure is controlled automatically by the AVS[®] 470 via the ViscoPump III module.

When using Ubbelohde viscometer, the design of the program ensures that the suspended spherical level will form prior to the start of the measurement.

Time recording extends up to 9999.99 seconds with a resolution of 0.01 s. The individual results of a measurement series can either be presented on the display, but it is also possible to document them on the printer which is available as an option.

The AVS[®] 370 can also be operated with ViscoPump II - modules, the older ViscoPump generation.

1.2 Capillary viscometry

Capillary viscometry is the most accurate method for the determination of the viscosity of liquids with a Newtonian flowing behaviour. The measurement as such consists in a time measurement. The time measured is that which a specific quantity of liquid requires to pass through a capillary having a defined width and length. Conventionally, this process is watched with the human eye, and the flow time is measured manually using a stop watch.

In the case of the AVS[®] 470, as with all viscometry measuring devices from SI Analytics[®], the liquid meniscus is captured on the measurement planes, either in an optoelectronic manner by means of light barriers, or else on a thermo-resistivity basis by thermistors.

1.3 Measurements principles

a) Optoelectronic sensing of the liquid meniscus

The near-infrared light which is generated in LEDs located in the upper section of the measurement stand AVS[®]/S is conducted through a glass-fibre light-conductor cable onto the measurement planes. The light shines through the viscometer before it arrives at another light-conductor cable located on the opposite side; inside this second cable, the light is conducted to a receiver in the upper section of the measurement stand. While the liquid meniscus passes through the measurement planes, the lens-like effect of the meniscus causes a short-term darkening of the light beam, followed by a magnification. This process generates a measurement signal which can be evaluated accurately.


Fig.1 Viscometer for optoelectronic measurements

b) Viscometer with thermistor sensors (TC viscometer)

In the case of TC viscometers, glass-coated thermistors serving as sensors are inserted on the level of the measurement planes. While the meniscus passes through the measurement planes, the differences in the thermal conductivity properties of air and liquid lead to a change in the heat balance. The thermistors of the TC viscometers are inserted hermetically tightly into the glass coating of the viscometer.


Fig. 2 Viscometer for optoelectronic measurements

2 Installation and Commissioning

2.1 Unpacking and setting up

i The installation of the AVS[®] 470 and the connection of the desired additional devices are generally conducted by a trained service technician.

⚠ **Please observe the operating voltage!** It is indicated on the type plate (back side of the device).

The device may be placed and operated on a plane surface.

i We recommend placing on the VZ 8571.

Up to two devices may be stacked.

2.2 Connecting the devices

2.2.1 Connecting cables

Description	Length	Connects:	to:
TZ 3086 Printer cable	1.5 m	AVS [®] 470	Printer Z 910

2.2.2 Deployable hose combinations


Hose combinations Type no.	Description	Application
VZ 5505	Silicone hose kit, oppressive, for Ubbelohde viscometers (3 legs), and Cannon-Fenske and Ostwald viscometers	Standard, but please observe: the sample can leak from the capillary tube during a malfunction
VZ 5505 + VZ 8526	Silicone hose kit, suctioning, for Ubbelohde viscometers (3 legs)	Standard, safer than oppressive operation, as the sample cannot leak from the capillary tube. Unsuitable for volatile samples.
VZ 8523	PTFE hose set, suction, for Ubbelohde viscometer (3 legs)	For aggressive samples that attack silicone, such as sulphuric acid. Hose lengths are designed to fit the AVS [®] 470 on the support console VZ 8571. All hose sets for suctioning operation can be combined with the suction set and the sample fill set.
VZ 8524	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (3 legs)	For aggressive samples whose vapors are absorbed by soda lime filters to protect the ViscoPump. Depending on the sample, an active carbon filter VZ 7216 can be used instead of the soda lime filter VZ 7215. Hose lengths are designed to fit the AVS [®] 470 on the support console VZ 8571.
VZ 7218 + VZ 8535	PTFE hose kit, suctioning, for Ubbelohde viscometers (4 legs)	For Ubbelohde viscometers with additional 4 th pipe for filling and cleaning Hose lengths are designed to fit the AVS [®] 470 on the support console VZ 8571.
VZ 8530	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (4 legs)	For Ubbelohde viscometers with additional 4 th pipe for filling and cleaning Hose lengths are designed to fit the AVS [®] 470 on the support console VZ 8571. For aggressive samples whose vapors are absorbed by soda lime filters to protect the ViscoPump.
VZ 5606	For TC viscometers (3 and 4 legs) with screw connections: Silicone hose fittings with connecting cable. For oppressive operation	For TC viscometers (3 and 4 legs) with 4 th pipe for filling and cleaning Typical applications are measurements of oils.

i The hose combinations are to be selected on the basis of the required application.


⚠ When using Micro Ubbelohde viscometers with TV sensors, the ignition temperature of the media to be measured has to be taken into account! It has to be higher than 250°C.

2.2.3 Suitable viscometer types, racks, and measurement stands


Viscosimeter (Type)	Rack (Type no.)	Measuring stand (Type)	
Ubbelohde (DIN) 532... 530... 501... 541... 545...	053 92	AVS [®] /S	AVS [®] /SK
Ubbelohde (ASTM) 525... 526... 527... 545...	053 92	AVS [®] /S	AVS [®] /SK
Mikro Ubbelohde 536... 537... 538...	053 92	AVS [®] /S	AVS [®] /SK
Ubbelohde dilute-solution viscometers 531...	---	AVS [®] /SK-V	
Cannon-Fenske-routine 513... 520...	---	AVS [®] /SK-CF	
Mikro-Ostwald 516... 517... 518...	053 97	AVS [®] /S	AVS [®] /SK
Ubbelohde (TC) 582... 583... 584... 587... 588... 589...	053 93	---	

DIN-Ubbelohde-viscometer
type 1 (DIN) resp. 2 (ASTM)


Micro-Ubbelohde-viscometer
type 3


Cannon-Fenske-Routine-viscometer
type 4

Micro-Ostwald-viscometer
type 5

Ubbelohde-viscometer
with TC sensor type 1

Micro-Ubbelohde-viscometer
with TC sensor type 3

Fig. 3 Suitable viscometer types

i The respective viscometers are assigned type numbers which are to be entered when programming the AVS 470 (📖 sections 2.3, 4.610 and 10)

2.3 Connecting the viscometers and other devices

The AVS[®] 470 allows the use of most various viscometer types:

DIN, ASTM, Ubbelohde and Micro Ubbelohde viscometers as well as Cannon-Fenske Routine, Micro TC and Micro Ostwald viscometer

Owing to careful manufacture and quality-assurance procedures, all viscometers from SI Analytics[®] meet the highest accuracy standards.

The K viscometer constant is determined individually by way of a calibration of each glass capillary viscometer. Owing to the use of high-quality measurement and testing equipment and the application of national standard gauges, SI Analytics[®] guarantees an precisely reproducible calibration.

For the flow time, a correction time can be determined based on the viscometer constant and the viscometer design, expressed by the type no. (Fig. 3). See also  chapters 4 and 10. This statement is true for both Ubbelohde viscometers of normal size as well as for micro viscometers.

i It is also possible to connect or control other devices (such as absorption traps, overflow guards etc.). Depending on the intended use of the AVS[®] 470, it may be highly recommendable to connect these devices, please refer to the items below.

The AVS[®] 470 is not intended to be connected to a computer, piston burettes or similar devices working according to the daisy-chain principle of SI Analytics[®]. If you wish to use a computer, please select the AVS[®] 370 Viscosity Measuring Unit from SI Analytics[®]; in this case, please contact SI Analytics[®] for the corresponding documentation.

In order to better follow the measurement process in the bath, the optionally available bath backlighting VZ 5405 can be connected to the AVS[®] 470. In addition, the adapter cable VZ 5408 is required for this purpose.

2.3.1 TC-Viscometer with thermistor sensors

Fill the viscometer and insert it in the thermostat bath. Connect the AVS[®] 470 Viscosity Measuring Unit and the TC viscometer using the hose/cable combination which comes with the device. To do so, place the device in the holders, then attach the quadruple plug of the cable to the viscometer and the ViscoPump III module VZ 8562 (first plug, then screw), subsequently, make the screwed connections in accordance with the numbers indicated on the hoses and the rack. In the case of “pressing” operation, the capillary tube remains open, for “sucking” operation the filling tube is to be left open. Please observe the colour codes (red = sucking, black = venting) when attaching the pneumatic screwed connections to the ViscoPump III module of the AVS[®] 470.

2.3.2 Viscometer with light-barrier sensing

Use the hose/cable combination to make an electrical and pneumatic connection between the AVS[®] 470 and the measurement stand. The plugs are firmly connected to the sockets by rotating the union sleeve. Please observe the colour codes (red = sucking, black = venting) when screwing the threaded pneumatic connections into the ViscoPump II module VZ 8561. Please insert the selected capillary viscometer into the fixating rack as is shown in Fig. 4, then fill it. Insert the fixating rack together with the viscometer into the measurement stand (with the cut-out at the bottom sheet pointing forwards). The cut-out will latch into the lug provided. Pressing the viscometer slightly against the fixating rack will latch it into the holding spring located on the measurement stand.

Fig. 4 Inserting or replacing a viscometer with light-barrier sensing

2.3.3 Connection of VZ 7215 absorption traps

i In the “sucking” mode (vacuum), volatile components can enter the ViscoPump III module. This is particularly problematic for corrosive solvents such as formic acid or dichloroacetic acid.

⚠ **For these cases, a hose fitting “suctioning” VZ 8524 or VZ 8527 must be used!**
(includes the absorption traps VZ 7215 and appropriate connecting hoses).

i The absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals. If sodium lime is used as an absorption agent with acidic solvents, the colour condition of the indicator is to be checked on a daily basis. As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.

⚠ **If such a colour shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a decolouration; this may then appear as “normal”, but it will definitely result in the destruction of the pneumatic system after some time!** This situation is explicitly excluded from the warranty coverage!

For non-corrosive solvents and oils, which contain volatile constituents, absorption traps with activated carbon filling are available. When using activated carbon as an absorption agent (e.g. with solvents or used mineral oils), a replacement should be made at intervals between 1 and 2 weeks; this depends on the load factor which, in turn, is a function of the volatility of the materials.

2.3.4 Connection of the VZ 8552 overflow guard

We urgently recommend the connection of the VZ 8552 overflow guard (available as an option) for the suction-mode operation of the ViscoPump III module. The connection of the VZ 8552 overflow guard (capacitive sensor for the safety bottle) excludes over-pumping in suction mode (contamination of the ViscoPump III module). The holder on the safety bottle accommodates the capacitive sensor.

When using the ViscoPump III module VZ 8561 (meniscus sensing by light barriers) the holder for the safety bottle is to be attached to the measurement stand, e.g. the AVS[®]/S. When using the ViscoPump III module VZ 8562 (thermo-resistive measurement) the holder for the safety bottle is to be attached to the “viscometer gallews” provided for the TC viscometer 5932.

i Should any liquid be over-pumped into the safety bottle, the safety sensor will trigger a stop. After emptying the safety bottle, the lateral LED on the capacitive sensor will go out. You may continue with the measurements.

The electrical connection of the VZ 8552 overflow guard is made using DIN plugs on the front side of the respective module of the ViscoPump III.

⚠ **The sensitivity of the capacitive sensor has to be adapted to the medium being used.**

To do so, please use the enclosed screw driver to adjust the lateral set screw in such a manner, that the capacitive sensor in the built-in condition (i.e. without medium) are just close from responding (i.e. the LED is on).

Viscosity depends on the temperature of the sample liquid. This means that the viscometers must always be thermostated during the measurement. The measurement temperature has to be kept constant in order to achieve an accurate result.

The transparent thermostats from SI Analytics[®] which were developed especially for capillary viscometry, meet the requirements imposed on precision and constancy. The CT 72/2, CT 72/4 thermostats, for instance, guarantee a temperature constancy of $\pm 0,02$ K at a command temperature in the range of 10° to 40 °C, and a maximum fluctuation of the ambient temperature of ± 3 K.

i As a rule of thumb, you may suppose that the temperature deviation, expressed in degrees, multiplied with a factor of 10 will correspond to the deviation from the result in terms of %. This means that a deviation of 0.05 °C corresponds to a possible error of 0.5 %.

In principle, two different transparent thermostats can be used on the AVS[®] 470:

For measurements at different temperatures, the CT 72/2, CT 72/4 transparent thermostats are available. These thermostats can be equipped with two or four viscometers including other equipment. For measurements up to 60 °C, the thermostat CT 72/P of acrylic can be used.

⚠ The bath body of the thermostat CT 72/P consists of acrylic (PMMA)! Acrylic can be damaged by a variety of organic solvents as well as by concentrated acids/bases. Therefore, a contact of the bath body with such substances should be avoided, as the acrylic may be damaged.

i Please read the separate operating manual of the transparent thermostats as well.

2.3.5 Flow coolers

⚠ As was mentioned above, viscosity measurement is highly dependent on temperature constancy. For reasons of control technology (self heating of the thermostat head), it is therefore necessary to use a CK 300/CK 310 flow cooler as a counter cooler at bath temperatures exceeding 40°C.

i Please read the separate operating manual of the flow coolers as well.

2.3.6 ViscoPump III module

The ViscoPump III modules control the entire measurement process, among other things the temperature pre-adaptation process of the samples in the viscometers, the process of pumping the liquid up into the storage containers of the viscometers, the measurement of the flow times etc.

To replace the ViscoPump III module, please proceed as follows:

- Switch off the AVS[®] 470 and unplug the power plug from the socket.
- Remove the pneumatic and electrical connections from the front panel of the ViscoPump III module to be replaced.
- Loosen the screws located at the corners of the front panel.
- Use the upper and lower insertion handles to leverage the ViscoPump III module out of its rearside plug connection.
- Pull the ViscoPump III module out of the AVS[®] 470.
- After inserting the new ViscoPump III module, please secure it again with the screws of the front panel. Re-establish the electrical and pneumatic connections.
- Insert the mains plug into the socket and switch on the AVS[®] 470.

⚠ After connecting, check whether the correct operating mode, "Suctioning" or "Pushing" is visibly set on the corresponding LED on the front panel. Changing the operating mode to the application is described in  section 4.6, Figure 23. Use the hose fittings designed for this purpose!

2.3.7 System enhancements

A number of functional units of the AVS[®] 470, such as viscometers and the ViscoPump III module, can be replaced. The system can be extended by an optional Z 910 printer.

⚠ Make sure that the AVS[®] 470 is always the first device you switch off!

Prior to replacing any functional unit, please be sure to the mains plug must be pulled out of the mains socket. **Caution: Liquid dripping off may be hazardous to the user!**

2.3.8 Troubleshooting

Check whether the AVS[®] 470 is switch on.

Trouble	Troubleshooting
Air bubbles in viscometer	<ul style="list-style-type: none"> • Is the filling quantity sufficient? <ul style="list-style-type: none"> - Check, fill viscometer if required. • Is the viscometer of properly connected? <ul style="list-style-type: none"> - in the case of pressing operation, please checked whether the filling tube is connected; if necessary connect properly. - for operation in suction mode, please check whether the capillary tube is connected; if necessary, connect properly. - please check whether the venting port is tightly connected; if necessary re-tighten its screwed connection.
Excessive pumping of measuring medium into the thermostat bath or into the safety bottle	<ul style="list-style-type: none"> • Is the viscometer properly connected? <ul style="list-style-type: none"> - for pressing operation? - for operation in suction mode? <p data-bbox="577 813 975 842"><u>When using AVS[®] measuring tripods</u></p> <ul style="list-style-type: none"> • check position of the rack within the stand • check the electrical connection from the viscometer to the ViscoPump type III module • Is the green LED on the measuring stand illuminated? <p data-bbox="577 983 884 1012"><u>When using TC viskometers</u></p> <ul style="list-style-type: none"> • Is the viscometer properly connected?

2.4 Initialization and software upgrade of the AVS® 470 measuring unit

2.4.1 Initialization

At the time of delivery, all memory values of the AVS® 470 Viscosity Measuring Unit are set to starting values (so-called default values). If it should prove necessary to restore this delivery status at any time, this can be achieved in the form of a so-called initialisation. This process is triggered by simultaneously pressing the “Up ↑” and “Down ↓” key for at least 2 seconds as shown in Fig. 5 after switching the device off and on again, initialisation is completed.

2.4.2 Software upgrade

The software of the AVS® 470 Viscosity Measuring Unit can be upgraded, but such an upgrade has to be done by specially trained service personnel. If such an upgrade should become necessary, please contact SI Analytics® (address can be found at the end of the present operating instructions) to take the required action.

2.5 Description of the front-panel elements


Fig. 5 Front panel

2.6 Measurement with the device

Setting the measurement parameters for a method

(Sample) description (optional)	< ▲ ▼ ↵ >
User (optional)	< ▲ ▼ ↵ >
Lot (optional)	< ▲ ▼ ↵ >
Measurement type [Abs., Saybolt, rel., blank value]	< ▲ ▼ ↵ >
Number of measurements [1 ... 99]	< ▲ ▼ ↵ >
Pre. temp. time [1...20 min]	< ▲ ▼ ↵ >

Input of the temperature pre-adaptation time in minutes. In the course of the temperature pre-adaptation, the liquid is permanently pumped upwards and flows through the viscometer, just as in the course of a measurement; this process is intended for a speedy temperature adaptation. The input of the temperature pre-adaptation time may be as high as 20 minutes.

Bath temperature	[°C]	< ▲ ▼ ↵ >
Max. aberration	[%]	< ▲ ▼ ↵ >
Viscometer ID	[1 digit]	< ▲ ▼ ↵ >
t ₀ time	[s]	< ▲ ▼ ↵ >
Constant	[mm ² /s ²]	< ▲ ▼ ↵ >

(For the further way of proceeding, see  section 4.6 Functional description)

2.7 Completing a measurement

Depending on the specific circumstances, there is a number of conditions which lead to the completion or cancellation of the measurement program in the AVS[®] 470 Viscosity Measuring Unit:

- Automatically, as soon as the specified number(s) of measurements including the repeated measurement has been performed; this is the normal case.
- Abortion of the respective measurement series by selecting “Stop”. In this case all device measurement values will be lost; this feature should only be used in emergency cases.
- Abortion of the respective measurement series by a time-out error. In this case all device measurement values will be lost, and the measurement is to be restarted.
- Mains failure: In the case of a power failure all device measurement values will be lost. After the return of the mains current the condition will be as it was after parameterisation, prior to the measurement. The set and stored parameters will remain preserved in the EPROM.

3 Data transfer

3.1 RS-232-C Interfaces

The AVS® 470 has an RS-232-C interface. This interface is used for documentation in connection with the optional data printer e.g. Z 910.

3.2 Interface configuration

All transmission parameters are firmly set to the following default values:

Parity:	none
Stop bits:	1
Data bits:	8

It is essential to set the printer to the same parameter settings. The optional Z 910 printer from SI Analytics® is ex-works set to these parameters.

The other possible settings include:

RS-Parameter:

Baud:	Bit:	Stopp:	Parity:
2400	7	2	No
4800	8	1	No
9600	7	1	Odd
	8	1	Odd
	7	1	Even
	8	1	Even

i.e. the baud rate applicable to the respective parameter sets can be set independently (see  section 4.6, figures 27 and 28).

4 The Operating Software of the AVS[®] 470

4.1 Introduction

The AVS[®] 470 Viscosity Measuring Unit including the ViscoPump II module are controlled by the operating software of the AVS[®] 470. The measurement values determined by the ViscoPump II module are received and evaluated by the software. The results may be output on a printer (Report) and shown on the display. The user can store the various measurement parameters of a method. Below please find a description of the functioning and operation of the software. The selection from the various options is done using the cursor, followed by a depression of the "Enter" key.

4.2 Hardware requirements

For the operation of the AVS[®] 470 we recommend the following minimum equipment:

- 1.) keyboard (TZ 2835)
- 2.) printer (RS-232-C) serial (e.g. Z 910)
- 3.) ViscoPump III module measurement plug-in unit (VZ 8561 or VZ 8562)

4.3 Operation

4.3.1 Operation using the PS2 keyboard

All the functions described in  section 4.6 can be called using the keyboard.

"Enter"	= To confirm input and continue with the programming cycle
"Esc"	= Escape, backward jump to the previous screen
"↓"	= To scroll backwards
"↑"	= To scroll forwards
"←"	= Move to the left
"→"	= Move to the right
"F1"	= Go to "start"
"F2"	= Go to "stop/reset"
"F3"	= selection on the main menu, refer to  section 4.6 Fig.4
"F4"	= Go to blank-value determination
"F5"	= Call measurement values
"F6"	= Result protocol; you may use the F6 function key to create a new printout of the result protocol as long as the method creating the protocol was not changed, and unless the AVS [®] 470 Viscosity Measuring Unit was switched off; if no printer is used, the protocol can be viewed on the display
"Pg Up"	= Increase LCD contrast
"Pg Dn"	= Decrease LCD contrast

For any input use the numeric and character keys.

4.3.2 Operation using the front membrane keyboard

The front membrane keyboard may be used for starting and stopping the program, or for restricted programming operations, i.e. only such programming steps can be modified or selected which do not require any data to be input, but just a selection to be made; please refer to  section 4.6.

"Enter"	= To confirm input and continue with the programming cycle
"Esc"	= Escape, backward jump to the previous screen
"↓"	= To scroll backwards
"↑"	= To scroll forwards
"Stop"	= Go to stop/reset
"Start"	= Start the program

4.4 General information

4.4.1 Selection of menu items

After selecting the required menu item using the selection bar = cursor in combination with the ↓ and ↑ arrow keys, the item will be highlighted in black. Pressing the “Enter” key will confirm your selection.

Input of values: Any values to be input are to be entered in the corresponding empty or default-value fields. After selecting the corresponding fields using the selection bar = cursor in combination with the ↓ and ↑ arrow keys, the respective field will be highlighted in black. If the upper or lower limits are ignored in the selection, you cannot proceed by simply pressing the “Enter”. The default value will be displayed.

Please note: If a numeric field (other than a floating-point field) allows the input of e.g. 3 digits, but only 1 or 2 of the digits to be input are of significance, a trailing “0” has to be added.

Example: The value 100 is to be changed to 30; the input would be 0 - 3 - 0; on the display the value will be shown as 030.

Storing the values or selection: The values entered or changed, or the selection made, respectively, are taken over and memorized only after an explicit confirmation. Return to previous screen: If you wish to return to the preceding screen without saving any changes, please press the “Esc” key.

4.4.2 Selection of the operating mode:

( Section 4.6, fig. 23) A selection can be made between “sucking” and “pressing” action. After power-up, the selected operating mode will be signalled in the form of the letter “P” or “S” on an indicator light (red LED) on the front panel.

 The user is responsible for using the appropriate hose sets and connecting them properly!

4.4.3 ViscoPump parameters:

( Section 4.6, fig. 23) ViscoPump parameters, you may set the “Ramp”, “Pump power”, “Meas. delay” and “Suction above N1” parameters in addition to the operating mode. Legend:

Ramp: This is the gradient angle applicable to the individual pump-pressure increases (dynamic pumping / suction); the default setting is 15 (fictitious measure).

Pump power: The max. achievable performance, the default value = 30%!

Meas. delay: A factor which may be added as a holding period to the period of time which was calculated from the flow time and which is to be observed between two subsequent measurements in order to shorten this time out of a system-inherent necessity, or to increase it, for instance, to allow the capillary to run empty.

Suction above N1: A factor which may be used to shorten or increase the period of time which was calculated from the flow time as required for sucking/pumping above the upper N1 level, if this seems to be appropriate out of a system-inherent or application-technological necessity.

The default values of these parameters are selected in such a manner that a large portion of all applications can be carried out without any problem. In special cases, for instance in the case of very short run times in micro viscometers, or in the case of very high toughness, or also for optimising the total throughput, these parameters may require an adaptation.

 **Please note that such an adaptation has to be made only in small steps and with due care. Any improper modification may cause damage to the ViscoPump III module of the AVS[®] 470 and lead to the loss of warranty.**

4.5 Notes on programming

Below you find a description of how to program the AVS[®] 470 Viscosity Measuring Unit. In principle, the settings are to be made as in the case of other automatic viscosity measurement devices and software packages from SI Analytics[®]. Users who are conversant with such devices will not encounter any problems with the AVS[®] 470 Viscosity Measuring Unit. Other users who are working for the first time with devices of this kind should at first try out the various settings to gather some experience on their influence on the measurements.

In the delivery state, all the settings are present in the form of so-called default or basic settings to which they will also be restored in the case of a deletion of the memory contents.

In the operating mode used to create a method or to set the parameters of the ViscoPump or the system parameters, you can use the “ESC” key to return to the preceding screen.

4.6 Functional description

After power up, the following screen will appear on the display:


Fig. 1 Power-up screen for 5-10 s
After the power-up screen has disappeared, either of the two screens below may appear:


Fig. 2 System ok,
press “Enter” to proceed to **Fig. 4**


Fig. 3 System is not ok,
follow the error removal instructions; until the elimination of the problem was confirmed.

If everything is ok, you you will be taken to **Fig. 2**,
As long as the error persists, you will be hold at **Fig. 3**

If the ViscoPump III is operated without a capacitive sensor for monitoring the safety bottle (recommended for suction working mode), the supplied dummy plug must be plugged into the “Safety bottle” socket of the ViscoPump. Otherwise, the error message “Safety bottle full” appears as shown.

method: absolute
ready!
start
create method
delete method
system paramete
ViscoPump parameter

Id: Sample XYZ
lot: 12075ADC
usr: Jonny Miller

temp equ. Time
122 s

meas. 1 out of 10
120 s

meas. 2
245.56 s

result corr.
average = 234.56 s
stand. dev. = 0.001
correction = 0.34 s
corr. ave. = 233.22
Abs Visc = 3.322 mm ² / s

Fig. 4. Main menu, selection

Display of the set method (e.g. absolute)

Start the set method

Fig. 5

Create/change the set method

Fig. 14

Delete the method

Fig. 32

Change/set the system parameters

Fig. 25

Change&/set the ViscoPump parameters Fig. 23

Fig. 5 Input after start

Id: The designation of a method; in the present case the field is empty (after a reset, or yet unused) or contains the last input

lot: The lot designation (as above)

usr.: User (as above)

No mandatory entry, press "Enter" to proceed

Fig. 6

If "pre. temp." (fig. 17) was set, this screen will display the remaining time of the running temperature pre-adaptation in terms of seconds.

Fig. 7 Measurement after temperature pre-adaptation

After temperature pre-adaptation, the set number of measurements (Fig. 17) will be preformed and displayed here on a running basis in steps of 1 s.

Fig. 8 Measurement value at the end of a measurement

Upon completion of a measurement, the respective result will remain on the display until the next measurement begins.

Fig. 9 Display of the measurement result

Example: Absolute measurement with correction

The present example shows the result of an absolute measurement in which the Hagenbach-Couette correction was applied. The average shown corresponds to the **(exemplary!)** amount, reduced by the number of correction seconds which are also shown.

Result n. corr.
Average = 234.56 s

SUS = 356
Temperature = 100° F

Result n. corr.
Average = 234.56 s

SFS = 234
Temperature = 250° F

Result n. corr.
Average corr. = 234.56 s

eta rel. = 1.23456
eta spec. = 0.23456
V.N. = 234.56 [ml / g]

Retrieve measurement values
Main menu

measurements [s]
meas. 1: 1234.67 *
main menu
press Enter

In the „Creation mode“, you can at any time return to the previous screen by pressing the “Esc” key

Create method

mode: absolut
Saybolt
relative
blank value

Fig. 10 Display of the measurement result
Example: Absolute measurement with calculation of SUS (Sayboldt Universal seconds) without Hagenbach-Couette correction (exemplary!)

Fig. 11 Display of the measurement result
Example: Absolute measurement with calculation of the SFS (Sayboldt Furol seconds) without Hagenbach-Couette correction (exemplary!).

Fig. 12 Display of the measurement result
Example: Relative measurement with correction. This example shows the result of a relative measurement with the application of the Hagenbach-Couette correction. The average shown represents the amount, reduced by the correction seconds which are also shown (**exemplary!**)
The VN (viscosity number) is determined using the concentration which was input in Fig. 19. As soon as you confirm the result with Enter, you will be taken to the “**Retrieve Measurements values**” mode.

Fig. 13 Retrieve measurement values
If the measurement mode was selected, you can retrieve the individual measurement values one by one at this point. If this is not needed, you can return to the starting screen (**Fig. 4**).

Fig. 14 Retrieve and scroll through measurement result
Use the ↑ and ↓ key to retrieve the individual measurement results. The “*” sign means that this value was used in the calculation. To finish, press the “Enter” key to return to **Fig. 4**.

Fig. 15 Creation mode following Fig. 14
Selection from 4 modes: With the “Enter” key If you select “absolute” and “relative”, you will proceed to **Fig. 16**. Selecting “Saybolt” will take you to **Fig. 18** Select “Blank value” to proceed to **Fig. 21**

Create method

temperature: 25.00
° F ° C

Create method

viscometer type:
constante:
of meas.:
pre temp. time:
delta % choice:
H.C. correktion:

Create method

temperature: 100° F
SFS SUS

Create method

eta rel: 1
eta spec. 1
V.N.: 0
conc. 0,250
dim.: [g/ml] [g/dl]

Create method

save parameters?

Yes No

Relative measurement:

blank value: 0,00 s
measure blank value ?

Measure blank value:

result corr.
average = 1234.56 s
H.C. corr. = 1.23 s
press Enter

Fig. 16 Creation mode following Fig. 15
Enter the temperature value for documentation and select the temperature scale (°C or °F), if “abs” or “rel” mode was selected – press the “Enter” key to proceed to Fig. 17 or to Fig. 19 in “relative” mode.

Fig. 17 Creation mode following Fig. 16
Select of the viscometer type: see  chapter 10
Input of the constants or guideline constants
Input of the number of measurements
Input of the temperature pre-adaptation time
Input of the maverick test:
“0”: no, number ≠ “0”: yes
Selection criterion in ± n.nn%: max. permissible deviation of transit times from mean value, if exceeded: substitute measurement
Selection of HC correction, 0 = no, 1 = yes
Press “Enter” to proceed to selection in Fig. 20

Fig. 18 Creation mode following Fig. 16
If Sayboldt calculation was selected in Fig. 15, this is the point to enter the working temperature in terms of °F; subsequently, you have to specify SUS or SFS for the calculation in the selection field below. Press “Enter” to proceed to Fig. 17

Fig. 19 Creation mode following Fig. 17
If “relative” was selected in Fig. 15:
Selection of the calculation types:
0 = no selection, 1 = selection.
Please note: eta rel. includes eta rel. and VN eta rel. and eta spec.
Input of the concentration, and selection of its unit; only one value is possible.
Press “Enter” to proceed to Fig. 17

Fig. 20 Creation mode, prompt following Fig. 17
This prompt is issued for data integrity reasons to avoid the automatic take-over of any erroneous input or changes.
Yes means: accept and store data
No means: discard changes
Press “Enter” to return to Fig. 4 or, in “Relative” mode, to Fig. 21 if “yes” was selected.

Fig. 21 “Relative” creation mode
If “Relative mode” or “Blank value” was selected in Fig. 15, you will be prompted here to know whether the blank value (t_0) is to be entered manually, or whether it is to be measured instead.
The measurement will be made using the parameters which were input in Fig. 17. If manual input is specified you will proceed to Fig. 4, Measurement, pressing “Enter” will take you to Fig. 22.

Fig. 22 “Relative” creation mode, measurement result
As soon as the measurement of a blank value is completed, the result will be displayed; by “Enter” you can accept the result as t_0 ; subsequently, you will be taken to Fig. 4.

ViscoPump parameter:

pressure	suction
ON	OFF
ramp:	30 %
pump power:	100 %
meas. delay:	1.0 f
time above N1:	1.0 f

ViscoPump parameter:

save parameters?

Yes No

System parameter

language:
RS parameter:
documentation:
date/time:
back:

System parameter

language: deutsch
 english
 francais
 espanol
 italiano

System parameter

RS parameter:
 2400 Baud
 4800 Baud
 9600 Baud

RS parameter:

Bit	Stopp	Parity
7	2	No
8	1	No
7	1	Odd
8	1	Odd
7	1	Even
8	1	Even

Fig. 23 Input mode for the ViscoPump parameters from Fig. 4

Pressure/Suction action: Adjustable working mode
 Ramp: This refers to the steepness of the pressure increase in terms of scale sections (1 - 50, default value = 15)
 Pump power: % of the programmed normal value
 Measurement delay: Between the individual measurements: Value between 0.5 and 2
 Above N1 Suction: Time of suction above the upper light barrier. Value between 0.5 and 2
 Pressing "Enter" will take you to **Fig. 24**.

Fig. 24 Input mode for the ViscoPump parameter; prompt following Fig. 23

This prompt is issued for data integrity reasons to avoid the automatic take-over of any erroneous input or changes.
 Yes means: accept and store data
 No means: discard changes
 Press "Enter" to return to **Fig. 4**.

Fig. 25 Input mode for the system parameters from Fig. 4

Selection of the language: Proceed to **Fig. 26**
 Selection of the RS parameters: Proceed to **Fig. 27**
 Selection of the documentation: Proceed to **Fig. 29**
 Set time and date: Proceed to **Fig. 30**
 Possibility of returning to **Fig. 4**
 Confirm, press the "Enter" key to jump to the function

Fig. 26 Input mode for the system parameters: Language

Select from:
 German, English, French, Spanish, and Italian
 Press the "Enter" key to return to **Fig. 25**

Fig. 27 Input mode for the system parameters: RS baud rate

Select the baud rate 2400, 4800, or 9600
 Used to adapt to the RS printer
 Press the "Enter" key to proceed to **Fig. 28**

Fig. 28 Input mode for the system parameters: more RS parameters

At this point you can select the word length, the stop bit, and the parity
 Used to adapt to the RS printer Press the "Enter" key to return to **Fig. 25**

System parameter

memory print out
protocol
documentation
 Yes No

System parameter

date: 21 12 03
 dd mm yy
time : 12 00 00
 hh mm ss

**System parameter
date & time
save parameters?
 Yes No**

**System parameter
delete method?
 Yes No**

**System parameter
method deleted
 press Enter**

**System report
is the capillary empty?
 press Enter**

**Error report
meas. time-out !
remove cause
 press Enter**

**Error report
safety sensor caused!
remove cause
 press Enter**

**Fig. 29 Input mode for the system parameters:
more RS parameters**

Pressing the “Enter” key will trigger an immediate memory printout, followed by the return to **Fig. 25**
Result protocol: **please refer to description of function key F6** in  4.3.1

Documentation yes / no, means that a printout is to be made after the end of a measurement series. This requires that a printer is connected; subsequently, press the “Enter” key to proceed to **Fig. 25**

**Fig. 30 Input mode for the system parameters:
date and time**

At this point you can set the internal clock. Confirm, press the “Enter” key to proceed to **Fig. 31**

**Fig. 31 Input mode for the system parameters:
confirmation**

This prompt is issued for data integrity reasons to avoid the automatic take-over of any erroneous input or changes.

Yes means: accept and store data

No means: discard changes

Press “Enter” to return to **Fig. 25**

**Fig. 32 Input mode for the system parameters:
delete confirmation**

If you confirm this question with “yes”, all system parameters will be reset to the starting valued (default values), followed by a move to **Fig. 33**; if you confirm with “no”, you will be taken to **Fig. 4**

Fig 33 System message: method deleted

If confirmed with the “Enter” key, you will be taken back to **Fig. 4**

Fig 34 System message: following a stop!

You are asked whether the capillary have run empty in order to ensure that no malfunction will be caused by bubbles or splashes in the case of a restart.

Fig. 35 Error message: measurement time-out

The cause of the time-out has to be removed. In most cases the measurement has to be restarted, since this error message is almost always caused by a fatal error, such as viscometer empty, connected incorrectly or not at all etc.

 Please refer to 2.3.8 Troubleshooting

Press the “Enter” key, and you will be taken back to **Fig. 4**

Fig. 36 Error message: safety sensor

In this case the safety bottle has to be discharged and cleaned, the cause of the overfilling has to be discovered and has to be removed.

Press the “Enter” key, and you will be taken back to **Fig. 4**

5 Maintenance and care of the device and the viscometers

Maintaining the proper functioning requires certain inspection and maintenance work.

Maintenance and service work includes:

- Visual check
- Interface function, ViscoPump III
- Once per quarter, the electrical contacts have to be inspected for corrosion, if the viscosity measuring unit is used in premises with an occasional occurrence of corrosive matters in their atmosphere.

Maintenance intervals

Normal operation	intervals for carrying out all work is 6 months
Under particular strain	the max. intervals for carrying out all maintenance work are 4 weeks
In case of disturbances	the work has to be carried out immediately

5.1 Maintenance work to be carried out

- Check the hoses, screwed connections for signs of visible damage, contamination, and leaks.
- Check the electrical plug contacts for corrosion and mechanical damage (on the device and on the cables).
- If necessary, the exterior of the casing of the viscosity measuring unit can be cleaned with a piece of cloth soaked with a household cleaning agent. The lower and rear sections have to be dry-treated.
- In no case must liquid penetrate into the interior of the lower section.
- Defective parts must be repaired or replaced with new ones.
- Defective glass parts must always be replaced.

5.2 Maintenance and care of the absorber traps VZ 7215

The VZ 7215 absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals.

 **The measurements will not work if there is an incorrect connection!** There is a risk of the sample leaking from the viscosimeter or being suctioned into the ViscoPump.

 **If sodium lime is used** acidic solvents **as an absorption agent**, the colour condition of the indicator is to be checked on a daily basis.

- As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.
- If such a colour shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a de-colouration; this may then appear as “normal”, but it will definitely result in the destruction of the pneumatic system after some time.
- This situation is explicitly excluded from the warranty coverage!

 When using **activated carbon as an absorption agent** (e.g. with solvents or used mineral oils), a replacement should be made monthly. This depends on the load factor which, in turn, is a function of the volatility of the materials.

5.3 Periods without operation

If the capillary viscometers are not used over a long period of time, the liquids contained in the system, in particular aggressive solutions, have to be drained. If the liquid is left in the system, one has to reckon that the solutions used will alter in the course of time and attack the glass, in particular the capillaries.

 **Cleaning agents should be matched to the previous samples or impurities!** In many cases, an aqueous cleaning agent (glass cleaners, detergents) or organic solvents (such as acetone or hydrocarbons) are sufficient.

⚠ Strong oxidizing cleaning agents such as chromic acid may only be used by trained personnel and must be suitably disposed for safety and environmental reasons - the current guidelines for handling hazardous materials must be observed.

⚠ In the last rinse cycle, the viscometer should be rinsed with a suitable solvent with a low boiling point (such as acetone), and dried by an air flow, which is preferably generated by underpressure (for example, water jet pump). The viscometer is dry and dust-free by this treatment and can thus be used for manual and automatic measurements.

5.4 Reproducibility of results

The measurement or analysis results depend on a variety of factors. Please check the plausibility of the measurement results or analysis results at regular intervals, and carry out the required reliability tests. In this regard, please adhere to the usual validation procedures and especially to the "Viscometers within quality assurance systems" chapter.

5.5 Viscometers within quality assurance systems

Recommendations for companies that have introduced a quality management system (OM system) according to DIN EN ISO 9001: In this quality assurance system, an inspection of the measuring equipment is planned. The intervals and required accuracy can be defined by each company according to its own requirements. The standard DIN/ISO 10 012, Part 1 serves as a guideline in this matter. We recommend regular inspection of the viscometers in defined intervals.

Inspection of the viscometer constants:

a) Calibration using comparative measurements with reference measuring standards

Comparative measurements must be performed with a viscometer (reference measuring standard) which was tested at the PTB (Federal German Physical-Technical Institute) and provided with a constant. During this comparative measurement, the viscometer to be inspected and the PTB - tested viscometer were placed simultaneously in the same thermostat bath. The test liquid tested, the viscosity of which must not be known exactly, is filled into both viscometers, tempered and the flow-through time then measured. The constants of the viscometers to be inspected are then calculated according to the following equation:

$$K = \frac{K_{PTB} \times t_{PTB}}{t}$$

K constant of the tested viscometer

K_{PTB} constant of the standard reference viscometer

t low time (HC) of the tested viscometer (corrected by Hagenbach-Couette)

t_{PTB} flow time (HC) of the standard reference viscometer (corrected by Hagenbach-Couette)

Within the quality management system in accordance with DIN EN ISO 9001, traceability of the measuring equipment to national measuring standards is demanded. This traceability can be achieved by inspecting the comparative viscometers (reference measuring standards) at regular intervals at the PTB. The time intervals are defined according to the specifications made in the quality assurance quality management system of the user.

b) Calibration of the capillary viscometer with normal oils of the PTB

During this calibration, a normal oil from the PTB with known viscosity is used as a reference measuring standard. The measurement is performed by means of flow-through measurement of the PTB normal oil in the viscometer to be inspected in a temperature bath, the temperature of which must correspond precisely to the test temperature of the PTB. In this case, it is extremely important to make sure that the temperature is absolutely correct. In case of temperature variation, this will always result in a constant for the viscometer that deviates from the constant applied. A temperature variation of 0.01 K, for instance, will result in a measuring error of up to 0.01 %. The calibration of the deviating temperature into the viscometer constant is not permitted.

c) Inspection by Xylem Analytics Germany with a quality certificate in accordance with DIN 55 350-18, 4.2.2

The inspection at the manufacturer is carried out by means of comparative measurements using viscometers as reference measuring standards that were tested at the PTB (corresponds to Item a).

i General Information on the stability of viscometer constants

Each inspection (even with a certificate) can guarantee the technical measuring direction only for a limited period of time. The constants of viscometers made of the borosilicate glass DURAN[®], however, can remain unchanged for long periods of time if the viscometers are kept away from altered influences. Especially extreme changes can be expected, for instance, during the use of liquids that attack glass, in particular hot caustic soda hydrated (NaOH) orduring glass-blowing repairs (even for apparently insignificant repairs).

Liquids whose components adhere to the glass wall also cause errors. In such cases, regular cleaning is required whereby the corrosive action cleaning agent on the glass must be eliminated. For this reason, we recommend that the user should write up a special processing instruction for all important measurements and include them in his quality management manual in accordance with DIN EN ISO 9001. In all cases the user is responsible for the correctness of his measuring and testing equipment and is not released from his responsibility for quality (cp. DIN 55 350, Part 18).

6 Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts. The breach of glass parts is also excluded. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

7 Storage and transportation

If the AVS[®] 370 has to be stored over some time, or to be dislocated, the use of the original packing will be the best protection of the devices. However, in many cases this packing will not be available anymore, so that one will have to compose an equivalent packaging system. Sealing the lower section in a foil is hereby recommended.

The devices should be stored in a room with a temperature between +10 and +40°C, and the (relative) humidity of the air should not exceed 70 %.

 If the interchangeable have to be stored over some time, or to be dislocated, the fluids inside the system, especially aggressive solution have to be removed.

8 Recycling and Disposal

Please observe the applicable local or national regulations concerning the disposal of “waste electrical and electronic equipment”.

In Germany, disposal of the device in household waste or via municipal collection points is not permitted. **Devices that are no longer used must be sent to the respective manufacturer for disposal with paid postage and the note “FOR DISPOSAL”!** The devices will then be disposed of at the manufacturer's expense.

Detailed information on “Recycling and disposal” can be found on our homepage in the customer support section under “WEEE disposal concept”. If you have any further questions about disposal, please contact the manufacturer directly (see the back of these operating manual).

The device and its packaging have been manufactured as far as possible from materials that can be disposed of in an environmentally friendly manner and recycled properly.

i This device contains batteries. Batteries must not be disposed of with household waste. The end user is obliged to dispose of used batteries (including non-hazardous batteries) at a collection point set up for this purpose or at the point of sale for proper recycling.

9 EC – Declaration of Conformity

The corresponding declaration of conformity of the device can be found on our homepage. It will also be made available to you on request.

10 Appendix viscometer type list

The following viscometers can be used in the AVS® 470 viscometer for evaluation with the Hagenbach-Couette correction:

Ubbelohde-Viskosimeter nach DIN	= type 1
Ubbelohde-Viskosimeter nach ASTM	= typ 2
Mikro-Ubbelohde-Viskosimeter	= typ 3
Cannon-Fenske Routine Viskosimeter	= typ 4
Mikro-Ostwald Viskosimeter	= typ 5

This type number is set in the menu (📖 section 4.6, fig.17).

The further differentiation is made by the entered constant or its table guide value.

11 Appendix spare parts

Spare parts can be found in our viscometry catalog.

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
Xylem Analytics Germany Sales GmbH & Co.KG
© 2022 Xylem, Inc. Version 220922 M 827 229 0