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

GA_Bezugselektroden_d_T814167_230224

Instruction manual for reference electrodes

The electrodes illustrated are example type

Xylem Analytics Germany GmbH
Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49(0)6131.894.5119
E-mail: Service-Sensors.si-analytics@xylem.com
www.XylemAnalytics.com

SI Analytics
a xylem brand

Preparation and general

The reference electrodes are delivered ready-to-use. For measuring, a cable (e.g. L1N) is needed. The irrigation cap located over the diaphragm must be removed for measurement. It contains electrolyte solution. The closure of the refill opening must be opened for measurement. The diaphragm must be immersed in the measurement solution. The electrode is ready to perform a measurement. Electrodes can be kept dry or soaked for of the different reference systems in different solutions: for Ag/AgCl reference system KCl 3.0 mol/l, for calomel reference system KCl 4.2 mol/l and for Hg/Hg₂SO₄ reference system K₂SO₄ 0.6 mol/l. Lacking electrolyte solution has to be refilled. The filling level of the electrolyte solution should always be at least 5 cm above the level of the measurement medium.

Electrodes with ground-joint diaphragm

The ground-joint diaphragm of this electrode is fixed with a shrink tubing when delivered. Remove the shrink tubing carefully with a sharp knife or scalpel for startup. To wet the ground joint surfaces, hold the electrode vertically with the refill opening open and carefully loosen the ground joint while turning it slightly until reference electrolyte solution emerges from the ground joint. Then put the ground joint sleeve back on while turning it slightly.

Measuring

Please also note the usage instructions of the measuring device during calibration and measurement. In order to minimize inaccuracies in the measurement results, electrodes used under extreme conditions or at the limits of the specified usage range should be checked more frequently. For measurement the reference electrode must be combined with a suitable indicator electrode. In combination with a pH electrode we recommend the use of our hot steam sterilized, certified buffer ampoules to DIN 19 266. Always use fresh buffer solution. In combination with a suitable metal electrode, we recommend checking with our redox test solutions, e.g. L 4648. In combination with a suitable ion-sensitive electrode, use calibration and test solutions according to the specifications in the operating manual of the electrode.

The reference electrode with double electrolyte system (B2220+ and B3920+) avoids the choice of a suitable bridge electrolyte, the interference of the measurement due to interfering ions.

For measurements in the medical and food-chemical fields, as well as in other cases where a physiological effect due to mercury is to be expected, electrodes with Ag/AgCl reference system should be used.

Storage and maintenance

Reference electrodes should be stored at a temperature of between 0° and 40 °C. Depending on storage conditions (temperature and air humidity), the irrigation fluid in the cap may dry out prematurely. In this case, the electrode must be irrigated for at least 24 hours in electrolyte solution before it is ready to use. The electrolyte for the reference electrode must be replenished or renewed occasionally. Crystals in the electrolyte area can be dissolved by warming in a water bath. The electrolyte solution should then be replaced.

Cleaning

Dirt and contamination at the membrane/Pt sensor and diaphragm will lead to measurement inaccuracies.

- **Coatings and deposits** can be removed with diluted mineral assets (e.g. hydrochloric acid 1:1), **Organic** contamination with the aid of suitable solvents, **Grease** with tenside solution and **Proteins** with hydrochloric pepsin solution (cleaning solution L 510).
- After cleaning, rinse off the electrodes with distilled water.
- Ceramic diaphragms blocked from outside can be restored to working order by carefully rubbing down with fine sandpaper or a diamond file.
- Platinum diaphragms must not be subjected to any mechanical treatment. Rinsing (e.g. vacuuming) may be followed by chemical cleaning (e.g. with diluted hydrochloric acid).
- Ground diaphragms are prepared for measurement by gently raising and then replacing the ground sleeve on the core. The refilling hole must be open during this process. **Caution:** More electrolyte flows out during this process to ensure full coverage of the ground surface.

Quality

Every electrode must meet the strict quality requirements of final testing. The durability depends mainly on the usage conditions. Extreme conditions include for example high or frequently fluctuating temperatures, strong acids and caustic solutions, protein and heavily contaminated solutions as sulphides, bromides and iodides. Hydrofluoric acid and hot phosphoric acid corrode glass.

Further information

Damaged sensors with calomel and Hg/Hg₂SO₄ reference system can deliver toxic mercury and mercury compounds in the environment. These sensors must be disposed in accordance with your local requirements.

Further information can be found in our meters and electrodes catalog.

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