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

Instruction manual for SI Analytics® Puncture knife Z 503

Xylem Analytics Germany GmbH

Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49(0)6131.66.5119
Fax: +49(0)6131.66.5001
Email: Service-Sensors.si-analytics@xylem.com
www.XylemAnalytics.com

SI Analytics
a xylem brand

Safety notes

! Read the safety notes and the operating manual of the used electrodes carefully in order to ensure a faultless use.

! Risk of injury!

- When removing the protective cap ensure that there is no hand contact with the blade of the puncture aid.
- Ensure during the assembly/disassembly of the puncture aid that there are no injuries through the blade of the puncture aid. Possibly wear cut protection gloves!
- Always put on the protective cap when not using the puncture blade!

After use of the puncture aid, always clean it thoroughly (see section "Cleaning")

Use as intended

The puncture knife together with the pH electrodes of the structural shape BlueLine 21pH serves for the direct pH puncture measurement in food (e.g. meat or cheese).

Specifications

Material:


Blade: Knife steel DIN 1.4112 and 57-59 HRC

Handle piece, protective cap, spherical cap, handle screws

Polybutylene terephthalate (PBT, white, food-compliant, FDA-compliant)

i All materials are suitable for cleaning in the dishwasher up to a temperature of 60 °C.

Handle screw versions

Commissioning

- Remove the watering cap from the electrode
- Possibly connect the cables to the electrode.
- Push the electrode through the handle piece.
- Introduce the electrode into the puncture knife.

i Observe the position of the diaphragm opening. **Orientation of the diaphragm parallel to the blade!**

- Screw the handle piece onto the puncture blade:
- Screw the cable through the slot of the of the appropriate handle screw (see section "Handle screw versions") and screw it into the handle piece.
- Guide the cable laterally out through the slot and screw the spherical cap on tightly.
- Connect the cable to the pH meter.

Measurement of the pH value

i For measuring the pH value and for calibration please observe the respective operating manual of the used pH electrode and of the pH meter.

- Remove the protective cap prior to the pH measurement.
- Plunge the puncture blade 35-50mm deep into the measurement good and read the pH value as soon as it is stable.

Storage of the installed pH electrode in the puncture aid

The storage of the installed pH electrode should only take place for short intervals in the protective cap of the puncture knife. When storing it for several hours or longer, please dismantle the electrode from the puncture aid, clean it an insert it into the electrode watering cap (filled with 3mol/l KCL).

Cleaning

The individual parts of the puncture aid can be cleaned with a customary dish washing detergent and a brush, afterwards with ethanol. Cleaning in a dishwasher is also possible.

Further information

You can find further information in our "Measuring devices and electrode catalog" or on our homepage. Technical changes reserved.

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