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
Instruction manual for generator electrode with and without diaphragm for coulometric measurement

Xylem Analytics Germany GmbH

Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49.(0)6131.66.5119
Fax: +49.(0)6131.66.5001
E-mail: Service-Sensors.si-analytics@xylem.com
www.XylemAnalytics.com

SI Analytics
a xylem brand

Preparation and general

The generator electrode with diaphragm must be dried before the first measurement and after every cleaning in the drying oven at **max. 50°C** for 30 min.

The generator electrode without diaphragm is ready to perform the measurement.

Storage and maintenance

The storage conditions for generator electrodes are uncritical. They are in the range of 0–50°C and 5–95% relative humidity. If often used, a storage in the titration vessel is recommended.

Change of catholyt: depending on application, and titration volume.

Change of anolyt: depending on sample consistence and water content.

Cleaning

Discoloration or precipitations on **anode** and **cathode** can be removed with conc. nitric acid. Oil contamination can be removed with alkaline tensid solutions or solvents.

Precipitations on the **diaphragm** will be softened some hours with conc. nitric acid. Therefore the electrode is fixed with e.g. a stand, placed in a beaker and filled with conc. nitric acid. Afterwards rinse it extensively with dist. water.

The cathode chamber is filled with ethanol and the solution should be spilled out through the diaphragm. This rinsing procedure should be done 2-3 times after every cleaning.

Please wear suitable protective equipment and consider the warning and safety instructions of the corresponding solutions.


Quality

Every electrode must meet the strict quality requirements of final testing. The durability depends mainly on the usage conditions.

Subject to technical amendment.

Application Hints for Coulometric KF Titration with Generator Electrode TZ 1753

In case of using a generator electrode with glass frit, different reagents have to be used: In the anode cell the anolyte for Iodine generation, in the cathode cell the catholyte for the related redox reaction (pic 1).

Pic. 1 Addition of different reagents in case of using the generator electrode with glass frit.

Fill in about 100 ml anolyte and up to 5 ml catholyte. Add as much catholyte as it is necessary to have one liquid level of catholyte and anolyte. Typically 5 ml ampouls are available. For different applications anolyte and catholyte can be adjusted to the sample matrix. In the case the sample contains Aldehydes or Ketons, special „K“ Reagents are used as anolyte as well as catholyte.


For optimizing the solubility of unpolar samples, unpolar solvents are added or special unpolar anolytes are used. In such cases it is possible that the total conductivity is reduced in the way that the electricity yield is not any longer complete and the result is higher than expected.

This happens often in cases, where Chloroform content exceeds a certain amount. In such cases we recommend to use ecological more harmless, modern Chloroform free reagents. In cases this is applicable, an correction factor can be used in the calculation or additional 5 ml catholyte in the anode cell improves the electric yield.


What is happening in detail?

In the coulometry the titration reagent is produced electrochemically out of the anolyte at the anode of the generation electrode. Iodide is oxidized to Iodine.

The indicator electrode indicates, whether the reaction is ready or not. As soon as there is a little amount of Iodine in surplus, a current is measured by the reversible redox pair. (pic.2). In the generator electrode the forming of Hydrogene takes place as the related redox reaction. The cathode cell is separated from the anode cell by a glass frit. (Pic. 3).

Pic. 2 Principle of the Coulometric titration

Pic. 3 The generator electrode