Lab 955

Contents

EN English	3
DE Deutsch.....	7
CS Český	11
DA Dansk	13
EL Ελληνικά	15
ES Español.....	17
FI Suomi	21
FR Français	23
HU Magyar	27
IT Italiano	29
JA 日本語	33
KO 한국어.....	35
NL Nederlands	37
NO Norsk	39
PL Polski.....	41
PT Português	43
RO Română	47
RU Русский	49
SK Slovensky	51
SV Svenska	53
TR Türkçe	55
ZH 中文.....	57

Copyright

© 2016 Xylem Analytics Germany GmbH
Printed in Germany.

Conductivity meter, Lab 955

Overview


The Lab 955 compact precision conductivity meter enables you to perform conductivity measurements quickly and reliably.

The Lab 955 handheld meter provides the maximum degree of operating comfort, reliability and measuring certainty for all applications.

Keypad


Keys are indicated by brackets <.>.

The key symbol (e.g. <ENTER>) generally indicates a short keystroke (under 2 sec) in this operating manual. A long keystroke (approx. 2 sec) is indicated by the underscore behind the key symbol (e.g. <ENTER_>).


<On/Off>:	Switches the meter on/off
<On/Off_>:	Resets calibration data
<M>:	Selects the measured parameter
<M_>:	Opens the setting menu for calibration and measurements
<CAL>:	Calls up the calibration procedure
<CAL_>:	Displays the calibration data
<▲>:	Increments values, scrolls
<▼>:	Decrements values, scrolls
<ENTER>:	Confirms entries
<ENTER_>:	Opens the setting menu for system settings

Display

Status display indicators	AR	Stability control (AutoRead) is active
	ARng	Automatic range switchover; meter measures with highest possible resolution
	Cal	Calibration
	LoBat	With battery operation: Batteries almost empty
	nLF	Nonlinear temperature compensation
	TP	Temperature measurement active
	Tref20	Reference temperature of 20 °C
	Tref25	Reference temperature of 25 °C
	TIME	Setting of calibration interval

Socket field connections

1	Conductivity measuring cell
2	Service interface
3	Connection for line adaptor

Initial commissioning

1	Insert the supplied batteries. When doing so make sure that the poles of the batteries are positioned correctly.
2	Connect the power pack if necessary.
3	Switch on the meter with the <On/Off> key.

Caution

Follow the safety instructions of the used sensors.

The operating manuals of the sensors are available on the supplied CD and on the Internet under www.WTW.com.

Determining/ setting up the cell constant [C]

Due to aging, the cell constant slightly changes. As a result, an inexact measured value is displayed. Calibration determines the current value of the cell constant and stores this value in the meter.

Thus, you should calibrate at regular intervals.

You can either determine the cell constant of the conductivity measuring cell in the range $0.800 \dots 0.880 \text{ cm}^{-1}$ by calibration in the control standard 0.01 mol/l KCl or adjust it manually in the range $0.800 \dots 1.200 \text{ cm}^{-1}$.

Determining the cell constant (calibration in control standard)

1	Press <CAL> repeatedly until <i>CAL CELL</i> is displayed.
2	Press <ENTER> or <CAL__> to confirm the selection of <i>CAL CELL</i> . The cell constant of the last calibration is displayed.
3	Immerse the measuring cell in the control standard solution, 0.01 mol/l KCl .
4	Start the calibration with <ENTER> . The determination of the cell constant with stability control starts. The [Auto] display flashes until a stable signal is available. The cell constant determined is displayed. The meter automatically stores the cell constant.
5	Switch to the measuring mode with <ENTER> . The determined cell constant is used.

Using the last calibrated cell constant

Prerequisite: No valid calibration exists.

1	Press <CAL> repeatedly until <i>USE CELL</i> is displayed.
2	Press <ENTER> or <CAL__> to confirm the selection of <i>USE CELL</i> .
3	If necessary, press <CAL> repeatedly until <i>CAL</i> and the last calibrated cell constant is displayed.
4	Confirm the selection with <ENTER> . The displayed cell constant is used. The meter switches to the measured value display.

Setting the cell constant manually

The cell constant to be set must either be taken from the operating manual of the measuring cell or is printed on the measuring cell.

1	Press <CAL> repeatedly until <i>USE CELL</i> is displayed.
2	Confirm the selection with <ENTER> or <CAL__> . The cell constant that was set last is displayed.
3	If necessary, press <CAL> repeatedly until a cell constant in the range 0.800 ... 1.200 cm ⁻¹ is displayed.
4	Set the cell constant to be used with <▲><▼> , e.g. 0.846 cm ⁻¹ .
5	Confirm the selection with <ENTER> . The new cell constant is used from now on. The meter switches to the measured value display.

Measuring the conductivity

6	Immerse the conductivity measuring cell in the test sample.
7	If necessary, scroll with <M> until the measured parameter \mathcal{X} with the unit mS/cm or μ S/cm is displayed.
8	Wait for a stable measured value. The [Auto] display indicator flashes as long as the measured value is not yet stable.

Calibration interval (*Int.C*)

The calibration interval reminds you to calibrate regularly. After the specified calibration interval (*Int.C*) has expired, the symbol for the calibration evaluation flashes. It is still possible to measure.

1	Open the menu for measurement settings with <M__> .
2	Confirm all settings with <ENTER> until <i>Int.C</i> is displayed.
3	Set the calibration interval with <▲><▼> .
4	Confirm the setting with <ENTER> .

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Xylem |'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water technology company.

We're a global team unified in a common purpose: creating advanced technology solutions to the world's water challenges. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. Our products and services move, treat, analyze, monitor and return water to the environment, in public utility, industrial, residential and commercial building services settings. Xylem also provides a leading portfolio of smart metering, network technologies and advanced analytics solutions for water, electric and gas utilities. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise with a strong focus on developing comprehensive, sustainable solutions.

For more information on how Xylem can help you, go to www.xylem.com.

SI Analytics

a xylem brand

Hersteller

(Manufacturer)

Xylem Analytics Germany GmbH

Am Achalaich 11
82362 Weilheim
Germany

SI Analytics

Tel. +49.(0)6131.66.5111

Fax. +49.(0)6131.66.5001

E-Mail: si-analytics@xylem.com

www.xylemanalytics.com

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(Service and Returns)

Xylem Analytics Germany Sales GmbH & Co. KG

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Fax. +49.(0)881.183.414

E-Mail: wtw.rma@xylem.com

SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
© 2021 Xylem, Inc.