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

Instruction manual for electrodes for measurement of pH values and redox potential

Also with integrated temperature sensor

The electrodes illustrated are example types

Xylem Analytics Germany GmbH
Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49(0)6131.894.5119
E-mail: Service-Sensors.si-analytics@xylem.com
www.XylemAnalytics.com

SI Analytics
a xylem brand

Preparation and general

The pH and redox combination electrodes are delivered ready to use. The irrigation cap located over membrane/sensor must be removed for measurement and calibration. It contains irrigation fluid. The electrode is ready to perform a measurement. Electrodes stored dry must be irrigated for 24 hours in irrigation fluid. Lacking potassium chloride solution $c(\text{KCl}) = 3 \text{ mol/l}$ in the electrolyte area of the reference system is refilled. The filling level of the electrolyte solution should always be at least 5 cm above the level of the measurement medium. The closure of the refill opening must be opened for calibration and measurement. The diaphragm must be immersed in the measurement solution. Refilling is not required for low-maintenance electrodes with gel filling or the REFERID® system. Irrigation with electrolyte solution is however particularly important for these electrodes.

Electrodes with ground-joint diaphragm

The ground-joint diaphragm of this electrode is fixed with a shrink tubing when delivered. Remove the shrink tubing carefully with a sharp knife or scalpel for startup. To wet the ground joint surfaces, hold the electrode vertically with the refill opening open and carefully loosen the ground joint while turning it slightly until reference electrolyte solution emerges from the ground joint. Then put the ground joint sleeve back on while turning it slightly.

Measuring the pH value

Please also note the usage instructions of the measuring device during calibration and measurement. In order to minimize inaccuracies in the measurement results, electrodes used under extreme conditions or at the limits of the specified usage range should be calibrated more frequently. For accurate calibration, we recommend the use of our hot steam sterilized, certified buffer ampoules to DIN 19 266. Always use fresh buffer solution.

Measuring the Redox voltage

For measurement purposes, the pH-meter is switched to the mV-range. In the case of metal single-rod measurement chains, an Ag/AgCl reference system is used as standard. Immerse the electrode measurement chain in the measurement solution and read off the measurement value immediately, without performing a calibration. Redox test solutions are available to check the results. The correct temperature of the reference system should be observed.

Storage and maintenance

Electrodes should be stored at a temperature of between 0° and 40 °C. Depending on storage conditions (temperature and air humidity), the irrigation fluid in the cap may dry out prematurely. In this case, the electrode must be irrigated for at least 24 hours in potassium chloride solution $c(\text{KCl}) = 3 \text{ mol/l}$ before it is ready for use.

In the case of pH single-rod measurement chains/Redox measurement chains and reference electrodes, the electrolyte must be topped up or replaced occasionally. Crystals in the electrolyte area can be dissolved by warming in a water bath. The electrolyte solution should then be replaced, and the correct potassium chloride concentration checked.

Cleaning

Dirt and contamination at the membrane/Pt sensor and diaphragm will lead to measurement inaccuracies.

- **Coatings and deposits** can be removed with diluted mineral assets (e.g. hydrochloric acid 1:1), **Organic** contamination with the aid of suitable solvents, **Grease** with tenside solution and **Proteins** with hydrochloric pepsin solution (cleaning solution L 510).
- After cleaning, rinse off the electrodes with distilled water, do not rub dry.
- Ceramic diaphragms blocked from outside can be restored to working order by carefully rubbing down with fine sandpaper or a diamond file. **Take care not to scratch the pH glass membrane!**
- Platinum diaphragms must not be subjected to any mechanical treatment. Rinsing (e.g. vacuuming) may be followed by chemical cleaning (e.g. with diluted hydrochloric acid).
- Ground diaphragms are prepared for measurement by gently raising and then replacing the ground sleeve on the core. The refilling hole must be open during this process. **Caution:** More electrolyte flows out during this process to ensure full coverage of the ground surface.
- The glass membrane can be cleaned by rubbing off with a fluff-free cloth soaked in ethanol.

Quality

Every electrode must meet the strict quality requirements of final testing. The durability depends mainly on the usage conditions. Extreme conditions include for example high or frequently fluctuating temperatures, strong acids and caustic solutions, protein and heavily contaminated solutions as sulphides, bromides and iodides. Hydrofluoric acid and hot phosphoric acid corrode glass.

Further information

Further information can be found in our meters and electrodes catalog.

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