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Short manual for SI Analytics®

Glass Capillary Viscometer

Thank you for choosing a product of Xylem Analytics Germany. This short manual contains the most important information for handling and cleaning the viscometers. Please observe the detailed operating manual of the individual viscometer types. You will find them on our homepage ([see QR Code](#)).

General and Preparation

SI Analytics® glass capillary viscometers are made of DURAN® borosilicate glass 3.3, the TC viscometers are also equipped with thermoelectric sensors – the different ways of cleaning the viscometers must be observed. The viscometers are delivered in ready-to-operate condition. In the case of calibrated viscometers, residues of the calibration or cleaning fluid may still be present despite factory cleaning. Before the first use, it is recommended to clean the viscometers with e.g. petroleum benzine and, if necessary, acetone if aqueous samples are to be measured afterwards.

Use/Application

Glass capillary viscometers are designed for use in thermostatic baths (e.g. type CT 72). Suitable holders or brackets must be used for this purpose. In the case of viscometers of the Ubbelohde type, brackets are a good protection against damage. The viscometers should preferably remain in the brackets, also for handling outside the bath, e.g. during cleaning and storage.

Filtration of samples

Viscometers with narrow capillaries are particularly sensitive to contamination: the flow times scatter when particles are present in the sample solutions. To remedy this situation, samples with particles must always be filtered in a suitable manner before being filled into the viscometer.

Cleaning of viscometer for optical detection – for the visual or automatic flow time measurement by means using IR light barriers

These viscometers have a particularly high chemical resistance.

Standard cleaning: After measuring aqueous solutions, water is suitable as a cleaning agent, possibly with the addition of surfactant. After measuring non-aqueous samples, rinse with suitable solvent. Organic solvents do not attack the glass.

Attention

Laboratory cleaners are alkaline and can attack glass. Therefore, use only diluted, at room temperature and not over long periods of time.

For **intensive cleaning**, all organic solvents or also acids such as conc. sulfuric acid can be used, depending on the contamination. Even oxidizing additives such as chromosulfuric acid, NoChromix® or piranha solution do not attack the viscometer.

Attention

For all aggressive cleaning fluids, the respective guidelines and safety precautions must be observed. Not suitable as cleaning agents due to glass attack are alkaline solutions, hydrofluoric acid and hot phosphoric acid.

Drying: Viscometers for optical detection can be dried in a drying cabinet at temperatures up to 150°C. Alternatively: Drying by rinsing with a highly volatile solvent (e.g. acetone) and drying in an air stream, preferably generated by negative pressure (e.g. water jet pump).

Cleaning of TC Viscometer – for automatic detection by thermoelectric sensors

These viscometers have built-in thermoelectric sensors and therefore have a limited chemical resistance.

These viscometers can be cleaned with hydrocarbons (e.g. benzene, toluene) and alcohols as well as aqueous surfactant solutions. Stronger organic solvents are also permitted for short periods.

Attention

Aggressive cleaning agents such as concentrated acids, especially with oxidizing additives, are not permitted. Exposure to strong solvents such as ketones (acetone) must be limited to a short time - do not leave standing for hours or overnight. Alkaline cleaning agents are also not permitted.

Drying: Drying by rinsing with a highly volatile solvent (e.g. acetone) and drying in an air stream, preferably generated by negative pressure (e.g. water jet pump).

Subject to technical changes.

Xylem Analytics Germany GmbH

Am Achalaich 11
82362 Weilheim
Germany

SI Analytics
Tel.: +49.(0)6131.66.5119
Fax: +49.(0)6131.66.5001
E-mail: Service-Instruments.si-analytics@xylem.com
www.XylemAnalytics.com

SI Analytics
a xylem brand