**BEDIENUNGSANLEITUNG**  
**OPERATING MANUAL**  
**MODE D'EMPLOI**  
**INSTRUCCIONES DE OPERACIÓN**

ba15352d04 12/2021

---

**Ca 800 DIN**

---

a xylem brand

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
a xylem brand

## View

1	Connection head with connection cable
2	Shaft
3	Exchange membrane cap
4	PVC membrane
5	O-ring
6	Inner junctions
7	Exterior junctions
8	Storage vessel for Ca 800


## Commissioning

### Combination electrode Ca 800

- 1 Remove the storage vessel.
- 2 Fill ELY/BR/503 bridge electrolyte into the filling opening so that the inner junctions are covered with bridge electrolyte.
- 3 Rinse the combination electrode with deionized water.
- 4 Wipe the shaft using a clean paper towel.

## Conditioning, calibration, measurement

### General information

When operating the electrode ensure that

- the filling opening (1) for the bridge electrolyte is open
- the inner junctions (3) are covered with bridge electrolyte
- no air bubbles are in the bridge electrolyte
- the depth of immersion is within the optimum range:

Minimum depth of immersion

The outer junction (4) must be covered

Maximum depth of immersion

Approx. 1 cm below the fluid level (2) of the bridge electrolyte

### Before measuring

- 1 Condition every ion sensitive electrode in standard solution before using it.

The conditioning time depends on the previous history of the ion sensitive electrode. To achieve a good responding behavior, the concentration of the standard solution should be similar to the concentration expected of the test sample. During commissioning and after long periods of storage: condition the electrode overnight.

Connect the electrode(s) to the meter and switch to the measured value display (see operating manual of the meter). The ion sensitive electrode is sufficiently conditioned when the measured value no longer changes considerably under otherwise constant conditions (temperature, composition of the solution).

- 2 Remove any air bubbles in the bridge electrolyte by slightly knocking against the shaft.
- 3 Calibrate according to the operating manual of the meter and the analysis specification.

## Sample preparation

Add 2 % ISA/Ca solution.

This sample conditioning solution creates optimum conditions for measuring. It provides a constant ionic strength and similar diffusion potentials at the reference electrode in standard solution and test sample.

If you would like to have more detailed information concerning sample preparation and measuring procedures, WTW provides a large number of application reports for various applications.

## Response times

The response time depends on the concentration range. It is

- several seconds at high concentrations,
- several minutes near the detection limit.

The measured value is stable if the value does not change by more than 0.1 mV within 30 seconds.

## Interferences

Interfering ions: 10 % error with the following concentration ratio (concentration ratio = interfering ion / measured ion):

Pb <sup>2+</sup>	Hg <sup>2+</sup>	H <sup>+</sup>	Sr <sup>2+</sup>	Fe <sup>2+</sup>	Cu <sup>2+</sup>	Ni <sup>2+</sup>	NH <sub>4</sub> <sup>+</sup>
0.01	4	4	6	20	40	50	200

Na <sup>+</sup>	Tris <sup>+</sup>	Li <sup>+</sup>	K <sup>+</sup>	Ba <sup>2+</sup>	Zn <sup>2+</sup>	Mg <sup>2+</sup>	
200	300	300	400	700	1000	1000	

## Aging

Please note that every (combination) electrode undergoes a natural aging process. The response time increases and the slope decreases with the age of the (combination) electrode. The following factors shorten the lifetime considerably:

- Incorrect storage
- Special measuring conditions (e.g. organic solutions, frequent measurement with high concentrations of interfering ions)
- High temperatures
- High changes in temperature

The warranty does not cover damage caused by measuring conditions and mechanical damage.

## Maintenance

- Combination electrodes:  
Refill any used up bridge electrolyte.
- Install an exchange membrane cap.

### NOTE

**Never pull the connection cable of the electrode.  
The cable might be damaged.**

## Storage

### Between two measurements

Condition the electrode in standard solution.

### Overnight

Rinse the combination electrode with deionized water and dry it. Then put on the storage vessel and store the electrode upright.

### Longer than overnight

Remove the bridge electrolyte, then rinse the combination electrode with deionized water and dab it dry with a clean paper towel. Store the combination electrode in the storage vessel in a dry place. For the following measurement, prepare the combination electrode as described in section COMMISSIONING.


Store the reference electrode according to the instructions in its operating manual.

## Recommended accessories

Description	Model	Order no.
Exchange membrane cap for calcium combination electrode Ca 800	Ca 800/AT	106656
Bridge electrolyte	ELY/BR/503	106 575
ISA sample conditioning solution for Ca <sup>2+</sup> measurement	ISA/Ca	140120
Standard solution 10 g/l Ca <sup>2+</sup>	ES/Ca	120200

## Calibration line of a calcium ion sensitive electrode

## What to do if ...

Measured value unstable	Cause	Remedy
	– Inner junctions not sufficiently wetted with bridge electrolyte (Ca 800)	– Fill up bridge electrolyte until the inner junctions are covered with bridge electrolyte
	– Junctions encrusted	– Leave the bridge electrolyte to react on the junctions for some hours until the crusts have dissolved.
	– Cable broken	– Exchange the electrode

Slope too low	Cause	Remedy
	– Conditioning time too short	– Extend the conditioning time
	– Standard solutions too old	– Use new standard solutions
	– Junctions encrusted	– Leave the bridge electrolyte to react on the junctions for some hours until the crusts have dissolved.
	– Electrode defective	– Exchange the electrode

## Technical data

<b>Measuring range</b>	0.02 ... 40,000 mg/l Ca <sup>2+</sup>	(5 x 10 <sup>-7</sup> ... 1 mol/L Ca <sup>2+</sup> )
<b>Reproducibility</b>	± 4 %	
<b>pH range</b>	2.5 ... 11	(see INTERFERENCES)
<b>Temperature range</b>	0 ... 40 °C	
<b>Membrane resistance</b>	1 to 4 MΩ	
<b>Length</b>	170 mm (including 50 mm connection head)	
<b>Diameter</b>	Shaft:	12 mm
	Connection head	16 mm

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**Enchufe** Enchufe DIN

# Xylem |'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water technology company.

We're a global team unified in a common purpose: creating advanced technology solutions to the world's water challenges. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. Our products and services move, treat, analyze, monitor and return water to the environment, in public utility, industrial, residential and commercial building services settings. Xylem also provides a leading portfolio of smart metering, network technologies and advanced analytics solutions for water, electric and gas utilities. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise with a strong focus on developing comprehensive, sustainable solutions.

**For more information on how Xylem can help you, go to [www.xylem.com](http://www.xylem.com).**

## **Service and Returns:**

Xylem Analytics Germany  
Sales GmbH & Co. KG  
WTW  
Am Achalaich 11  
82362 Weilheim  
Germany

Tel.: +49 881 183-325  
Fax: +49 881 183-414  
E-Mail [wtw.rma@xylem.com](mailto:wtw.rma@xylem.com)  
Internet: [www.xylemanalytics.com](http://www.xylemanalytics.com)

Xylem Analytics Germany GmbH  
Am Achalaich 11  
82362 Weilheim  
Germany

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**Copyright** © 2021, Xylem Analytics Germany GmbH