

Water determination according to Karl Fischer in organic liquids

Date: 09.08.2013 Page 1 of 8



Use

This method is a general description how to determine the water content in organic liquids such as methanol, isopropyl alcohol, glycol and other similar liquids.

#### **Appliances**

- Titrator: TL 7500 KF (TL 7500 KF 05/10/20) or TL 7750 with KF accessories
- µl syringe or one-way insulin syringe with needle (> 0.8 mm)
- 4 digit balance

#### **Electrodes**

KF 1100 with fixed cable

#### Reagents

Titration agent: 1 component titrant or 2 component Titrant with 2 or 5 mg/ml concentration

Solvent: dry methanol, combi solvent or a two component solvent/solvent E

#### **Description**

The setup of the instrument is described in the operating manual.

#### Determination of the titrant concentration

The determination of the titrant concentration (titer determination, standardization) is necessary before the sample titration and we recommend a weekly titer determination.

Fill the titrant in the exchange unit and the solvent in the titration vessel as described in the manual. Do not forget to enter the concentration of the titrant (e.g. 2 or 5) into the WA unit. Otherwise the drift determination is wrong. Load the standard method (default method) Titer 1-Comp. (liq. std.) or Titer 2-Comp. (liq. std.) as user method and start the method.

The KF titrator starts automatically the conditioning routine. The solvent and also the titration vessel contents also water which has to be titrated before. When the conditioning is ready you see "Conditioning ready" on the screen.

The conditioning remains active until the actual titration is started by pressing <F1/START>.

- Open the glass ampoule with the water standard (10 mg/g) and rinse the 10 ml syringe 1-2 times with 0.5 1 ml of the standard.
- Fill the complete content of the liquid standard then into the syringe and wipe then the needle with a tissue.
- Place the syringe with needle on a 4 digit balance. Maybe use a beaker, place this beaker on the balance and put the syringe then into the beaker during the weighing process.
- Press the tara button.
- Press Start/F1 to interrupt the conditioning. You are prompted immediately to add the sample.
- Please add between 0.5 and 1 g standard for a 2 mg/ml titrant and 1.5 2 g for a 5 mg/ml titrant through the septum,
- Place the syringe again on the balance and read the absolute weight form the balance (ignore the minus)
- Enter the sample ID (optional)
- Enter the sample weight and confirm with ENTER/F1

Date: 09.08.2013 Page 2 of 8


## **Application**

a xylem brand

The titrator starts then automatically the titration. The result will appear on the screen at the end of the titration. Repeat the titer determination two times (total 3 times). The average value of the titrant concentration in mg/ml will be stored automatically in the exchange unit and used for the calculation of the sample water concentration.

## GLP documentation





#### Method data

Method name:	Titer 2-Comp. (liq. st.)	Titration duration:	1 m 43 s
End date:	27.05.13	End time:	16:03:32

#### Titration data

Sample ID:	10.0 merck	Weight:	0.90070 g
Start µA:	-0.052 μA	End µA:	22.873 μΑ
Start drift:	59 μg/min	End drift:	43 µg/min
Consumption:	4.624 ml	Drift correction:	Off
Tibour	1.0470		

Titer: 1.9479 mg/ml

Mean value: --- RSD: ---

## Calculation formula

riter:	(W*F2)/((EP-B)*M*F1) -> WA	MOI (M):	1.00000

 Weight (W):
 0.90070 g (m)
 Factor 2 (F2):
 10.0000

 Blank value (B):
 0.0000 ml
 Factor 1 (F1):
 1.0000

Statistics: 1 from 3

Date: 09.08.2013 Page 3 of 8

# **SI Analytics**

05/27/13 13:26:56

05/27/13 15:57:51

GLP

Created at:

Last modification:

Documentation:

a xylem brand

### Method data overall view

Method name: Titer 2-Comp. (liq. st.)
Method type: Automatic titration

Titration mode: KF Linear steps: 0.004 ml Fixed delay time: 0 s

 $\begin{array}{lll} \text{Extraction time:} & 0 \text{ s} \\ \text{Pretitration:} & \text{Off} \\ \text{Endpoint:} & 20.0 \ \mu\text{A} \\ \text{delta endpoint:} & 15.0 \ \mu\text{A} \\ \end{array}$ 

Endpoint delay: 10 s
Drift: 50 µg/min
Drift correction: Off

Min. titration time: 60 s Max. titration time: 600 s

Polarization voltage: 100 mV

## Dosing parameter

Dosing speed: 60.00 % Filling speed: 30 s Maximum dosing volume: 50.00 ml

## Unit values

 Unit size:
 10ml

 Unit ID:
 10035454

 Reagent:
 Composite 2

 Batch ID:
 no entry

 Concentration [mg/ml]:
 1.94360

Determined at: 05/14/13 23:46:49

Expire date: -Opened/compounded: -Test according ISO 8655: 08/21/12

Last modification: 05/27/13 15:56:13

Date: 09.08.2013 Page 4 of 8

# SI Analytics

## **Application**

a xylem brand

#### **Sample Titration**

The procedure for organic liquid samples is the sample as for the liquid standard. The only difference is the used method. Please load the default method "Sample 1 Comp." or "Sample 2 Comp." for the sample titration. Use also a suitable syringe (1, 2, 5 or 10 ml) with needle for the liquid sample.

- Rinse syringe 1-2 times with 0.5 1 ml of the sample. Use between 0.2 10 g of the sample depending on the expected water concentration (Please refer to the enclosed table)
- Place the syringe with needle on a 4 digit balance. Maybe use a beaker, place this beaker on the balance and put the syringe then into the beaker during the weighing process.
- Press the tara button.
- Press Start/F1 to interrupt the conditioning. You are prompted immediately to add the sample.
- · Add the sample trough the septum
- Place the syringe again on the balance and read the absolute weight form the balance (ignore the minus)
- Enter the sample ID (optional)
- Enter the sample weight and confirm with ENTER/F1

The titrator starts then automatically the titration. The result will appear on the screen at the end of the titration.

## Recommended sample amounts:

expected water amount	Titre 2 mg/ml	Titre 5 mg/ml
10 ppm – 100 ppm	> 10 g	
100 – 500 ppm	2- 10 g	
500 ppm – 1000 ppm (0.1 %)	1- 5g	5 – 10 g
0.1 % - 10 %	0.5 – 1 g	1 – 2 g
1 % - 10 %		0.1 – 2 g
10 % - 100 %		max. 0.1 g

Date: 09.08.2013 Page 5 of 8

## **Application**


a xylem brand

#### Methods

## Sample titration (page 1)

## GLP documentation

### Titration graph



## Method data

Method name: CombiSolvent 1 Titration duration: 2 m 1 s
End date: 06.09.12 End time: 16:34:09

#### Titration data

Ethanol Weight: Sample ID: 0.65630 g End µA: Start µA: 20.565 μΑ 23.324 μΑ Start drift: 45 µg/min End drift: 99 µg/min Consumption: 3.015 ml Drift correction: Off

result: 2.379 %

### Calculation formula

result: (EP-B)\*T\*M\*F1/(W\*F2)

Mol (M): 1.00000

Blank value (B): 0.0000 ml Titre (T): 5.17870000 (a)

 Factor 1 (F1):
 0.1000
 Weight (W):
 man

 Factor 2 (F2):
 1.0000
 Statistics:
 Off

Date: 09.08.2013 Page 6 of 8

## a xylem brand

### sample titration (page 2):

## Method data overall view

Method name: CombiSolvent 1 Created at: 09/06/12 16:25:47

Method type: Automatic titration Last modification: 09/06/12 16:32:18

Documentation:

Titration mode: KF
Linear steps: 0.005 ml
Fixed delay time: 1 s

 $\begin{array}{lll} \text{Extraction time:} & 10 \text{ s} \\ \text{Pretitration:} & \text{Off} \\ \text{Endpoint:} & 20.0 \ \mu\text{A} \\ \text{delta endpoint:} & 4.0 \ \mu\text{A} \\ \end{array}$ 

Endpoint delay: 10 s Drift: 100 µg/min

Drift correction: Off

Min. titration time: 30 s Max. titration time: 600 s

Polarization voltage: 100 mV

## Dosing parameter

Dosing speed: 60 % Filling speed: 30 s

Maximum dosing volume: 50.00 ml

Unit values

Unit size: 10ml Unit ID: 10035463

Reagent:


Batch ID: no entry
Concentration [mg/ml]: 5.17870
Determined at: -Expire date: -Opened/compounded: --

Test according ISO 8655: 08/27/12

Last modification: 09/06/12 16:15:24

Date: 09.08.2013 Page 7 of 8

## **Application**



a **xylem** brand

### **Hints**

If you have any questions concerning the application, you are welcome to contact us.

Literature

SI Analytics GmbH Hattenbergstr. 10 55122 Mainz Germany Phone: +49 (0) 6131 / 66 – 5062

+49 (0) 6131 / 66 – 5118 +49 (0) 6131 / 66 – 5001

Fax: +49 (0) 6131 / 66 – 5001 E-Mail: titration@si-analytics.com Homepage: www.si-analytics.com

Date: 09.08.2013 Page 8 of 8