

CR 2200

THERMOREACTOR



Contents

CR 2200 - Contents

1	Ove	Overview5				
	1.1	Components of the thermoreactor	6			
2	Safe	ety	7			
	2.1	Safety information	7			
		2.1.1 Safety information in this operating manual	7			
		2.1.2 Safety signs on the product				
		2.1.3 Further documents providing safety information				
	2.2	Safe operation				
		2.2.1 Authorized use				
		2.2.2 Requirements for safe operation				
	0.0	2.2.3 Unauthorized use				
	2.3	Personal protective equipment	٠. و			
3	Con	nmissioning	. 10			
	3.1	General handling instructions	. 10			
	3.2	Scope of delivery	. 10			
	3.3	Unpacking and positioning				
	3.4	Initial commissioning	. 11			
4	Basic principles of operation13					
	4.1	Operating and display elements				
		4.1.1 Keys				
		4.1.2 Display	. 15			
		4.1.3 Control lamps (LEDs)	. 15			
	4.2	Operating modes	. 15			
5	Ope	ration	. 17			
	5.1	Inserting and removing the reaction cells				
	5.2	Starting a temperature program				
	5.3	Stopping a temperature program				
		Settings				
	0	5.4.1 Setting the display contrast				
		5.4.2 Setting the timer for the reaction time				
6	Mai	ntenance, cleaning, disposal	. 22			
	6.1	Maintenance				
	6.2	Exchanging the fuses				
	6.3	Cleaning the enclosure				
	6.4	Cleaning the thermoblock of spilled cell contents				
	6.5	Disposal				
7	Wha	ıt to do if	25			
•	44110	it to ao il	. 20			
8	Tec	nnical Data	. 26			

Contents CR 2200

9	Replacement parts and accessories / options	7
10	Index	29

4

CR 2200 Overview

1 Overview

The CR 2200 thermoreactor is a dry temperature control device for use in the laboratory. It facilitates and safeguards the digestion with reaction cells.

8 temperature programs are permanently set.

- 1: 148 °C for 120 minutes
- 2: 120 °C for 30 minutes
- 3: 120 °C for 60 minutes
- 4: 120 °C for 120 minutes
- 5: 100 °C for 60 minutes
- 6: 148 °C for 20 minutes
- 7: 150 °C for 120 minutes
- 8: 100 °C for 30 minutes

The thermoreactor accommodates 12 reaction cells with an outer diameter of 16 mm.






You will find information on accessories in the WTW catalog LABORATORY AND FIELD INFORMATION or on the Internet.

Overview CR 2200

1.1 Components of the thermoreactor



1	Protective cover
2	Thermoblock with cell shafts
3	Sample display
4	Keypad

CR 2200 Safety

2 Safety

2.1 Safety information

2.1.1 Safety information in this operating manual

This operating manual provides important information on the safe operation of the instrument. Read this operating manual thoroughly and make yourself familiar with the instrument before putting it into operation or working with it. The operating manual must be kept in the vicinity of the instrument so you can always find the information you need.

Important safety instructions are highlighted in this operating manual. They are indicated by the warning symbol (triangle) in the left column. The signal word (e.g. "CAUTION") indicates the level of danger:



WARNING

indicates a possibly dangerous situation that can lead to serious (irreversible) injury or death if the safety instruction is not followed.



CAUTION

indicates a possibly dangerous situation that can lead to slight (reversible) injury if the safety instruction is not followed.

Note

indicates a situation where goods might be damaged if the actions mentioned are not taken.

2.1.2 Safety signs on the product

Note all labels, information signs and safety symbols on the product. A warning symbol (triangle) without text refers to safety information in this operating manual.

2.1.3 Further documents providing safety information

The following documents provide additional information, which you should observe for your safety when working with the thermoreactor:

- Analysis specifications for the cell tests
- Safety datasheets for the cell tests

Safety CR 2200

2.2 Safe operation

2.2.1 Authorized use

The authorized use is exclusively the carrying out of reactions for photometric determinations, with the aid of commercial test sets, in round cells with a diameter of 16 ± 0.2 mm. Only the operation and running of the product according to the instructions given in this operating manual is authorized (see chapter 8 TECHNICAL DATA). Any other use is considered unauthorized.

2.2.2 Requirements for safe operation

Note the following points for safe operation:

- The product may only be operated according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be operated under the environmental conditions mentioned in this operating manual.
- The round cells to be used including the screw caps must be suitable and tested for the intended application (temperature resistance, chemical stability against the reaction mixture, pressure resistance, tightness).
- The reaction mixture in the round cells must not be inflammable.
- For operation, the cell shafts have to be clean and free of obstacles.
- Round cells may only be inserted and heated with the cap screwed on.
- The instrument may only be opened if this is explicitly described in this operating manual (example: inserting and removing the cells).
- The connection cable and fuses may only be replaced with original WTW replacement parts.

2.2.3 Unauthorized use

The product must not be put into operation if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storing conditions, see chapter 8 TECHNICAL DATA).

CR 2200 Safety

2.3 Personal protective equipment

For your safety, wear the following personal protective equipment to be protected against any (residual) risk while carrying out certain activities:

- Protective goggles and protective clothing with long sleeves as a protection against any splashes of hot and corrosive liquids
- Chemical-resistant gloves as a protection against contact with possibly noxious substances
- Heat protection gloves to touch hot parts
- Safety shoes for protection against heavy falling objects

Activity	Protective goggles	Protective clothing with long sleeves	Chemical resistant gloves	Heat protection gloves	Safety shoes
Transport, positioning					*
Inserting or removing cells or the external temperature sensor	1	1		(if the instrument is hot)	
Cleaning after a cell has broken or its content has been spilled	1	1	1		

Commissioning CR 2200

3 Commissioning

3.1 General handling instructions

Always protect the instrument against conditions that might damage its mechanical and electronic components. Observe the following points especially:

- For operation and storage, the temperature and air humidity must be within the limits specified in chapter 8 TECHNICAL DATA.
- The following impacts on the instrument have to be avoided in any situation:
 - Extreme dust, moisture and wetness
 - Intensive light and heat
 - Corrosive or solvent-containing vapors.
- Any spilled or spattered liquid and broken glass in the thermoblock have to be removed immediately (see section 6.4).
- The protection cover should always be closed when the instrument is not operating.
- The thermoblock must be empty while the instrument is being transported.

3.2 Scope of delivery

- Thermoreactor CR 2200
- Power cable
- Operating manual

3.3 Unpacking and positioning



CAUTION

The instrument is quite heavy. Therefore, there is a risk of injury for the lower extremities if it falls down. Wear safety shoes during the unpacking and positioning of the instrument.



CAUTION

Fire hazard! The thermoreactor may only be operated on a surface of nonflammable material. The ventilation slots in the bottom and rear panel of the instrument must always be free. Air circulation in the area of the ventilation slots must not be impeded.

CR 2200 Commissioning

Note

Please make sure to keep the original packing including the inner packages so the instrument will be optimally protected against strong impacts if it has to be transported. The original packing is also required for the appropriate return transport if the instrument has to be repaired. Please note that the warranty does not cover any damage caused by inappropriate transport.

Note

If the thermoreactor was transported from a cold environment to a warm environment, the formation of condensate can lead to the faulty functioning of the instrument. Wait until the thermoreactor has adapted to the changed environmental conditions before putting it back into operation.


Setting up the thermoreactor

- 1 Place the thermoreactor firmly on a flat, nonflammable surface.
- 2 Make sure that there is enough space between the thermoreactor and other instruments or devices that are heat sensitive.

3.4 Initial commissioning

Adjusting the line voltage

3 Check whether the arrow on the back of the housing points to the line voltage stated on the fuse holder that is provided by the power supply system (115 or 230 V).




- 4 If the wrong line voltage is set, perform steps 5 to 7.
- 5 Pull out the fuse holder.

Commissioning CR 2200

- Turn the fuse holder so that the arrow on the housing points to the line voltage (115 or 230 V) provided by the power supply system.
- 7 Push the fusion holder completely in.

Connecting the power cable



- 1 Connection for power cable
- 2 Fusion holder




CAUTION

Each time before commissioning, check the outside of the power cable for intactness. A damaged cable must be removed and replaced by an original power cable. Original power cables are listed in chapter 9 Replacement parts and accessories / Options. Exclusively use original power cables suitable for the country-specific power supply system. Do not use any line adapters (fire hazard).



If the power cable is not suitable for the power supply system of your country, please contact your supplier for an exchange.

- 8 Connect the power cable to the socket 1 on the thermoreactor.
- 9 Connect the power cable to an easily available power socket.
 The thermoreactor is now in the Standby mode.
 The display shows the name of the instrument.



The thermoreactor is ready for operation.

4 Basic principles of operation


This chapter provides you with basic information on how to operate the thermoreactor.

4.1 Operating and display elements


Using the six keys of the keypad (see section 4.1.1) you control the thermoreactor.

Temperature values, available temperature programs or settings can be viewed in the display (see section 4.1.2).

The control lamps above the operating panel are assigned to the thermoblock. Their color (red or green) and their state (flashing or illuminated) show the current operating state of the thermoreactor (see section 4.1.3).




4.1.1 Keys



Key	Meaning
(4)	On/off switch
RUN	 Making or confirming a selection or Starting the timer for the reaction time (active temperature program).
RUN ENTER (U)	Keep depressed and simultaneously press Changing to the SETUP menu from the Standby mode
(4)	 Changing between the temperature setting and the reaction time setting (in the SETUPmenu) or Canceling the active temperature program
(A) (V)	 Selecting the temperature program (program selection) Or Changing settings and switching between settings (SETUP) Starting the scrolling through settings by keeping the key depressed

4.1.2 Display

Example: Program selection



4.1.3 Control lamps (LEDs)

The control lamps above the operating panel are assigned to the thermoblock and indicate the current operating state.

LED	flashes	is illuminated
green		Program selection
red	active temperature program: heating period or cooling period	active temperature program: reaction temperature reached
red and green	active temperature program: program canceling selected	

If the control lamps are off the thermoreactor is in the Standby mode.

4.2 Operating modes

The thermoreactor has three operating modes:

Standby

The display shows the model and version number of the thermoreactor

With the key combination and o you move to the *SETUP* menu. There you can:

- Setting the display contrast (CONTRAST:0 to 9, see section 5.4.1)
- Activating a manual confirmation before the timer for the reaction time is started (START TIMER:MAN. or AUTO, see section 5.4.2)

- Program selection
 - After switching on with (b) the display shows the current temperature value for the thermoblock.
 - In the second display line there are the temperature programs to be selected. The control lamps above the display for the thermoblock light up green.
- Active temperature program

The display shows the current temperature value for the thermoblock. The control lamps above the display light up or flash red.

CR 2200 Operation

5 Operation

5.1 Inserting and removing the reaction cells



WARNING

Always wear protective goggles and protective clothing with long sleeves as a protection against any splashes of hot and corrosive liquids. When dealing with cell tests always follow the relevant safety data sheets and analysis specifications.



CAUTION

Very hot - burn hazard. Wear heat protection gloves to remove the cells when they are hot. Place hot cells in a cell rack.

- Open the protective cover.
 Insert or remove the reaction cells from the cell shafts without jamming.
 Close the protection cover.
- i

When cold reaction cells are inserted in the preheated thermoblock it can temporarily cool down by approx. 3 °C.

5.2 Starting a temperature program

1 Switch on the thermoreactor with (b).

Select a temperature program with (a) (v).
You can select from 8 predefined temperature programs.

3 Start the displayed temperature program with ...
The control lamps of the thermoblock flash red.
The nominal reaction time (in hours and minutes) appears on the display.



The reaction temperature is reached when the temperature in the thermoblock is in a range of \pm 1 °C around the adjusted temperature for two minutes constantly. The control lamps of the thermoblock will then light up red.

Operation CR 2200

Operation with automatic timer

If the start of the timer for the reaction time has been set to automatic in the *SETUP* menu (*START TIMER:AUTO* see section 5.4.2), the reaction time automatically starts after the reaction temperature has been achieved. The reaction temperature is kept constant during the reaction time.

After the reaction time has expired the control lamps flash red.

The temperature program is finished.

The thermoreactor is in the program selection mode.

As soon as the thermoblock has cooled down to under 50 °C, the control lamps switch themselves off.

Operation with manual timer

If the start of the timer for the reaction time has been set to manual in the *SETUP* menu (*START TIMER:MAN*. see section 5.4.2), an *S* is displayed in front of the nominal reaction time. With this setting the thermoreactor controls the temperature until the timer for the reaction time is started by pressing (a).

4 Start the timer for the reaction time with . The S in front of the reaction time disappears.

The reaction temperature is kept constant during the reaction time. The control lamps of the thermoblock light up red.

After the reaction time has expired the control lamps flash red. An audio signal sounds in addition.

Confirm the end of the reaction time for the thermoblock with ...

The temperature program is finished.

The audio signal is finished.

The thermoreactor is in the program selection mode.

As soon as the thermoblock has cooled down to under 50 °C, the control lamps switch themselves off.

5.3 Stopping a temperature program

You can terminate a running program at any time.

1 Using **◆ •**), terminate the running temperature program. The control lamp for the thermoblock flashes red/green. The safety query *STOP?* is displayed.

CR 2200 Operation

54°C STOP?

2 Using confirm the safety query STOP?

The temperature program is finished.

The control lamps of the thermoblock will then light up green.

Using **④ ▶**, exit STOP?.

The query *STOP?* disappears from the display. The temperature program is continued.



While STOP? is displayed the temperature program goes on running. As soon as a section of the temperature program is finished (e.g. after the end of the heating period or after the end of the reaction time), the STOP? display is overwritten.

5.4 Settings

5.4.1 Setting the display contrast

The display contrast can be set in 10 steps.

- 1 Switch to the Standby mode.
- 2 Keep the key depressed while pressing to switch to the SETUP menu.

SETUP and, in the second line, an editable parameter are displayed.

3 Using (a) (v), select CONTRAST.

SETUP CONTRAST:5

4 Using edit the contrast setting.

The editing is marked on the display by *.

SETUP *
CONTRAST:5

Operation CR 2200

5	Using (A) (T), set the contrast from 0 to 9.
6	Using , confirm the changes. The marking on the display (*) disappears.
7	Using (b) leave the SETUP menu. The changes are stored. The thermoreactor is switched on (operation mode: program selection).

CR 2200 Operation

5.4.2 Setting the timer for the reaction time

After the start of a temperature program the thermoblock starts heating up. Depending on the setting, the timer for the reaction time starts automatically after the reaction temperature has been reached or only after confirmation by keypressing.

With the setting *START TIMER:AUTO* the timer for the reaction time starts immediately after the reaction temperature has been reached.

With the setting *START TIMER:MAN*. the timer for the reaction time only starts after confirmation by keypressing.

- Switch to the Standby mode.
 Keep the key depressed while pressing to switch to the SETUP menu.
 SETUP and, in the second line, an editable parameter or a temperature program are displayed.
 Using (), select START TIMER.
- SETUP START TIMER:MAN.

 - 5 Using (A) (7), select MAN. or AUTO.
 - 6 Using , confirm the changes.
 The marking on the display (*) disappears.
 - Using (a) leave the SETUP menu.
 The changes are stored.
 The thermoreactor is switched on (operation mode: program selection).

6 Maintenance, cleaning, disposal

6.1 Maintenance

The thermoreactor CR 2200 is maintenance free. Check the ventilation slots on the rear panel of the instrument at regular intervals. Remove any accumulations of dust in the slots, e.g. with a vacuum cleaner. Do not use compressed air!



CAUTION

Fire hazard! The ventilation slots in the bottom and rear panel of the instrument must always be free. Air circulation in the area of the ventilation slots must not be impeded.


6.2 Exchanging the fuses



CAUTION

A defective fuse may only be replaced by an original fuse with UL approval. The original fuse is listed in chapter 9 REPLACEMENT PARTS AND ACCESSORIES / OPTIONS.

1 Disconnect the line power cable from the thermoreactor.



- 2 Pull out the fuse holder.
- 3 Exchange defective fuses.
- Turn the fuse holder so that the arrow on the housing points to the line voltage (115 or 230 V) provided by the power supply system.
- 5 Push the fusion holder completely in.

6.3 Cleaning the enclosure

Wipe the thermoreactor with a damp cloth.

Note

The housing is made of synthetic material. Thus, avoid contact with acetone or detergents that contain solvents. Remove any splashes immediately.

6.4 Cleaning the thermoblock of spilled cell contents

If liquid penetrated a thermoblock (e.g. from a leaking or broken cell), clean the thermoblock as follows:



WARNING

Cells can contain poisonous or corrosive substances. If the content was released, heed the safety datasheet of the cell test and take the required protective measures (protective goggles, protective gloves etc.).



CAUTION

Danger of cuts due to broken glass. Remove any broken glass with a pair of tweezers completely. Do not pick any broken glass with your bare fingers.



CAUTION

Burn hazard on the hot thermoblock. Prior to cleaning, allow the thermoreactor to cool down.

1	Switch off the thermoreactor and disconnect the power plug.
2	Allow the thermoreactor to cool down.
3	Unscrew the cover plate on top of the thermoblock.
4	Remove any broken glass with a pair of tweezers completely.
5	Clean the cover plate, block surfaces and borings with a damp cloth.
6	Screw on the cover plate again.



CAUTION

Fire hazard. If any flammable detergents are used, completely remove all detergent residues prior to the next commissioning.



Discoloration that remains on the thermoblock and cover plate does not affect the functioning of the thermoreactor.

6.5 Disposal

Dispose of the thermoreactor as electronic waste at an appropriate collection point. It is illegal to dispose of the thermoreactor in household refuse.

CR 2200 What to do if...

7 What to do if...

There is nothing on the	Cause	Remedy
display	The power supply is interrupted	Check power cable and connections
		 Exchange the fuses
		 Have repaired by service department
Bars are displayed instead of the	Cause	Remedy
temperature (-°C)	In the program selection mode: Internal temperature sensor defective	 Have repaired by service department
Cell contents spilled /	Cause	Remedy
thermoblock contaminated	e.g. leaking cell	- see section 6.4

Technical Data CR 2200

8 Technical Data

Reactor type	Dry temperature control device with safety cover
Cell shafts	12 cell shafts for reaction cells, diameter: 16 ± 0.2 mm
Reaction time setting	20 min, 30 min, 60 min, 120 min (via fixed programs)
Temperature setting	100 °C, 120 °C, 148 °C, 150 °C via fixed programs
Controlling accuracy	± 1 °C ± 1 digit
Temperature stability	± 0.5 K
Overtemperature protection	190 °C ± 5 °C
Heating time (with empty thermoblock) from 25 °C to	100 °C approx. 5 min 120 °C approx. 7 min 148 °C approx. 10 min
Temperature of the enclosure at an environmental temperature of 25 °C	< 30 °C with a block temperature of 148 °C
Housing materials	Upper part: PC/ABS Blend Lower part: PA66 Protective cover: PC Cover plate: PBT GF30
Dimensions	D x W x H: 312 x 255 x 185 mm
Weight:	3 kg
Power supply	
	115/230 VAC 50/60 Hz, allowed voltage tolerance: ±15 % Maximum power consumption: 280 W Fuses 2 x 3,15 AT
Protective class	voltage tolerance: ±15 % Maximum power consumption: 280 W
Protective class Overvoltage category	voltage tolerance: ±15 % Maximum power consumption: 280 W Fuses 2 x 3,15 AT I according to DIN VDE 0700
	voltage tolerance: ±15 % Maximum power consumption: 280 W Fuses 2 x 3,15 AT I according to DIN VDE 0700 part 1

Climatic class	2 according to VDI/VDE 3540 Relative humidity: Yearly mean: < 75 % 30 days /year: 95 % Other days: 85 % Light dew: yes Site altitude: max. 5000 m above sea level
EMC	EN 61326-1 FCC 47 CFR Part 15
Test certificates	cETLus, CE
Safety standards	EN 61010-1 EN 61010-2-010 UL 61010-1 CAN/CSA C22.2#61010-1 CAN/CSA C22.2#61010-2-010

9 Replacement parts and accessories / options

Replacement parts

Designation	Accessory
REP3-3000 1004	Original fuse set (2 items; 3.15 A) with UL approval
REP3-3000 1006	Original power cable, type EU
REP3-3000 1007	Original power cable, type UK
REP3-3000 1008	Original power cable, type USA/Canada
REP3-3000 1009	Original power cable, type Australia

CR 2200 Index

10 Index

Commissioning	10
Components of the thermoreactor Control lamps (LEDs)	6
F	
Fuses	22
K	
Keys	14
M	
Maintenance	22
0	
Operating and display elements Operating modes	13 15
Р	
Program selection	16
S	
Safety	
Sample displayScope of delivery	
Setting the display contrast	19
StandbyStarting a temperature	15
program	17
Т	
Temperature program active Trouble shooting	16 25
W	
What to do if	25

Index CR 2200

What can Xylem do for you?

We're a global team unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to xyleminc.com.



Service address:

Xylem Analytics Germany Sales GmbH & Co. KG WTW Dr.-Karl-Slevogt-Str. 1 82362 Weilheim Germany

Tel.: +49 881 183-325 Fax: +49 881 183-414

E-Mail wtw.rma@xyleminc.com

Internet: www.WTW.com



Xylem Analytics Germany GmbH Dr.-Karl-Slevogt-Str. 1 82362 Weilheim Germany