



a xylem brand

The Photometry Dictionary

TIPS AND TRICKS FOR PHOTOMETRIC MEASUREMENTS -
FROM OUR CUSTOMERS' MAGAZINE

WATERWORLD

More AQA support for the photoLab[®] Series

wtw.com/en/photometry

photometry-compendium.com

More AQA support for the photoLab® 6000 Series

AQA 3 - The Matrix Check

Non-plausible and incorrect measuring results are a phenomenon that can occur in samples depending on the sample composition/matrix. The cause of these is often interference presented by reaction-relevant sample content substances: For example, they simulate a higher concentration due to the similar load or ion size or they remove the measuring parameter by forming complexes and thus show smaller concentration values. Classic examples from the sewer system are the interference of the ammonium determination by potassium or by nitrate measurement by chloride.

Following good laboratory practices (GLP), interferences are first detectable by the determination of a suitable standard solution, as the nominal value will not be found correctly. For the cause study search for interfering forces, the MatrixCheck is a suitable aid; it works with a so-called buildup or dilution of the sample. For this buildup, defined small volumes of a known concentration are added, sometimes in the course of several steps. As the concentration is known, you can observe, whether and how it is reflected in the measuring result. The result can also be indicated as a so-called recovery rate of the added quantity. In most cases, a recovery of < 90% and >110% points to an interference.

For a convenient and sample matrix check, photoLab® 6100 VIS and photoLab® 6600 UV-VIS offer the new function AQA 3 - matrix check for the supplementation and dilution with recovery rate:

With the commercial multi-parameter standards CombiCheck® for programmed test sets, the respective instruction values are already given via the display and menu guide; for your own methods or when using individual standards, there is also the possibility of a step-by-step buildup (max. 3).

The results can be saved and put out as a result log for the documentation.

Observation

The photoLab® 6600 UV-VIS and photoLab® 6100 VIS spectrophotometers have meanwhile been replaced by the successor models photoLab® 7600 UV-VIS and photoLab® 7100 VIS respectively. These models offer the same AQA 3 functionality.

MatrixCheck				16.04.07 9:52	
Method		1: C3/25			
Sample concentration		45 mg/l CSB			
Sample [ml]	Standard [ml]	Target value [mg/l]	nominal [mg/l]		
10	0.5	62	58	94 % ✓	
10	1	77	71		
10	1.5	91	77		
Back		Measure		Finish	

Example of a printout:

photoLab 6600	09130512	1.30-WTW-1.60	Administrator		
MatrixCheck	OK				
Protocol ID	7				
Method:	1: C3/25 COD				
Sample concentration	45 mg/l COD				
Standard ID	COD 1500				
Standard concentration	400 mg/l COD				
Sample ml	Standard ml	Target value mg/l	Actual value mg/l		
10	0.5	62	58	94%OK	
10	1	77	71	92%OK	
10	1.5	91	77	85%OK	