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
GA\_BlueLine\_pH\_8188023\_150209

## Operating instructions for pH combination electrodes from SI Analytics

### Putting into operation

The pH combination electrode are delivered ready for use. If a irrigation cap is located over the membrane and the diaphragm, please remove it. It contains potassium chloride (KCl) solution 3 mol/l. A KCl solution shortage in the electrolyte chamber of the reference system should be topped up. Topping up is unnecessary in the case of low-maintenance electrodes. With these electrodes, irrigation in electrolyte solution is particularly important.

### Calibration and measuring

For calibrating and measuring, the gate of the refilling hole has to be opened. To do so, please move the slide gate which is integrated in the electrode head into position 1. The diaphragm has to immerse into the measurement solution. Low-maintenance pH combination electrodes have a gel electrolyte, therefore they are not equipped with such a closure slide. For accurate calibration, we recommend the use of our hot-steam sterilised, certified buffer ampoules according to DIN 19 266. For detailed information, please refer to the operating instructions of the pH meter. Used buffers should not be poured back. Use fresh buffers only!

### Maintenance

- Soil on membranes and diaphragms will lead to measurement deviations. Calcareous deposits may be removed using diluted mineral acids (for instance diluted hydrochloric acid), for organic contamination please use a suitable solvent. Fats are best removed using tenside solutions, and proteins should be dissolved using hydrochloric pepsin solution (e.g. L 510)  
- After measuring and cleaning, please rinse the pH combination electrodes with distilled water. Do not rub them dry, and please do not use them as stirrers.  
- To refill and store the pH combination electrode, KCl 3 mol/l (for instance, L 300) is to be used.  
- When putting on the watering cap of the liquid electrolyte electrodes, please open the slide feed first and close it afterwards!  
- Extreme pH values and temperatures will shorten the life time of the pH combination electrodes and restrict the accuracy of the measurements.  
- Hydrofluoric acid, hot phosphoric acid and strong alkaline solutions will destroy the pH glass membrane. For further information, please refer to our website or our specialist literature.

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### SI Analytics GmbH