**Mikroprozessor Konduktometer**  
***Microprocessor Conductivity Meter***

BA51101de

**LF 538**

**Standard-Leitfähigkeitsmeßzelle**  
***Standard Conductivity Cell***

**TetraCon<sup>®</sup> 325**

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***Please read this information carefully before putting the device into service!***

*This device has been built and tested according to the IEC 1010 safety standards for electronic measuring instruments and has left our works in a condition complying with all the requirements of technical safety.*

*The perfect functioning and operational safety of the instrument can only be ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present operating instruction.*

- *Before connecting the plug-in power supply unit to the electricity supply network, it must be ensured that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).*
- *The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of these operating instructions.*
- *When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.*
- *Balancing, maintenance and repair work must only be carried out by a suitably qualified technician authorized by us.*
- *If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.*
- *The safety of the user may be endangered, e.g., if the instrument*
  - *shows visible damage,*
  - *no longer operates as specified,*
  - *has been stored over a longer period under unsuitable conditions,*
  - *has been subjected to difficult conditions during transport.*
- *If in doubt, the instrument should as a rule be sent back to the manufacturer's - "Wissenschaftlich-Technische-Werkstätten GmbH" - for repair and maintenance.*

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|                                                   |    |
|---------------------------------------------------|----|
| Measurement of total dissolved solids (TDS) ..... | 22 |
|---------------------------------------------------|----|

|                                           |    |
|-------------------------------------------|----|
| Selecting the reference temperature ..... | 23 |
|-------------------------------------------|----|

|                                                                         |    |
|-------------------------------------------------------------------------|----|
| Switch-over of measuring range selection (automatically/manually) ..... | 24 |
|-------------------------------------------------------------------------|----|

|                                         |           |
|-----------------------------------------|-----------|
| <b>Memory for measured values .....</b> | <b>26</b> |
|-----------------------------------------|-----------|

|               |    |
|---------------|----|
| Storage ..... | 26 |
|---------------|----|

|               |    |
|---------------|----|
| Readout ..... | 27 |
|---------------|----|

|               |    |
|---------------|----|
| Erasure ..... | 28 |
|---------------|----|

|                            |           |
|----------------------------|-----------|
| <b>Analog output .....</b> | <b>29</b> |
|----------------------------|-----------|

|                               |           |
|-------------------------------|-----------|
| <b>Serial interface .....</b> | <b>30</b> |
|-------------------------------|-----------|

|                                                 |    |
|-------------------------------------------------|----|
| Automatic transmission of measured values ..... | 31 |
|-------------------------------------------------|----|

|                                                 |    |
|-------------------------------------------------|----|
| Measured value transmitted by keypressing ..... | 31 |
|-------------------------------------------------|----|

|                         |    |
|-------------------------|----|
| Readout of memory ..... | 32 |
|-------------------------|----|

|                                             |           |
|---------------------------------------------|-----------|
| <b>Platinizing the measuring cell .....</b> | <b>33</b> |
|---------------------------------------------|-----------|

|                             |           |
|-----------------------------|-----------|
| <b>RESET function .....</b> | <b>34</b> |
|-----------------------------|-----------|

|                                     |           |
|-------------------------------------|-----------|
| <b>Maintenance / Cleaning .....</b> | <b>35</b> |
|-------------------------------------|-----------|

|                   |    |
|-------------------|----|
| Maintenance ..... | 35 |
|-------------------|----|

|               |    |
|---------------|----|
| Cleaning..... | 35 |
|---------------|----|

|                              |           |
|------------------------------|-----------|
| <b>Troubleshooting .....</b> | <b>36</b> |
|------------------------------|-----------|

|                            |           |
|----------------------------|-----------|
| <b>Technical data.....</b> | <b>41</b> |
|----------------------------|-----------|

|                                                            |           |
|------------------------------------------------------------|-----------|
| <b>Measuring cells and accessories for the LF 538.....</b> | <b>47</b> |
| Conductivity cells and temperature sensors .....           | 47        |
| Recommended accessories.....                               | 48        |
| Spares.....                                                | 49        |


**Instruction manual**

|                                                      |           |
|------------------------------------------------------|-----------|
| <b>Standard conductivity cell TetraCon® 325.....</b> | <b>50</b> |
| <b>Construction.....</b>                             | <b>50</b> |
| <b>Operation / Maintenance.....</b>                  | <b>51</b> |
| Cleaning.....                                        | 51        |
| Aging of the TetraCon® 325 cell .....                | 52        |
| Disposal .....                                       | 52        |
| <b>Plug arrangement.....</b>                         | <b>53</b> |
| <b>Technical data .....</b>                          | <b>55</b> |

|                                                                                                           |           |
|-----------------------------------------------------------------------------------------------------------|-----------|
| <b>Application report: "Control of conductivity systems<br/>acc. to ISO 9000, GLP requirements" .....</b> | <b>60</b> |
|-----------------------------------------------------------------------------------------------------------|-----------|

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


*For power supply use only the Friemann and Wolf Line Adapter, model FW1199.*

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
*Before the instrument is plugged in and switched on, it has to be controlled if the working voltage of the instrument and line voltage are corresponding!*

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*Connect power-supply unit with electric circuit and the instrument.*

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



*If a linear temperature function has been selected, the following appears:*

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*Display of the temperature coefficient in case of selection of linear temperature function, e.g. 2.18 % K.*

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

bis / until

LF-Messung läuft kontinuierlich.  
*Conductivity measurement runs continuously.*

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*No valid measured value available.*

#### Reason:

- Temperature out of specified measuring range.
- No or defective measuring cell connected.
- Measuring cell without temperature sensor connected (e.g. KLE 325). Such cells can only be used together with a TFK 325 as an extra temperature sensor.

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



*Fixed cell constant C=0.1 $\text{cm}^{-1}$  selected.*

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
*Switch back to freely adjustable cell constant.*

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*Back to measurement.  
Display of measured value and  
temperature of test sample.*

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*Further information see application report no. 395054.*


**Temperaturfunktion TC einstellen**  
**Setting of the temperature function TC**

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\* *Correction of the intrinsic conductivity of water by the ultrapure water function works also with measurements with linear temperature coefficient.*

• International oceanographic tables.

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Statusanzeige: nLF

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
*Temperature function **nLF** selected.*

*Status indication: nLF*

*Indication of measured value and **nLF**.*

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

*No temperature compensation selected. Status indication: Lin  
Display of measured value and temperature coefficient 0.00 %/K.*

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
*Linear temperature function selected. Status indication: Lin  
Display of measured value and currently set temperature coefficient.*

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*Back to measurement.  
Display of measured value and temperature of sample solution.*

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
*Further information see application report no. 1194055.*

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*Rinse measuring cell and immerse it into sample solution.*

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*Salinity measurement runs continuously.*

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*Further information see application report no. 993109.*

### Messen des Filtrattrockenrückstands (TDS) Measurement of total dissolved solids (TDS)

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

*Measurement of total dissolved solids*

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
*Set determined TDS factor.  
TDS measurement runs continuously.*

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*Further information see application report no. 995084.*

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

Switch over reference temperature  
TREF.  
Press key TC during switching on:  
reference temperature  
e.g. 25°C set.

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
Switch over reference temperature  
TREF.  
Press key TC during switching on:  
reference temperature  
e.g. 20°C set.

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
**Automatic** selection of measuring range ARng ("AutoRange") **on**.  
Press key during switching on.  
Status indication: **ARng**  
Instrument measures automatically in measuring range with highest possible resolution.

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

**Automatic** selection of measuring range ARng ("AutoRange") **off**.  
Press key during switching on.

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**Manual** selection of measuring range by repeated pressing the key .

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
*General information on temperature compensation at the conductivity measurement:  
Application report no. LF 1194055*

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Memory capacity: Max. 50 pairs of measuring values (conductivity + temperature)

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
*Storage is only possible in operation mode LF, Sal or TDS.*

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
*38th storage location.  
Mem: storage function switched on.*

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
*After occupying of the 50th storage location the instrument displays Full Sto before each storage.*

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
*RUN/ENTER overwrites beginning at 1 an already occupied storage location.*

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
*Each other key prevents the storage.*

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
*Readout is only possible in LF, Sal or TDS operation mode.*

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
*903 µS/cm: stored value.  
 04: storage location.  
 Mem: storage function.*

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*After approx. 2 sec.:  
 903 µS/cm: stored value.  
 21.2°C: stored temperature value.  
 Mem: memory function.*

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*Select storage location (for readout only). For fast run keep key depressed.*

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*If no value is stored, no reaction results from RCL.*

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*Depress key STO during switching on.*

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*Clr Sto is not displayed if memory is already empty.*

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
*Every other key cancels the erasure process.*

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*Connect analog output to recorder with cable.*

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
**Resolution:** 1 mV corresponds to the specific resolution of the measuring range.

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*Connect serial output of the instrument with serial port PC/Printer.*

PC: Cable AK 325/B, O. no. 902 836  
Printer: Cable AK 325/S, O. no. 902 837

Connections: 1 CTS  
2 plug codification  
3 ground  
4 TxD

*Connection of 2 and 3 activates serial output.*

*Adjust transmission settings at the PC/Printer.*


Baud rate: 4800  
Parity (PC only): none  
Handshake: RTS/CTS Handshake  
Data bits (PC only): 8  
Stop bit (PC only): 1

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*Serial interface ON.  
Switch over from OFF to ON:  
Keep key X depressed, switch on  
instrument.*

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
*The instrument automatically transmits the data via the RS 232 interface. Renewal of  
measured value approx. every 3 up to 5 s.*

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*Serial interface OFF.  
Switch over from ON to OFF:  
Keep key X depressed, switch on  
instrument.*

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
*Instrument transmits value pair  
e.g. 903 µS/cm and 21.2°C by  
keypressing via RS 232 interface.*

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
Readout of memory selected.  
(Indication: Mem)

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Start data transmission. Readout from oldest to newest value.

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Data are transmitted in ASCII-format, the data sets are separated with CR, LF.

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#### Example for screen terminal program:

| Datei      | Bearbeiten | Einstellungen | Telefon | Übertragung |
|------------|------------|---------------|---------|-------------|
| 1902       | µS/cm      | 25.0°C        |         |             |
| Lin TREF25 |            |               |         |             |
| Sa1        | 000.9      | 25.0°C        |         |             |
| nLF        |            |               |         |             |
| TDS        | 1901       | mg/l          | 25.0°C  |             |
| nLF TREF25 |            |               |         |             |
| 01         | 1902       | µS/cm         | 25.0°C  |             |
| Lin TREF25 |            |               |         |             |
| 02         | 190.1      | µS/cm         | 25.0°C  |             |
| Lin TREF25 |            |               |         |             |
| 03         | 099.8      | mS/cm         | 25.0°C  |             |
| nLF TREF25 |            |               |         |             |
| 04         | 089.2      | mS/cm         | 25.0°C  |             |
| nLF TREF20 |            |               |         |             |
| 05         | Sa1        | 005.7         | 25.0°C  |             |
| nLF        |            |               |         |             |

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
Further information see application report no. 895235.

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*Connect measuring cell to platinizing output.*

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*Display shows **PLA 1** for negative platinizing current (starting condition).*

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*Positive platinizing current is switched on. Display shows **PLA 2**.*

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*Negative platinizing current is switched on. Display shows **PLA 1**.*

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
*Further information see application report no. 1195243.*

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
Set back:

- cell constant to  $0.475 \text{ cm}^{-1}$  (TetraCon 325).
- temperature compensation to nLF
- reference temperature to  $T_{REF} 25$

Anzeige:

Zellenkonstante  $0,475 \text{ cm}^{-1}$ ,  
 $T_{REF} 25$  und nLF.

Display:

Cell constant  $0.475 \text{ cm}^{-1}$ ,  $T_{REF} 25$   
and nLF.

## **Wartung / Maintenance**

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*The instrument is maintenance-free!*

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*Cleaning with aqueous soap solutions.*

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| Error message                          | Cause                                                 | Action                              |
|----------------------------------------|-------------------------------------------------------|-------------------------------------|
| <b>FL</b><br>Indication range exceeded | The measured value is not within the measuring range. | Select next-higher measuring range. |
| <b>to</b>                              | Timeout of the serial interface.                      | Check connected instrument.         |

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Additional measuring range with cell constant  
0.1  $\text{cm}^{-1}$ :

Range

Resolution

|                                 |                       |
|---------------------------------|-----------------------|
| 0.00 ... 19.99 $\mu\text{S/cm}$ | 0.01 $\mu\text{S/cm}$ |
|---------------------------------|-----------------------|

Accuracy:  
(Temp. of test sample  
=  $T_{REF}$ )

$\pm 0.5\%$  of measured value  $\pm 1$  digit at operating  
temperature  $-10 \dots +55^\circ\text{C}$

Cell constant

Setting ranges: 0.1  $\text{cm}^{-1}$  fixed  
0.45 ... 1.30  $\text{cm}^{-1}$  adjustable

Reference temperature

$20^\circ\text{C}$  or  $25^\circ\text{C}$  selectable.

Temperature  
compensation

automatically, 3 modes selectable:

1. Linear compensation ( $\overline{\text{Lin}}$ ),  
coefficient adjustable between 0.5 %/K and  
3.00 %/K in the range of  $-5 \dots 90^\circ\text{C}$ .
2. Compensation of natural water ( $\overline{\text{nLF}}$ )  
acc. to EN 27888 (replaces DIN 38404),  
Ultrapure water function for neutral water;  
Coefficient fixed.
3. No compensation.

Accuracy:  
(at operation temp.  
 $-10^\circ\text{C} \dots +55^\circ\text{C}$ )

1. Linear compensation ( $\overline{\text{Lin}}$ ):  
0.5 % of measured value  $\pm 1$  digit for  
TC:  $\alpha = 0.5 \dots 2.5\%$ /K and test sample  
temperature  $+5^\circ\text{C} \dots 75^\circ\text{C}$  and  
TC:  $\alpha = 2.5 \dots 3.00\%$ /K and test sample  
temperature  $10^\circ\text{C} \dots 75^\circ\text{C}$ .
2. Non-linear compensation ( $\overline{\text{nLF}}$ ):  
0.5 % of measured value  $\pm 1$  digit for test  
sample temperature  $0 \dots 35^\circ\text{C}$ .  
For test sample temperature  $35 \dots 50^\circ\text{C}$   
expanded nLF-function acc. to WTW-  
measurements.
3. No compensation:  
0.5 % of measured value  $\pm 1$  digit for test  
sample temperature  $-5 \dots +99.9^\circ\text{C}$ .

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                                                                                                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><u>Salinity</u></p> <p>Range: 0.0 ... 70.0<br/>Resolution: 0.1<br/>Accuracy: after adjustment (via cell constant) to 35.0:</p> <ul style="list-style-type: none"> <li>· <math>\pm 0.1</math> in the range 0.0 ... 42.0 at test sample temperature 5 ... 25°C.</li> <li>· <math>\pm 0.2</math> in the range of 0.0 ... 42.0 at test sample temperature 0 ... 5°C resp. 25 ... 30°C</li> <li>· better than 3 % of measured value in the range 42.0 ... 70.0.</li> </ul> <p><u>TDS</u></p> <p>Range: 0 ... 1999 mg/l<br/>Resolution: 1 mg/l<br/>Accuracy: <math>\pm 2 \% \pm 1</math> digit at test sample temperature 15°C ... 35°C.<br/>Factor adjustable between 0.40 ... 1.00.</p> <p><u>Temperature</u></p> <p>Range: -5 ... 99.9°C<br/>Resolution: 0.1 K<br/>Accuracy: 0.1 K <math>\pm 1</math> digit at operating temperature -10°C ... +55°C</p> | <p>according to IOT chart<br/>source: "UNESCO Technical Papers in marine sciences 36"</p>                                                                                                                                                                        |
| <p><b>Measuring rate</b></p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | <p>Renewal of measured value: approx. every 1.2 sec<br/>With active serial interface: approx. every 3.5 sec.</p>                                                                                                                                                 |
| <p><b>Data storage</b></p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | <p>Ring memory<br/>Number of values: 50 pairs<br/>Output: Display, RS 232 interface<br/>Data protection: any long</p>                                                                                                                                            |
| <p><b>Display</b></p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | <p>WTW-specific LCD-display<br/>Height: 35 mm<br/>Width: 60 mm</p>                                                                                                                                                                                               |
| <p><b>Keyboard</b></p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | <p>Field with 10 keys</p>                                                                                                                                                                                                                                        |
| <p><b>Inputs</b></p> <p>Temperature sensor</p> <p>Socket for power supply</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | <p>8 pin WTW special socket</p> <p>2 x 4 mm jack<br/>TFK 325 as separated temperature sensor (only together with conductivity cells without integrated temperature sensor)</p> <p>2 pin special socket for line adaptor<br/>Voltage: 12 V<br/>Current: 70 mA</p> |

| <b>Outputs</b>                            |                                                                                                                                                                                                                                                                                                                                                       |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>Platinizing output</i>                 | 8 pin WTW special socket<br>Impressed direct current 6 mA<br>Polarity can be chosen with key<br>Accuracy: + 1 mA with operating temperature<br>-10°C ... +55°C                                                                                                                                                                                        |
| <i>Analog output</i>                      | 2 x 4 mm jack for recorder,<br>not galvanically isolated<br>Signal range: 0 ... 1.999 V for range voltage<br>0 ... 1999 digit<br><br>Accuracy:<br>$\pm 0.5 \%$ of displayed value $\pm 0.1$ mV<br>at ambient temperature -10°C ... +55°C<br><br>Source resistance:<br>$< 5 \text{ W}$ (Current limit<br>max. 0.2 mA output current)                   |
| <i>Digital interface</i><br><i>RS 232</i> | 4 pin special socket for<br>RS 232 cable AK 325/B or<br>Printer cable AK 325/S<br>Serial interface, not galvanically isolated<br>Data output, no separate excitation<br>Baud rate: 4800<br>Data bits: 8<br>Stop bit: 1<br>Parity: none<br>Handshake: RTS/CTS Handshake<br>Max. cable length: 15 m<br>Output format: ASCII-strings<br>separator CR, LF |
| <b>Environmental conditions</b>           |                                                                                                                                                                                                                                                                                                                                                       |
| <i>Ambient temperature</i>                | Operation: + 5°C ... + 40°C<br>Storage: - 25°C ... + 65°C                                                                                                                                                                                                                                                                                             |
| <i>Relative humidity</i>                  | Yearly mean: < 75 %<br>30 days/year: < 95 %<br>Remaining days: < 85 %                                                                                                                                                                                                                                                                                 |
| <b>EMC</b>                                |                                                                                                                                                                                                                                                                                                                                                       |
| <i>Interference transmission</i>          | Generic standard EN 50081-1<br>FCC class A                                                                                                                                                                                                                                                                                                            |
| <i>Disturbance</i>                        | Generic standard EN 50082-2<br>Namur recommendation (standard requirements)                                                                                                                                                                                                                                                                           |
| <b>Power supply</b>                       |                                                                                                                                                                                                                                                                                                                                                       |
|                                           | Line adapter<br><br>230 VAC: Friwo FW1199/11.7864,<br>Friwo Part. No. 1762613<br>120 VAC: Friwo FW1199/11.7880,<br>Friwo Part. No. 1794043                                                                                                                                                                                                            |

### **Data storage and operation mode in case of power failure**

*Settings and measuring values remain stored.*

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■ EMV-Prüfzentrum ■ EMI/EMC-Testcenter ■

## CERTIFICATE

### to Electromagnetic Compatibility

to Test Report No. 52501-5403

**EUT:** LF 538  
**Description:** LF-Meter  
**Applicant:** WTW GmbH  
**Regulations:** prEN 50082-2:1/1993  
 Namur Recommendation (standard requirements) 12/1990

#### Testresult:

The tested sample is in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

#### RFI Emissions

EN 50081-1:1992 See test report no. 52501-5142-1 for OXI 538  
 FCC rules part 15/B See test report no. 52501-5142-2 for OXI 538

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
| <b>Contamination</b>                         | <b>Cleaning solution</b>                                 | <b>Reaction time at room temperature</b> |
|----------------------------------------------|----------------------------------------------------------|------------------------------------------|
| Water-soluble substances                     | Deionized water                                          | any                                      |
| Grease and oil<br><br>at heavy contamination | Warm water and household cleaning solution<br><br>Spirit | any<br><br>max. 5 minutes                |
| Lime and hydroxide coatings                  | Acetic acid (10 %)                                       | any                                      |

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*Fundamentally, the conductivity measuring cell does not age. The cell life is considerably shortened or the cell damaged by excessive temperatures or special measuring solutions (e.g. strong acid and lye solutions, organic solvents). We give no warranty for defects and mechanical damage caused by the measuring medium.*

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
*We recommend disposal as electronic waste.*

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
- number of electrodes 4
- electrode material *graphite*

**Shaft material***Epoxy, black***Dimensions**

- shaft diameter 15.3 mm
- shaft length 120 mm
- diameter of connection head 21.7 mm
- total length of the sensor 162.5 mm
- cable length 1.5 m (*cable is exchangeable by WTW-Service*)

**Immersion depth**

*min. 36 mm*  
*max. total length incl. cable (J < 80°C)*  
*max. 120 mm (J < 100°C)*

**Pressure resistance***2 bar***Connections**

- cable, plug *protection system IP67, 8-pin plug*
- measuring cell *protection system IP68 (20 m)*

**Cell constant***0.475 cm<sup>-1</sup> ± 1.5 %***Range of application***1 µS/cm .. 2 S/cm***Temperature measurement**

- thermistor type *integrated NTC (30 kW / 25°C)*
- thermistor material *graphite*
- temperature range *-5 .. +100°C*
- response time of thermistor *t<sub>99</sub> < 20 sec.*
- sensor accuracy *± 0,2 K*

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## Control of conductivity systems acc. to ISO 9000; GLP - requirements

---

Replaces reports 53023,  
293054, 594054 and 994054

As part of the introduction of a quality system according to standard series DIN ISO 9000, a growing number of users find themselves obliged to carry out test equipment checks for the measuring systems applied. In this context, the system of rules of the GLP (Good Laboratory Practice) results in almost identical requirements.

To enable you to meet the ISO 9000 and GLP requirements, it is necessary to establish the following conditions:

- The measuring systems applied must be calibrated at regular intervals in order to ensure the accuracy and comparability of the measuring data acquired.
- The completion of the specified calibrations must always be recorded in writing: type and identity number of measuring instrument and measuring cell, date, signature of the person carrying out the calibration.

It is advisable to have a separate log form or log book for each measuring system (instrument with measuring cell) where the calibrations are entered consecutively. This documentation should also state the calibration intervals, so that the user can see at any time when the next calibration is due.

Consequently, it will no longer be necessary to attach calibration labels on the instrument and the measuring cell.

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---

LF 1194054

- Calibration intervals and instructions for carrying out calibrations are to be specified in writing:

The appropriate calibration interval depends on the conditions of application and must be fixed by the operator on his own responsibility according to his experience.

The procedure for carrying out the calibration is laid down in standard EN 27888. It is therefore sufficient to refer to this standard as instructions for calibration.

If desired, we can also check the instrument and measuring cell for you at our works and issue a corresponding test certificate. You will nevertheless still need to carry out regular calibration.

## Calibration

### Standard solutions

Use potassium chloride solutions according to the specifications of standard EN 27888

| Concentration<br>in mol/l | Conductivity at 25°C<br>in mS/m                  µS/cm |       |
|---------------------------|--------------------------------------------------------|-------|
| 0,0005                    | 7,4                                                    | 74    |
| 0,001                     | 14,7                                                   | 147   |
| 0,005                     | 72                                                     | 720   |
| 0,01                      | 141                                                    | 1413  |
| 0,02                      | 277                                                    | 2770  |
| 0,05                      | 670                                                    | 6700  |
| 0,1                       | 1290                                                   | 12900 |
| 0,2                       | 2480                                                   | 24800 |

Disturbing gases in water for preparing potassium chloride solutions with concentrations  $c < 0.01$  mol/l must be removed by degassing the water with nitrogen or heating it up. The conductivity should be less than 2 µS/cm. The solutions have to be freshly prepared before use.

Solutions with low potassium chloride concentrations ( $c < 0.01$  mol/l) are very unstable. During calibration one can expect larger deviations of the measured values from the theoretical values.

## Calibration procedure

- Select a reference temperature of 25°C on the measuring vessel.
- Rinse the measuring cell and the measuring vessel several times with some of the standard solution.
- Fill the measuring vessel with standard solution.
- Stabilize the temperature of the standard solution and the measuring solution at 25°C. For high-accuracy measurements it is necessary to observe a tolerance of  $\pm 0.1$  K, while  $\pm 0.5$  K will be sufficient for routine measurements.
- Compare the measured value with the theoretical value of the standard solution.

With WTW plastic cells, the measured value should not deviate from the theoretical value by more than 1.5 % (2 % in the case of older measuring cells).

Undesirably low values are caused by measuring cells:

- with dirty electrodes. Clean the electrodes.
- with platinized electrodes if the surface has dried out. Replatinize the electrodes.
- with polarized electrodes. Cells with a cell constant of  $k < 0.5/\text{cm}$  already polarize if the conductivity amounts to  $> 200 \mu\text{S}/\text{cm}$ . Use a standard solution with lower conductivity.

## Adjustment

The adjustment of a conductivity meter is effected by setting of the cell constant. This may, e.g., be necessary with a new cell or after platinizing of the electrodes. The conductometer must have a potentiometer or a corresponding function for setting of the cell constant.

With new measuring cells it will often be sufficient to enter the value stated on the cell cable or the operating instructions.

For high-accuracy measurements involving cells with a constant of  $k > 0.1/\text{cm}$  it is better to carry out a determination of the cell constant.

## Procedure

Use a potassium chloride solution  $c = 0.01 \text{ mol/l}$  as standard solution.

- Select a reference temperature of  $25^{\circ}\text{C}$  on the measuring vessel.
- Rinse the measuring cell and the measuring vessel several times with some of the standard solution.
- Fill the measuring vessel with standard solution.
- Stabilize the temperature of the standard solution and the measuring solution at  $25^{\circ}\text{C}$ . For high-accuracy measurements it is necessary to observe a tolerance of  $\pm 0.1 \text{ K}$ , while  $\pm 0.5 \text{ K}$  will be sufficient for routine measurements.
- Set the cell constant until the instrument shows a conductivity of  $1413 \mu\text{S/cm}$ .

### Note

It is risky to try compensating errors due to dirty, dried out or polarized electrodes by adjusting the cell constant.

### Literature

EN 27888, Determination of Electric Conductivity, European Standard, 9/93