

Operating manual

D. O. probe TriOxmatic[®] 700IN

Accuracy when	The use of advanced technology and the high	
going to press	quality standard of our instruments are the result of a	
	continuous development.	
	This may result in differences between this operating	
	manual and your instrument. Also, we cannot guarantee that	
	there are absolutely no errors in this manual. Therefore, we	
	are sure you will understand that we cannot accept any legal	
	claims resulting from the data, figures or descriptions.	


Copyright © Weilheim 2000, WTW GmbH Reprinting - even as excerpts - is only allowed with the explicit written authorization of WTW GmbH, Weilheim. Printed in Germany.


Contents

1	Overview		
	1.1	Overall view of the D. O. probe TriOxmatic [®] 700IN .5	
	1.2	Construction of the D. O. probe $TriOxmatic^{ entropy}$ 700IN . 6	
2	Saf	ety7	
	2.1	Authorized use	
	2.2	General safety instructions	
3	Со	mmissioning11	
	3.1	Scope of delivery11	
	3.2	Commissioning / Readiness for measuring 11	
	3.3	Recommended fields of application11	
	3.4	Calibration11	
	3.5	Measuring12	
	3.6	Exterior cleaning and function check	
	3.7	Cleaning with normal operation	
		(e.g. domestic sewage)14	
	3.8	Storing	
4	Mai	intenance, cleaning, disposal, spare parts 17	
	4.1	Replacing the electrolyte solution and the	
	4.0	membrane head	
	4.2	Cleaning the electrodes	
		4.2.2 Cleaning the counter electrode in the	
		glass beaker	
	4.3	5	
		Checking whether the probe is zero current free 32	
	4.5	Disposal	
	4.6	Maintenance and spare parts	
5	Wh	at to do if	
6	Red	commended accessories	
7	Тес	hnical data41	
	7.1	General data	
	7.2	Measuring conditions	
	7.3	-	
	-	Electrical data	

1 Overview

1.1 Overall view of the D. O. probe TriOxmatic[®] 700IN





1.2 Construction of the D. O. probe TriOxmatic[®] 700IN

1	Protection cap (for transport and storage)
2	Protection hood
3	Membrane head, WP 600
4	Shaft
5	Closing head
6	Protection ring for G1 thread
7	Working electrode (platinium cathode)
8	Isolator
9	Counter electrode (silver anode)

2 Safety

This operating manual contains basic instructions that you must follow during the commissioning, operation and maintenance of the probe. Consequently, all responsible personnel must read this operating manual before working with the probe. The operating manual must always be available within the vicinity of the probe.

Target groupThe D. O. probe was designed for use with online
instrumentation - mostly in waste water treatment plants.
Thus, we assume that, as a result of their professional
training and experience, the operators will know the
necessary safety precautions to take when handling
chemicals.

Safety instructions In the various chapters of this operating manual, the following safety instructions point out the different degrees of danger:



Warning

indicates safety instructions which must be adhered to precisely to rule out possible serious hazards to persons.



Caution

indicates safety instructions which must be adhered to precisely to avoid slight injuries, damage to your equipment or environmental harm.

Other Notes



Note

indicates notes that draw your attention to special features.



Note

indicates cross references to other documents, e.g. operating manuals.

2.1 Authorized use

The technical specifications as given in chapter 7 TECHNICAL DATA must be observed. Only the operation and running of the probe according to the instructions given in this operating manual is authorized.

Any other use is considered **unauthorized**.

2.2 General safety instructions

The probe left the factory in a safe and secure technical condition.

Function and operating safety

The smooth functioning and operational safety of the probe can only be guaranteed if the generally applicable safety measures and the specific safety instructions in this operating manual are followed during operation.

The smooth functioning and operational safety of the probe can only be guaranteed under the environmental and electrical operating conditions that are specified in chapter 7 TECHNICAL DATA.

The temperature specified for using the probe (chapter 7 TECHNICAL DATA) has to be observed. Protect your probe especially against frost or overheating.



Caution

The probe is only allowed to be opened by personnel authorized by WTW.

Safe operation	If safe operation is no longer possible, the probe must be taken out of service and secured against inadvertent operation!
	Safe operation is no longer possible if the probe
	 has been damaged in transport
	 has been stored under adverse conditions for a lengthy period of time
	 is visibly damaged
	 no longer operates as described in this manual.

If you are in any doubt, please contact the supplier of the probe.

Obligations of the
purchaserThe purchaser of this probe must ensure that the following
laws and guidelines are observed when using dangerous
substances:

- EEC directives for protective labor legislation
- National protective labor legislation
- Safety regulations
- Safety datasheets of the chemical manufacturers.

Safety

3 Commissioning

3.1 Scope of delivery

- TriOxmatic[®] 700IN with accessory kit ZBK 600 (with 1.5 m cable length: protection ring for the G1 thread included)
- The probes are filled with electrolyte and are equipped with a protection hood and a protection cap.

3.2 Commissioning / Readiness for measuring

1	Connect the probe to the measuring instrument.
2	Take off the gray protection cap.
3	Der Sensor ist nach einer Verzögerungszeit von 2 Minuten meßbereit.
4	Calibrate.



Note

Recommendation for precision measurements:

Recalibrate the next day or let the probe polarize for the night.

3.3 Recommended fields of application

On-line measurements in rivers, lakes and waste water and in activated sludge basins of waste water treatment plants.

3.4 Calibration



Note

For calibrating, please refer to the operating manual of the measuring instrument.

3.5 Measuring

Provide the necessary minimum flow (chapter 7 TECHNICAL DATA), e.g.:

- the flow velocity of the water is sufficient
- pull the probe slowly through the water by hand
- the turbulence in the activated sludge basin is sufficient.

3.6 Exterior cleaning and function check

The function check is the simplest way of letting you know whether the probe needs to be cleaned and calibrated.


Note

The function check can either be performed in water vaporsaturated air or in air-saturated water.

In the case of air temperatures under 5 °C we recommend to perform the function check not in air but in air-saturated water that has a higher temperature.

Function check in water vaporsaturated air Proceed as follows:


Take the probe out of the measuring solution and position it approx. 2 cm above a water surface, best in a narrow bucket or similar container with water. When doing so, no liquid should be brought on the membrane.





Wait at least 15 min for adjustment.

5 Switch the measuring instrument to the % saturation operating mode. Read the measured value and determine whether it lies within the range of precision required by the user.





Note

If the measured value is within the accuracy tolerance required by the user it is not necessary to clean or recalibrate the probe.

Example:

Required accuracy: 5 % Measured value in air: 97 % saturation \Rightarrow no cleaning or recalibrating required (the measuring error is within the required accuracy tolerance)

Note:

If the temperature difference between the measuring solution and the air above the water surface is too high, an adjusting time of more than 15 min may be necessary.

Function check in air-saturated water

In the case of air temperatures under 5 °C do not perform the function check in air but in air-saturated water that has a higher temperature. You obtain air-saturated water by pouring water several times in and out of two vessels so that it sparkles.

Then immserse the probe in a vessel with air-saturated water and perform the function check as described above.

3.7 Cleaning with normal operation (e.g. domestic sewage)

Cleaning and calibrating is required:


- approx. every 4 to 6 weeks
- if the measured values are too low

Contamination	Cleaning procedure
Sludge and removable sediment	Warm tapwater with washing- up liquid, a smooth cloth or a smooth sponge
Salt and/or lime sediment	Acetic acid (w = 25 %), a smooth cloth or a smooth sponge


Thoroughly rinse with deionized water.

Exterior cleaning

1 Screw off the protection hood.



2 **Carefully** clean the probe shaft and the membrane.








Caution Clean the membrane very carefully!

The membrane may be damaged by sharp-edged objects or when lying on stones.


3 Clean the protection hood as well.









5 Rinse the probe with clear tapwater.



6 Carefully dry the membrane using a lint-free paper towel.





Note

probe.

For calibrating the probe refer to the relevant operating manual of the industrial monitor.

Wait at least 15 minutes for adjustment, then calibrate the



Note

Maintenance of the probe is required (see chapter 4 MAINTENANCE, CLEANING, DISPOSAL, SPARE PARTS) if:

- the membrane is damaged (leakage display)
- the probe cannot be calibrated any more.

3.8 Storing

Store the probe with the gray protection cap put on.

Temperature: -5 °C to +50 °C Position: any

4 Maintenance, cleaning, disposal, spare parts

4.1 Replacing the electrolyte solution and the membrane head

WTW delivers the probe ready for operation

An exchange is necessary if:

- the membrane is damaged (leakage display)
- the membrane is strongly contaminated; the probe cannot be calibrated
- the electrolyte solution is depleted.



Warning

Contact with the measuring solution can endanger the user! Take protective measures suitable for the kind of measuring solution (protective clothes, eye protectors etc.).




Warning

The ELY/A electrolyte solution is strongly alkaline!

- Keep out of the reach of children!
- Wear suitable gloves and eye/face protection!
- After contact with eyes and skin rinse immediately with plenty of water!
- After contact with the eyes seek medical advice!

ELY/A




R: 36/38 S: 2-26-27-37/39




Caution

For all maintenance work, disconnect the probe from the instrument!


1 Disconnect the probe from the instrument.




2 Screw off the membrane head.




3 Thoroughly rinse the probe head with deionized water.



4 Carefully wipe and dry the counter electrode using a paper towel.




5 Moisten the new membrane with some drops of the ELY/A electrolyte solution.





Warning Observe the warning on the ELY/A bottle! 6 Fill the membrane head with electrolyte solution.







Note

Always use a new membrane head!


7 Remove air bubbles by careful knocking.




8 Rinse the probe head with electrolyte solution.




9 Screw the WP 600 membrane head on the shaft.
 Hold the probe inclined. Excess electrolyte solution is pressed out at the vent.



10 The filling should be free of air bubbles. Small air bubbles do not interfere.



- 11 Refill only in case of:
 - Big air bubbles
 - Luftblasen an der Platinkathode
 - Measuring with high water pressure
 - Probe cannot be calibrated or is sensitive to knocking.





Note

After being connected to the instrument for approx. 60 minutes the probe is ready for operation, then **it is essential to recalibrate the probe**.

Recommendation:

For precision measurements, calibrate the probe the following day once again.

4.2 Cleaning the electrodes




Note


The counter electrode always shows a coloring. This is required for the operation of the probe. It is no contamination.

Cleaning is only required if the slope is too high or too low (the probe cannot be calibrated), and if this cannot be remedied by replacing the membrane head and the electrolyte solution.


1 Disconnect the probe from the instrument.



2 Screw off the membrane head.




3 Thoroughly rinse the probe head with deionized water.



4.2.1 Cleaning of the platinium working electrode

Moisten the cathode and the SF 300 grinding foil with deionized water.

1 Using the **wet** SF 300 grinding foil, polish impurities off the cathode under slight pressure.






Caution


Do not use normal grinding paper or glass fibre pen!

2 Thoroughly rinse the probe head with deionized water.




4.2.2 Cleaning the counter electrode in the glass beaker


1 Screw off the membrane head.




2 Using a lint-free paper towel, wipe the counter electrode and carefully remove any loose coating.



3 Fill the glass beaker (150 mL, high shape) with approx. 25 mL RL-AG/Oxi cleaning solution.



4 Put the probe in the stand. Immerse the electrode unit into the RL-AG/Oxi cleaning solution up to a little bit above the counter electrode.






Note

At this point the cleaning procedure deviates from the common maintenance procedure for WTW D. O. probes! This procedure applies for the TriOxmatic[®] 700IN only.




Allow the cleaning solution to react for 1 hour.

5 Thoroughly rinse the probe head and the electrode unit with deionized water several times.




6 Immerse the probe head and the electrode unit into deionized water.





Soak in deionized water for at least one hour.

7 Carefully shake off water drops.






Note

See section 4.3 REFILLING THE MEMBRANE HEAD.

4.3 Refilling the membrane head

1 Moisten the new membrane with some drops of the ELY/A electrolyte solution.






Warning

Observe the warning on the ELY/A bottle!

2 Fill the membrane head with electrolyte solution.







Note

Always use a new membrane head!


3 Remove air bubbles by careful knocking.



4 Rinse the probe head with electrolyte solution.



5 Screw the WP 600 membrane head on the shaft. Hold the probe inclined. Excess electrolyte solution is pressed out at the vent.



6 Connect the probe to the measuring instrument again. The probe is ready to operate after being connected to the instrument for approx. 60 minutes.



Note

After each refilling, it is essential to recalibrate the probe.

Recommendation:

For precision measurements, calibrate the probe the following day once again.

4.4	Checking whether the probe is zero current
	free

The probe is zero current free. Therefore, it is not possible to calibrate the probe in air free of oxygen.

It is possible, however, to check whether the probe is zero current free according to DIN EN 25814/ISO 5814. This is required in

Test solution 1 g/L aqueous sodium sulfite solution, Na₂SO₃ (adding 1 mg/L of a cobaltous salt accelerates the removal of oxygen from the solution.)

case of a suspected malfunction.

The probe should have been in operation for at least one hour before the check.

Course of the check

- Perform an exterior cleaning of the probe.
- Switch the measuring instrument to the % saturation operating mode.
- Immerse the probe into the test solution.
- Wait for one hour.



Note

Leave the probe in the test solution for one hour.

Test criterion

The probe is OK when the measuring instrument displays < 2 % oxygen saturation.

If the measured value is > 2 %:

- Reinigen der Platin-Arbeitselektrode und Wechsel von Elektrolyt und Membrankopf.
- If necessary, also clean the counter electrode (see section 4.2.1 CLEANING OF THE PLATINIUM WORKING ELECTRODE to section 4.3 REFILLING THE MEMBRANE HEAD)

4.5 Disposal

We recommend to dispose of the probe as electronic waste.

4.6 Maintenance and spare parts

	Model	Order no.
Exchange membrane heads (2 pieces)	WP 600/2	202 850
Electrolyte solution (1 bottle with 50 mL)	ELY/A	205 212
Cleaning solution for the silver counter electrode (1 bottle with 50 mL)	RL-Ag/Oxi	205 200
Grinding foil	SF 300	203 680
Accessory kit	ZBK 600	202 620

5 What to do if ...

The measuring instrument displays 0.0 mg/L or 0 % Sat O₂ (the probe is in air)

Cause	Remedy	
 Measuring instrument and probe are not connected to each other 	 Check the connection between measuring instrument and probe 	
 No electrolyte in the membrane head 	 Replace the WP 600 membrane head 	
 Cable defective 	 Send probe to WTW 	
 Electronics defective 	 Send probe to WTW 	

Probe cannot be	Cause	Remedy
calibrated	 Membrane head contaminated 	 Perform an exterior cleaning of the probe according to the operating manual. Then wait for at least 15 minutes and recalibrate. If there is contamination that cannot be removed: Replace the electrolyte solution and the membrane head

Probe still cannot be	Cause	Remedy
calibrated after the electrolyte and the membrane head have been replaced	 The electrodes are contaminated or the probe is poisoned 	 Clean the electrodes according to the operating manual-

The probe has been mechanically damaged	Cause	Remedy Send probe to WTW
Measured values too low	Cause Membrane contaminated 	Remedy Perform an exterior
	 Probe has not been calibrated for a lengthy period of time 	cleaning and calibration according to the operating manual
	 Membran nicht an Platin- Arbeitselektrode anliegend 	 Replace the membrane head and calibrate according to the operating manual
	 Shunt resistance caused by moisture in the probe, plug or cable 	 To check, disconnect the probe, screw off the membrane head, thoroughly rinse the electrode unit with distilled water and completely dry it Connect the probe to the instrument and switch to the highest measuring range (0 - 600 %): If the display is > 1 digit: There is a shunt. Send the probe to WTW for repair.

The measured value strongly fluctuates	Cause	Remedy
	 Membran nicht an Platin- Arbeitselektrode anliegend 	 Replace the membrane head and calibrate according to the operating manual
	 There is a tiny hole in the membrane, directly in front of the working electrode 	 Replace the membrane head and recalibrate according to the operating manual, if necessary clean the electrodes
Measured values are	Cause	Remedy
too high	 The has not yet completely polarized 	 Wait for the polarization to be completed (2 h)
	 Probe has not been calibrated for a lengthy period of time 	 Perform an exterior cleaning and calibration according to the operating manual
	 Shunt resistance caused by moisture in the probe, plug or cable 	 To check, disconnect the probe, screw off the membrane head, thoroughly rinse the electrode unit with distilled water and completely dry it Connect the probe to the instrument and switch to the highest measuring range (0 - 600 %): If the display is > 1 digit: There is a shunt. Send the probe to WTW for repair.

The leakage display appears	Cause	Remedy	
	 The membrane head is loose 	 Screw the membrane head tight. 	
	 There is a hole in the membrane 	 Replace the membrane head and calibrate according to the operating manual 	

Wrong temperature	Cause	Remedy
display	 Defective temperature sensor 	 Send probe to WTW

6 Recommended accessories

Model	Order no.	
EBST 700-DU	203 716	
EBS 700-DU	203 715	
D 702	203 738	
EH/U 170	109 320	
UA 55	109 260	
UA-SO	109 261 109 263 V	
EH/F 170-1.5	109 272	
EH/F 170-2.5	109 273	
S 200	108 540	
	EBST 700-DU EBS 700-DU D 702 EH/U 170 UA 55 UA130 UA-SO EH/F 170-1.5 EH/F 170-2.5	

Electrode weight GWA

109 232

(only for floater mounting or extremely strong current and/or turbulence)



Note

For further holders and simulators refer to the WTW-catalog.

	7	Technical data	
	7.1	General data	
Probe types	700I 700I 700I	xmatic [®] N-1.5 N-7 N-15 N-SO	Cable length 1.5 m, 7 m, 15 m customized lengths of up to 100 m
Measuring principle	Membrane covered amperometric probe with a 2 electrodes system potentiostatically operated; electrode with integrated preamplifier		
Electrolyte	ELY,	/Α	
Output signal	Voltage signal, proportional to the O ₂ partial pressure (corrected by temperature and water vapor partial pressure)		
Temperature measurement	via integrated NTC		
Temperature compensation	IMT		

Material	Membrane head:	POM
	Membrane:	Fluoroplastics
	Probe head:	POM
	Isolator:	Epoxy, PEEK
	Closing head:	POM
	Shaft:	VA steel 1.4571
	Protection hood:	POM
	Cable screw joint:	VA steel 1.4571
	Cable coating:	PU (Polyurethane)

Dimensions	Shaft length	193 mm
	Shaft diameter	40 mm
	Membrane thickness	50 µm
	Cable length	1.5 m, 7 m, 15 m customized lengths up to 100 m

Weight (1.5 m 675 g cable included)

7.2 Measuring conditions

Measuring range	0 - 60 mg/L O ₂ 0 - 600 % saturation (≈ 0 - 1200 mbar pO ₂)		
Polarization time	bei Elektrolytwechsel:at least. 60 minbei kurzen Polaristionsunterbrechungen:2 min		
Temperature range	•	0 °C + 50 °C - 5 °C + 50 °C	
Max. admissible overpressure	10 bar		
Immersion depth	min. 10 cm; max. 100 m depth		
Connection	7-pin plugConnectionProbe - cable:IP 68 (100 m)Cable - plug:IP 65Plug - instrument:IP 65		
Operating position	any, if the minimum flow is granted		
Flow to the probe	\geq 0.5 cm/s for 5 % measuring accuracy \geq 5 cm/s for < 1 % measuring accuracy		
Fields of application	Water and waste water control		

7.3 Characteristics when delivered

Slope	- 3.75 mV/mbar; mean value at 20 °C	
Zero signal	< 0.2 % of the saturation value	
Response time at 20 °C	t_{90} (90 % of the final displayed value after)< 180 s	
Own consumption at 20 °C	0.0059 µg h⁻¹ (mg/L)⁻¹ at 20 °C	
Drift	< 1 % per month with permanent polarization	
Operational lifetime	approx. 5 years per electrolyte filling	

7.4 Electrical data

Pin assignment



Top view of the plug

Pin	Assignment	Voltage	Max. current
1	Ub+	+ 6.5 V + 8.5 V	< 1mA
2	0V		
3	NTC (30kOhm)	< 3.5 V	< 0.15mA
4	NTC (30kOhm)	< 3.5 V	< 0.15mA
5	Ub-	- 6.5 V 8.5 V	< 1 mA
6	Reg/Leak	- 8.5 V + 8.5 V	
7	(output signal) pO ₂ (output signal)	0 V + 8.5 V	



Warning

- All voltages must be protective low voltages or safety extra low voltages without shock-hazard according to IEC 1010 or UL3101.
- All currents to the probe must not exceed 8 A even in case of a short-circuit.