(English version see page 2)
(Version français voir page 3)
(Versión en español vea página 4)

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Readme File for Firmware Update of Turb 430

General information:

The program "Firmware Update Turb 430" enables you to **update the firmware of the Turb 430** to the latest version with the aid of a PC. To do so, you need a free port on your PC and a suitable interface cable (see detailed operating manual on CD-ROM or Internet).

Program installation:

With the installation program "Turb430_Vx_yy.exe" you install the firmware update program on your PC.

Requirements:

- PC with current operating system Microsoft Windows
- Administrator rights for the installation.

Program start:

In the Windows start menu, execute the program "Update" in the "Xylem Analytics/Update_Turb430" directory.

The program automatically selects the first free interface. The selected interface is displayed on the left-hand side of the status bar at the bottom edge of the window. Via the Language menu, you can change the language (German or English).

Note:

For Microsoft Windows 7, administrator rights are required.

Firmware update:

Using a suitable cable (see detailed operating manual on CD-ROM or Internet), connect the Turb 430 to be programmed to the interface of the PC displayed in the status bar.

Make sure the instrument to be programmed is switched on.

Click the OK button to start the update process.

The program automatically recognizes the instrument type connected.

Proceed by following the instructions of the program.

During the programming process, a corresponding message and a continuation display (in %) are displayed.

The programming process takes about 8 minutes. After successful programming, a final message appears. The programming process is finished now. Disconnect the instrument from the PC. The instrument is ready for operation. After switching Off/On you can check whether the new software version was taken over in the start display of the instrument.

Note:

Instrument functions can change due to a firmware update.

The current operating manual is available on the Internet under www.xylemanalytics.com.

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